

OXIDATION-REDUCTION REACTIONS OF  
MORPHINE IN DIFFERENT MEDIA

By

CHINDA WONGWIECHINTANA

Master of Science  
Oklahoma State University  
Stillwater, Oklahoma  
1977

Bachelor of Science  
Mahidol University  
Bangkok, Thailand  
1974

Submitted to the Faculty of the  
Graduate College of the  
Oklahoma State University  
in partial fulfillment of  
the requirements for  
the Degree of  
DOCTOR OF PHILOSOPHY  
December, 1983

Thesis  
1983D  
W8720  
cop. 2



OXIDATION-REDUCTION REACTIONS OF  
MORPHINE IN DIFFERENT MEDIA

Thesis Approved:

*W. C. Purditt*

Thesis Adviser

*Warren T. Ford*

*Tom E. Moore*

*L. Herbert Breenece*

*Norman D. Durham*

Dean of the Graduate College

## ACKNOWLEDGMENTS

I wish to express my sincere appreciation to my major adviser, Dr. Neil Purdie, for his invaluable guidance, suggestions, patience and encouragement throughout my graduate work. Special thanks to Dr. E. M. Holt for her guidance, effort, patience, and endless amount of time spent with me throughout my X-ray study. Without their help this study would have been impossible.

I also wish to thank the members of my faculty committee, Dr. W. T. Ford, Dr. T. E. Moore and Dr. L. H. Bruneau for their understanding, encouragement, and helpful suggestions. A special note of thanks is given to Dr. K. D. Berlin for his counsel and advice in the NMR study.

A special thanks is extended to Mr. Norman Perreira and Mr. Stan Sigle for mass spectrometry and NMR services, to all the faculty and staff of the chemistry department and to my fellow graduate students especially Mr. Dennis Fine, Dr. John Bowen for their friendship and encouragement and Dr. Nasar Pourahmady for his guidance.

I wish to thank my parents, Mr. and Mrs. Wong, and my brothers and sisters for their understanding and moral support throughout my study and my life.

Finally, I would like to thank the Chemistry Department, the National Science Foundation and Dow Chemical Company for the assistantships and fellowships during my study.

## TABLE OF CONTENTS

Chapter	Page
I. INTRODUCTION . . . . .	1
Statement of the Problem . . . . .	5
II. BACKGROUND AND THEORY . . . . .	6
Circular Dichroism . . . . .	6
Nuclear Magnetic Resonance . . . . .	12
Information Available from NMR Spectroscopy . . . . .	20
X-ray Diffraction . . . . .	24
III. EXPERIMENTAL . . . . .	30
Instruments . . . . .	31
Chemicals . . . . .	32
Experimental Procedures . . . . .	35
Oxidation Reaction of Morphine in Basic Solution . . . . .	35
Oxidation Reaction of Morphine by Metal Ions in Aqueous Solution . . . . .	35
Oxidation of Pseudomorphine . . . . .	36
Oxidation of Morphine and Apomorphine in Strong Acid Solution . . . . .	36
X-ray Structural Analysis . . . . .	38
Solid State Circular Dichroism Study of Morphine Methyliodide . . . . .	40
IV. EXPERIMENTAL RESULTS . . . . .	42
Oxidation of Morphine in Basic Solution . . . . .	42
NMR Study of Pseudomorphine . . . . .	42
Oxidation of Pseudomorphine in Base . . . . .	53
Oxidation of Morphine and Pseudomorphine by Metal Ions in Neutral Solution . . . . .	66
Oxidation of Morphine and Apomorphine in Strong Acid Solutions . . . . .	71
Morphine Sulfate . . . . .	71
Apomorphine Hydrochloride . . . . .	81
Thin-layer Chromatography . . . . .	90
X-ray Structural Analysis . . . . .	90

Chapter	Page
Structural Comparison of Morphine Sulfate and Morphine Methyl-iodide . . . . .	90
The Structure of Compound I . . . . .	109
Solid State Circular Dichroism of Morphine Methyliodide . . . . .	109
V. DISCUSSION . . . . .	126
Oxidation of Morphine in Basic Solution . . . . .	126
Oxidation of Morphine and Pseudomorphine by Inorganic Ions in Aqueous Solutions . . . . .	131
Oxidation of Morphine Sulfate and Apomorphine Hydrochloride in Concentrated Sulfuric Acid Solution . . . . .	134
X-ray Structural Analysis . . . . .	137
Structural Comparison of Morphine Sulfate and Morphine Methyl-iodide . . . . .	137
The Structure of Compound I . . . . .	144
Conclusion and Suggestions for Further Work . . . . .	146
BIBLIOGRAPHY . . . . .	150
APPENDICES . . . . .	153
APPENDIX A - X-RAY DATA OF MORPHINE SULFATE . . . . .	154
APPENDIX B - X-RAY DATA OF MORPHINE METHYL- IODIDE . . . . .	178
APPENDIX C - X-RAY DATA OF COMPOUND I . . . . .	195

# LIST OF TABLES

Table		Page
I.	Reagents Used in Acid Oxidation Reaction of Morphine . . . . .	33
II.	Reagents Used in Neutral and Basic Reactions of Morphine . . . . .	34
III.	List of Color Test Reagents . . . . .	37
IV.	$^{13}\text{C}$ NMR Chemical Shifts of Morphine Sulfate and Pseudomorphine in DMSO ( $d_6$ ) . . . . .	50
V.	$^{13}\text{C}$ NMR Chemical Shifts of Phenol and $o,o'$ -Biphenol in DMSO ( $d_6$ ) . . . . .	62
VI.	$^{13}\text{C}$ NMR Chemical Shifts of Anisole and 3,3'-Bianisole in DMSO ( $d_6$ ) . . . . .	63
VII.	CD Data of Morphine, Pseudomorphine and HOP in Distilled Water . . . . .	65
VIII.	UV-Visible Data of Morphine Sulfate in Concentrated Sulfuric Acid . . . . .	75
IX.	Color Change Data of Apomorphine in Concentrated Sulfuric Acid . . . . .	81
X.	Crystal Data for $(\text{C}_{17}\text{H}_{20}\text{NO}_3^+)_2\text{SO}_4 \cdot 3\text{H}_2\text{O}$ . . . . .	92
XI.	Crystal Data for $(\text{C}_{18}\text{H}_{22}\text{NO}_3^+)\text{I}^- \cdot \text{H}_2\text{O}$ . . . . .	93
XII.	Nonhydrogen Positional Parameters for $(\text{C}_{17}\text{H}_{20}\text{NO}_3^+)_2\text{SO}_4$ . . . . .	94
XIII.	Hydrogen Positional Parameters for $(\text{C}_{17}\text{H}_{20}\text{NO}_3^+)_2\text{SO}_4$ . . . . .	95
XIV.	Nonhydrogen Atom Anisotropic Thermal Parameters for $(\text{C}_{17}\text{H}_{20}\text{NO}_3^+)_2\text{SO}_4$ . . . . .	96
XV.	Hydrogen Isotropic Thermal Parameters for $(\text{C}_{17}\text{H}_{20}\text{NO}_3^+)_2\text{SO}_4$ . . . . .	98



Table	Page
XVI. Nonhydrogen Positional Parameters for (C <sub>18</sub> H <sub>22</sub> NO <sub>3</sub> ) <sup>+</sup> I <sup>-</sup> . . . . .	99
XVII. Nonhydrogen Atom Anisotropic Thermal Parameters for (C <sub>18</sub> H <sub>22</sub> NO <sub>3</sub> ) <sup>+</sup> I <sup>-</sup> . . . . .	100
XVIII. Hydrogen Positional Parameters for (C <sub>18</sub> H <sub>22</sub> NO <sub>3</sub> ) <sup>+</sup> I <sup>-</sup> . . . . .	101
XIX. Bond Lengths(A <sup>o</sup> ) and Angles( <sup>o</sup> ) for (C <sub>17</sub> H <sub>20</sub> NO <sub>3</sub> <sup>+</sup> ) <sub>2</sub> SO <sub>4</sub> <sup>=</sup> . . . . .	102
XX. Bond Lengths(A <sup>o</sup> ) and Angles( <sup>o</sup> ) for (C <sub>18</sub> H <sub>22</sub> NO <sub>3</sub> ) <sup>+</sup> I <sup>-</sup> . . . . .	104
XXI. Torsional Angles of Morphine and Morphine Salts . . . . .	106
XXII. Crystal Data for C <sub>17</sub> H <sub>13</sub> NO <sub>3</sub> . . . . .	114
XXIII. Nonhydrogen Positional Parameters for C <sub>17</sub> H <sub>13</sub> NO <sub>3</sub> . . . . .	115
XXIV. Hydrogen Positional Parameters for C <sub>17</sub> H <sub>13</sub> NO <sub>3</sub> . . . . .	116
XXV. Nonhydrogen Atom Anisotropic Thermal Parameters for C <sub>17</sub> H <sub>13</sub> NO <sub>3</sub> . . . . .	117
XXVI. Hydrogen Isotropic Thermal Parameters for C <sub>17</sub> H <sub>13</sub> NO <sub>3</sub> . . . . .	118
XXVII. Bond Lengths(A <sup>o</sup> ) and Angles( <sup>o</sup> ) for C <sub>17</sub> H <sub>13</sub> NO <sub>3</sub> . . . . .	119
XXVIII. SSCD and CD Data of Morphine and Its Salts . . . . .	125
XXIX. Reduction Potential Data of Metal Ions in the Neutral Redox Reactions . . . . .	133

## LIST OF FIGURES

Figure	Page
1. End View of Electric Field Vector $E$ as the Resultant of Two Rotating Vectors, $E_L$ and $E_R$ . . . . .	8
2. End View of Electric Vectors on Passage of Plane-Polarized Light Through an Optically Active Material . . . . .	8
3. End View of Electric Field Vectors on Passage of Plane-polarized Light Through an Optically Active Material . . . . .	11
4. Energy Levels and Transitions for a Nucleus ( $I=1/2$ ) in a Magnetic Field . . . . .	15
5. Proton Processing in a Magnetic Field ( $H_0$ ) . . . . .	16
6. Oscillator Generates Rotating Component of Magnetic Field $H_1$ . . . . .	17
7. X-ray Spectra with Characteristic Peaks from a Molybdenum Anode X-ray Tube Operated at 35 kV . . . . .	26
8. The Bragg Scattering Condition . . . . .	27
9. $^1H$ NMR Spectrum of Pseudomorphine in DMSO ( $d_6$ ) . . . . .	43
10. $^1H$ NMR Spectrum of Morphine Sulfate in DMSO ( $d_6$ ) . . . . .	44
11. $^{13}C$ NMR Spectra of Pseudomorphine in DMSO ( $d_6$ ): a) Fully Decoupled, and b) Off-resonance . . . . .	46
12. The Unsaturated Region Expansion of the Pseudomorphine $^{13}C$ Off-resonance NMR Spectrum . . . . .	47
13. $^{13}C$ NMR (Fully Decoupled) Spectrum of Morphine Sulfate in DMSO ( $d_6$ ) . . . . .	48

Figure	Page
14. $^{13}\text{C}$ NMR (Off-resonance) Spectrum of Morphine Sulfate in DMSO (d6) . . . . .	49
15. $^{13}\text{C}$ NMR (Off-resonance) Spectrum of Codeine in DMSO (d6) . . . . .	51
16. $^{13}\text{C}$ NMR (Fully Decoupled) Spectrum of Codeine in DMSO (d6) . . . . .	52
17. $^{13}\text{C}$ NMR (Fully Decoupled) Spectrum of Phenol in DMSO (d6) . . . . .	54
18. $^{13}\text{C}$ NMR (Off-resonance) Spectrum of Phenol in DMSO (d6) . . . . .	55
19. $^{13}\text{C}$ NMR (Fully Decoupled) Spectrum of o,o'-Biphenol in DMSO (d6) . . . . .	56
20. $^{13}\text{C}$ NMR (Off-resonance) Spectrum of o,o'-Biphenol in DMSO (d6) . . . . .	57
21. $^{13}\text{C}$ NMR (Fully Decoupled) Spectrum of Anisole in DMSO (d6) . . . . .	58
22. $^{13}\text{C}$ NMR (Off-resonance) Spectrum of Anisole in DMSO (d6) . . . . .	59
23. $^{13}\text{C}$ NMR (Fully Decoupled) Spectrum of 3,3'-Bianisole in DMSO (d6) . . . . .	60
24. $^{13}\text{C}$ NMR (Off-resonance) Spectrum of 3,3'-Bianisole in DMSO (d6) . . . . .	61
25. CD Spectra of a) Morphine in Distilled Water, b) Pseudomorphine in Distilled Water, and c) Oxidation Product of Pseudomorphine in Distilled Water . . . . .	64
26. Structures of Morphine, Apomorphine and Pseudomorphine . . . . .	67
27. CD Spectra of Solutions from $\text{Ag}^+$ -morphine at Mole-ratio 3:1 at Different Times . . . . .	68
28. CD Spectra of Solutions from $\text{Ag}^+$ -Morphine Reaction at Different Times . . . . .	69
29. CD Spectra of Solutions from Morphine in Excess of $\text{K}_3\text{Fe}(\text{CN})_6$ in Distilled Water . . . . .	70

Figure	Page
30. CD Spectra of Solutions Containing 1:10 Mole Ratio of Pseudomorphine to a) $\text{Cu}^{+2}$ 2 Days After Mixing; b) $\text{Ag}^{+}$ 2 Days After Mixing; c) $\text{KMnO}_4$ After Mixing; d) $\text{K}_3\text{Fe}(\text{CN})_6$ After Mixing . . . . .	72
31. UV Spectra of Morphine Sulfate in Concentrated Sulfuric Acid.	
32. Visible Spectra of Morphine Sulfate in Concentrated Sulfuric Acid . . . . .	74
33. CD Spectra of Morphine Sulfate in Concentrated Sulfuric Acid . . . . .	76
34. UV-visible Spectra of the Chloroform Extracts of the Solutions Containing Morphine Sulfate with a) Froehde Reagent, b) Mandelin Reagent, and c) Mecke Reagent . . . . .	78
35. UV-visible Spectrum of the Chloroform Extract of the Solution Containing Morphine Sulfate and Concentrated Sulfuric Acid . . . . .	79
36. Mass Spectrum of the Chloroform Extracted Solution from Morphine Sulfate-Froehde Reactions . . . . .	80
37. UV Spectrum of Apomorphine Hydrochloride in Concentrated Sulfuric Acid (After Mixing) . . . .	82
38. UV Spectrum of Apomorphine Hydrochloride in Distilled Water . . . . .	83
39. UV Spectrum of Apomorphine Hydrochloride in Concentrated Sulfuric Acid After 2 Weeks . . . .	84
40. Visible Spectrum of Apomorphine Hydrochloride in Concentrated Sulfuric Acid After 2 Weeks . . .	85
41. CD Spectra of Apomorphine Hydrochloride in Concentrated Sulfuric Acid . . . . .	87
42. CD Spectra of Apomorphine Hydrochloride in Distilled Water . . . . .	88
43. UV-visible Spectra of the Chloroform Extracts of a) Apomorphine hydrochloride with Concentrated Sulfuric Acid, b) with Froehde Reagent . . . . .	89

Figure	Page
44. An ORTEP Projection of Morphine Sulfate . . . . .	107
45. An ORTEP Projection of Morphine Methyliodide . . . . .	108
46. An ORTEP Projection of Hydrogen Bonding Network of Morphine Sulfate . . . . .	110
47. An ORTEP Projection of Hydrogen Bonding Network of Morphine Methyliodide . . . . .	111
48. An ORTEP Projection of Hydrogen Bonding Network of Morphine Hydrate . . . . .	112
49. Mass Spectrum of Compound I . . . . .	113
50. An ORTEP Projection of Compound I . . . . .	120
51. An ORTEP View of Hydrogen Bonding of Compound I . . . . .	121
52. SSCD Spectrum of Morphine (...), Morphine Sulfate (—) and Morphine Methyliodide (---) . . . . .	123
53. CD Spectra of Morphine (...), Morphine Sulfate (—), and Morphine Methyliodide (---) in Distilled Water . . . . .	124
54. The Possible Structures of Pseudomorphine . . . . .	128
55. Structures of Phenol, o,o'-Biphenol, Anisole and m,m'-Bianisole . . . . .	130
56. Postulated Mechanism of Formation of Compound I . . . . .	138
57. The Numbering System of Compound I . . . . .	145

## CHAPTER I

### INTRODUCTION

Morphine, the naturally occurring alkaloid or nitrogenous base, is found principally in opium, which is the dried sap of the unripe seed capsule of Papaver somniferum and is present to about 10% by weight of the total alkaloid extract. Other important alkaloids are codeine, thebaine and papaverine which together make up 30% by weight (1). Morphine was isolated from opium being the first "vegetable base" reported by Serturmer in 1805 (2). Methods of isolation, structural determination, synthesis and physical properties of morphine have been summarized by Bentley and Holmes (3, 4).

Morphine is very readily oxidized, reducing gold and silver salts to the metals (5). Mild oxidation of morphine by a variety of reagents, especially in alkaline solutions, and including atmospheric oxygen in ammonia (6) condenses together two molecules of morphine to produce a dimer referred to most commonly as pseudomorphine. The molecule is also less frequently referred to as oxymorphine. oxydimorphine, dehydromorphine and pormine (7).

Pseudomorphine was first described in 1835 by Pelletier (8) in an investigation of the alkaloids

obtained from opium, and its characteristic properties were studied. It was noted that pseudomorphine is physiologically inert and undergoes further oxidation easily. In opium it is present to only 0.02% by weight (9) and is most likely produced from the partial oxidation of morphine. No evidence for further oxidation in the plant is available. There is no adequate proof that pseudomorphine is a metabolite of morphine in human organisms (7). The chemistry and the structure are not well known. It has been studied only in a desultory way and over a long period of years. However, it was reported as a first oxidation product of morphine (10). It was proposed by Bentley and Dyke (11) that the structure of the product was a 2,2'-dimer. Small and Turnbull (12) on the other hand suggested that it was a 1,1'-dimer. This ambiguity has not been resolved. Only two physical methods were applied to the identification of pseudomorphine, namely UV and fluorescence spectroscopy (11, 13).

Morphine is classified as a controlled or illegal substance. Its abuse is of international concern. Rapid identification is critical in pharmaceutical, clinical and forensic science applications. The formation of pseudomorphine has been suggested as a confirmatory test for the identification of morphine (14). Besides the formation of pseudomorphine, there are several other physical and chemical methods employed in the identification and/or determination of morphine e.g. UV-visible (15), IR (16),

and NMR spectroscopy (17, 18, 19), chromatography (20, 21, 22), circular dichroism (23), microcrystalline tests (24) and color spot tests (25). Color spot tests are the area of interest to this study.

Color Tests are used as a preliminary screen to place an anonymous substance into a particular classification group which expedites its identification. None of the tests is specific for any drug in that members of the same group will react similarly with the test reagent to produce a consistent color. The color test method was introduced by Dragendorff (26) in 1868 and was developed for the identification of drugs by Umberger (27). Most of the color test reagents used for opium alkaloids require concentrated sulfuric acid as the reaction solvent and some oxidative metal salts. It was reported by Fulton in 1928 (28) that morphine develops a series of colors in the sulfuric acid solution itself. In the same paper, he also postulated that the product was not apomorphine as most authorities had intimated (p. 24):

The statements in the textbooks as to the compound formed on heating morphine in sulfuric acid seem to be nothing more than poor guesses.

These very empirical color tests have been used although little is known about the chemistry involved including the mechanisms and the identities of the products formed. Possible interfering reactions have been identified; however, in 1964, Schieser (29) described an ESR study of a number of color test solutions for alkaloids,



and specifically for morphine sulfate. His results suggested that the reactions involved free radicals and the changes in color with time were due to the formation of different free-radical anions which appeared concurrently with each color. He also suggested that the site of radical formation in the alkaloid is the aromatic ring moiety. No reaction mechanism was proposed. In 1975, Ahlers and Auterhoff (30) used UV-visible, IR spectroscopy and thin-layer chromatography to identify the products of the morphine after reaction with the color test reagents developed by Froehde, Mandelin and Erdmann. They concluded that the product is an o-quinone of apomorphine. Again the chemistry and the complex mechanisms of these reactions are as yet unexplained.

The metals in the inorganically based color reagents are for the most part in their highest oxidation state, so reduction of these complex salts by the alkaloid is implied. Accordingly, this study covers the oxidation of morphine under two different conditions: the mild conditions pertaining to the alkaline medium and the more severe conditions which prevail in concentrated sulfuric acid. A variety of experimental and instrumental techniques have been applied to the study of these processes with a view to revealing additional information on the products so that the complexities of the reactions might be better interpreted.

## Statement of the Problem

This research was undertaken to study the redox reactions of morphine in basic, neutral and strong acid solutions. The chemistry and mechanisms of the reactions and the structures of the products in these reactions are of interest.

The particular physical methods used in the study are circular dichroism (CD), nuclear magnetic resonance (NMR) and X-ray crystallography. The background and theory of these techniques are discussed in the following chapter. Details of the reactions and structural information on the products are presented in subsequent chapters.

## CHAPTER II

### BACKGROUND AND THEORY

#### Circular Dichroism

Circular dichroism (CD) is a useful spectroscopic technique when applied to structural determinations of compounds having one or more optically active sites. Besides being optically active the compounds which are CD active must also contain a chromophore to absorb electromagnetic radiation. As such CD is simply a modification of absorption spectrophotometry. The basic difference is that the incident light is circularly polarized and the experimental parameter measured is the difference in absorption of the left and right circular components. It is most amenable to studies of molecules in solution, but has been applied to a few solids supported in KBr pellets.

To understand the basic theory of CD, one should begin with the discussion of the properties of light. Ordinary light behaves as if it consisted of a large number of electromagnetic waves oscillating in all possible orientations around the direction of propagation. The electric and magnetic fields associated with the propagating

beam oscillate at right angles to each other and perpendicular to the beam direction. By passing the beam through special devices known as polarizer elements, only waves oscillating in a particular plane are transmitted. Devices are available for the infrared, visible and ultraviolet ranges of the electromagnetic spectrum. The transmitted beam is commonly referred to as plane-polarized light. Considering only the electric vector  $E$  of the plane-polarized beam, it can be represented as consisting of two circularly polarized components moving in vacuo at the same speed. The vectors are imagined to rotate in opposite directions at the same angular velocity determined by the refractive index ( $n$ ) of the medium. Their vectors can be represented by  $E_L$  (left rotation) and  $E_R$  (right rotation) as shown in Figure 1. So long as the angular velocities are equal the vector sum is  $E$  oscillating in the direction of the original plane of polarization. If the angular velocities differ,  $\omega \neq \omega'$ , then the vector sum is no longer directed along the original plane but is rotated out of that plane by an angle  $\alpha$ .

For Figure 2  $\omega > \omega'$  and rotation is to the right. The only reason angular velocities would differ is that the medium through which the beam passes would have a different refractive index for each polarization, i.e.  $n_L \neq n_R$ . When this occurs the medium is said to be circular birefringent.  $E$  lies in the direction of the




Figure 1. End View of Electric Field Vector  $E$  as the Resultant of Two Rotating Vectors,  $E_L$  and  $E_R$




Figure 2. End View of Electric Vectors on Passage of Plane-Polarized Light Through an Optically Active Material

diagonal line of the square whose sides are  $E_L$  and  $E_R$ .

Since  $n_L$  and  $n_R$  are wavelength dependent,  $\alpha$  varies as a function of wavelength. The spectrum of  $\alpha$  plotted versus the wavelength  $\lambda$  is called Optical Rotatory Dispersion (ORD). The relationship was described mathematically by Fresnel (32) according to equation (2.1).

$$\alpha = \pi/\lambda(n_L - n_R) \quad (2.1)$$

where:  $\alpha$  has units of radians per unit length and  $\lambda$  is the incident light wavelength.

Media which are circularly birefringent are optically active. For solutions the activity is a consequence of dissolved chiral molecules. The property  $\alpha$  is, as expected, dependent upon the solution concentration. To normalize  $\alpha$  for comparison among solutions of different concentrations, a quantity known as the specific rotation  $[\alpha]$  is defined, equation (2.2).

$$[\alpha] = 1800 * \alpha / c * \pi \quad (2.2)$$

where:  $c$  is the concentration in grams per milliliter of a solution, and  $1800/\pi$  is the conversion factor needed to give the specific rotation in degrees per decimeter.

To allow for comparisons among different compounds in solutions of different concentrations, the more general molecular rotation  $[\Phi]$  term is defined, equation (2.3).

$$[\Phi] = [\alpha] * M/100 \quad (2.3)$$

where:  $M$  is the molecular weight of the optically active substance.

If the optically active medium contains a compound which absorbs energy from the beam then not only the speed of rotation will change but also the compound will absorb left and right circularly polarized light to different extents. Therefore in the spectral regions where optically active absorption bands are present, the length of vector  $E_R$  is no longer equal to  $E_L$ , and their resultant  $E$  no longer oscillates along the circumference of a circle. Instead, the resultant vector  $E$  now traces out an ellipse, as shown in Figure 3.

The emitted light is said to be elliptically polarized and the medium is said to exhibit circular dichroism. In terms of quantum mechanics, this means that the transition probabilities are different for left and right circularly polarized light (33). The resulting ellipse is characterized by the angle of ellipticity,  $\theta$ , which is given by the following equation (2.4).

$$\theta = \pi * (A_L - A_R) / \lambda \quad (2.4)$$

where:  $A_L$  and  $A_R$  are the absorbances for the left and right circularly polarized components, respectively. The  $A$  values are related to the molar absorbances,  $\epsilon$ , and the molar concentration  $C$  of the absorbing substance by equation (2.5).

$$A = \epsilon * b * d \quad (2.5)$$

By analogy with specific rotation, the specific ellipticity  $[\theta]$  is defined by equation (2.6).




Figure 3. End View of Electric Field Vectors on Passage of Plane-polarized Light Through an Optically Active Material. The wavelength of light is within an absorption band.



$$[\theta] = \theta/lc \quad (2.6)$$

where:  $\theta$  is measured in degrees,  $l$  is the path length in decimeters and  $c$  is the concentration in grams per milliliter of solution.

Similarly, the molar ellipticity  $[\Psi]$  is defined by equation (2.7).

$$[\Psi] = [\theta]M/100 = \theta/Cd = 3000 (\epsilon_L - \epsilon_R) \quad (2.7)$$

where:  $d$  is the path length in centimeters and  $C$  is the molar concentration of absorbing species.

CD spectra are obtained by measuring the ellipticity  $\theta$  as a function of wavelength  $\lambda$ . The spectra may show positive and/or negative peak(s). The wavelengths of the maxima are referred to as  $\lambda_{\max}^+$  or  $\lambda_{\max}^-$ . These maxima and the wavelength(s) of crossover point(s) ( $\lambda^0$ ) are characteristics of each CD active compound and can be used for qualitative distinction. All CD curves, by definition, are Cotton effect curves which are named in honor of the French physicist, Aime Cotton, who discovered this phenomenon in 1896 (34).

Since CD is a modified form of absorption spectroscopy, the method is not only useful in compound identification but also in quantitative analysis, kinetic studies, and stereochemical conformational assignments.

### Nuclear Magnetic Resonance

Nuclear Magnetic Resonance spectroscopy (NMR) is a

branch of chemical spectroscopy which concerns radio-frequency (rf) induced transition between quantized energy states of magnetic nuclei that have been oriented by magnetic fields.

A fundamental property of an atomic nucleus besides mass and charge is the nuclear spin quantum number or nuclear spin ( $I$ ) which has values  $0, 1/2, 1, 3/2, \dots$  etc. in units of  $h/2$ , where  $h$  is Planck's constant. The actual value of the nuclear spin of any given nucleus depends on the mass number and atomic number. The empirical rules relating the mass number and atomic number to nuclear spin are:

1. If both the mass number and the atomic number are even then  $I = 0$ .
2. If the mass number is odd and the atomic number is odd or even then  $I$  has half-integral values.
3. If the mass number is even and the atomic number is odd then  $I$  has integral values.

Both protons and neutrons have nuclear spin  $= 1/2$ , an unpaired nuclear spin leads to a nuclear moment ( $\mu$ ) given by the equation (2.8).

$$\mu = \gamma I h / 2\pi \quad (2.8)$$

where:  $\gamma$  is the magnetogyric ratio which is constant for each particular nucleus.

In general, there are  $2I + 1$  possible orientations or states of the nucleus derived from the value of the magnetic quantum number ( $m_x$ ).

$$m_x = I, (I-1), \dots, -I$$

The most commonly encountered nucleus in NMR ( $^1\text{H}$  and  $^{13}\text{C}$ ) has a value of  $I$  equal to  $1/2$ , therefore  $m_x = \pm 1/2$ . These states have the same energy in the absence of a magnetic field. When a uniform magnetic field of strength  $H_0$  is applied, the states are no longer degenerate as shown in Figure 4.

Nuclei for which  $I \neq 0$  behave as minute bar magnets whose axes coincide with the axis of spin. When these nuclei are made to spin in a powerful uniform magnetic field ( $H_0$ ), they do not flip over to align their magnetic moments with the field but instead align themselves at some angle  $\theta$  with respect to the direction of the field and their spin axes undergo precessions about the  $H_0$  axis as shown in Figure 5. The precession frequency ( $\omega_0$ ) is measured in radians per second. The angular velocity  $\omega_0$  is independent of  $\theta$  but is proportional to the field strength  $H_0$  as shown in equation (2.9).

$$\omega_0 = \gamma H_0 \quad (2.9)$$

When a second small magnetic field  $H_1$  (represented by the horizontal vector in Figure 6) is applied at an angle perpendicular to  $H_0$  and is rotated about the  $H_0$  axis with an angular frequency  $\nu_0$ , then the spinning nuclei will absorb energy and flip into a higher energy level, if  $\nu_0$  is exactly equal to the precession frequency  $\omega_0$ . The rotating nuclear magnetic moment is in resonance with the




Figure 4. Energy Levels and Transitions for a Nucleus ( $I=1/2$ ) in a Magnetic Field




Figure 5. Proton Precession in a Magnetic Field ( $H_0$ )




Figure 6. Oscillator Generates Rotating Component of Magnetic Field  $H_1$

field  $H_1$  and the absorption energy ( $\Delta E$ ) is given by

$$\Delta E = h\nu_0 = \mu H_0 / I \quad (2.10)$$

The selection rules for NMR which establish the resonance conditions are:

1.  $\Delta m_x = \pm 1$
2.  $\Delta E = \gamma h H_0 / 2\pi$
3.  $\nu_0 = \gamma H_0 / 2\pi = \omega_0 / 2\pi$

Let us consider a nucleus of spin  $1/2$  which has two possible orientations of nuclear spin  $\alpha$  and  $\beta$ . The coefficients for the absorption and emission of energy between these states are equal for NMR and therefore there would be no net energy transfer from the radiation source to the sample if the population of the two states were equal. However, the sample is in thermal equilibrium and the Boltzmann distribution of energies is maintained.

$$N_\beta / N_\alpha = \exp(-h\nu_0 / kT) \quad (2.11)$$

If  $h\nu_0 \ll kT$ , then

$$N_\beta / N_\alpha = 1 - h\nu_0 / kT \quad (2.12)$$

Whenever  $N_\beta < N_\alpha$  a net absorption occurs. In the rf spectral region (5-400 MHz), values of  $h\nu_0$  are on the order of  $1E-2$  cal. Therefore the excess population of spins in the lower state is  $1E-5$  for ambient temperatures. For this reason NMR is less sensitive than either IR or UV spectroscopy. However the coefficient of absorption is a constant for any given nucleus. Therefore the observed

NMR signal intensity is directly proportional to the number of nuclei producing it. Obviously, the populations of the spin states are perturbed on the absorption of radiation. The system relaxes to the original Boltzmann distribution by two mechanisms. These are spin-spin relaxation and spin-lattice relaxation.

In the spin-spin relaxation mechanism, a nucleus in its upper spin state transfers its energy to a neighboring nucleus of the same isotope by a mutual exchange of spins. This relaxation process does not change the relative spin state population and therefore does nothing to maintain the absorption condition, but it does shorten the lifetime of a given nucleus in the higher state and therefore produces line broadening.

In spin-lattice relaxation, energy from the nuclear spin at higher energy is converted into thermal energy. The high-energy spin state must be properly oriented with respect to the lattice molecules, such that energy can be transferred to the lattice molecules, giving them extra translational or rotational energy. This process is directly responsible for maintaining the unequal distribution of spin states. As a result there is always an excess of nuclei in the lower energy state, and a continuous absorption of energy from the rf source by the sample can occur.



Information Available from  
NMR Spectroscopy

Each nucleus in a molecule has a different electronic environment. The magnetic field associated with the electronic movement around a nucleus interacts with the applied magnetic field  $H_0$ , which successfully shields the applied magnetic field  $H^0$ . The extent to which  $H_0$  is reduced depends upon the shielding constant ( $\sigma$ ) which is different for different nuclei. The effective magnetic field at the nucleus ( $H_{\text{eff}}$ ) is given by equation (2.13).

$$H_{\text{eff}} = H (1 - \sigma) \quad (2.13)$$

$\sigma$  is a characteristic property for each nucleus in a given electronic environment and has a nondimensional value from zero to one.

For the above reason, different resonance conditions prevail for different nuclei. Each resonance can be excited in turn either by varying the frequency  $\nu_0$  at constant  $H_0$  or varying  $H_0$  at constant  $\nu_0$ . Data are usually obtained in the form of a spectrum in which resonance intensities are plotted as a function of  $\nu_0$  at constant  $H_0$ .

Since  $\nu_0$  and  $H_0$  are mutually dependent parameters,  $\nu_0$  is not a general characteristic quantity since it will vary with  $H_0$  for any given nucleus. To overcome this inconvenient problem, a chemical shift ( $\delta$ ) term is introduced. This term is defined to be independent of  $H_0$  and  $\nu_0$  and depends only upon the molecular environment. Nuclei

which have identical environments will have the same chemical shift. The chemical shift is defined as

$$\delta = \Delta\nu * 1E6 / \text{oscillator frequency (ppm)}$$

where,  $\Delta\nu$  is the difference between the radio-frequency at which energy is absorbed by a reference nucleus and the resonance frequency for the nucleus of interest. By using a difference measurement made under identical experimental conditions it is not necessary to determine  $H_{\text{off}}$ , a term which is difficult to obtain with any accuracy. In practice the relative resonance frequencies between the sample and the reference are measured, and  $\delta$  for each nucleus is obtained from its displacement on the frequency scale from the standard which is arbitrarily set at zero ppm. Since the oscillator frequency is in units of MHz but  $\Delta\nu$  is in units of Hz, the fraction is multiplied by  $1E6$  to give numbers for  $\delta$  which are more convenient to use.

Reference compounds must be unreactive and have resonance frequencies in a range which do not overlap or coincide with those for the compounds of interest. Common reference materials are tetramethylsilane (TMS) for non-aqueous solvents and sodium 2,2-dimethyl-2-silapentane-5-sulphonate (DSS) for aqueous solutions.

By convention, NMR spectra are reported in such a fashion that the field increases from left to right. In a comparison of nuclei, the more shielded nuclei are upfield (on the right of the spectra) and the less shielded nuclei are downfield (on the left of spectra).

Chemical shift data are used to distinguish among the types of nuclei. The area under each peak at a given chemical shift is used to determine the number of identical nuclei which contribute to that peak. Determination requires a peak integration step which is a standard electronic capability of all spectrometers. Integrations under  $^1\text{H}$  peaks can be done with  $\pm 2\%$  accuracy. For  $^{13}\text{C}$  however, integrations are less reliable because of the low natural abundance of the isotope.

In high resolution NMR such as Fourier transform (FT) NMR, spin-spin coupling data are the other useful experimental parameters. Spin-spin coupling is a consequence of interactions between nuclear spins of nuclei which are either adjacent or separated by several intervening bonds. A narrow splitting of energy levels allows for several transitions to occur, which produce groups of lines around the principal resonance band referred to as multiplets. The magnitude of the interaction between two nuclei in a particular environment is expressed in terms of a spin-spin coupling constant in Hz ( $J_{\text{AB}}$ ) which is independent of the applied magnetic field. However the analysis of spin-spin multiplets is applicable only under first-order analysis conditions. These conditions are:

1. The magnitude of the chemical shift difference between the nuclei or the group of nuclei must be larger than the spin-spin coupling constant between them.
2. The coupling must involve groups of nuclei that

are magnetically equivalent, not just chemically equivalent.

When first-order analysis is applicable, then

1. A nucleus or a group of nuclei coupled to a set of  $n$  nuclei with spin  $I$  will have the principal resonance split to  $2nI+1$  lines.

2. The relative intensities of the  $2nI+1$  lines are given by the coefficients of the terms in the binomial expansion of  $(x+1)^n$ .

3. The  $2nI+1$  lines are separated by equivalent coupling constants.

4. The coupling between a magnetically equivalent set of nuclei does not affect the spectrum.

Since each nucleus in a molecule is capable of coupling with every other nucleus, a spectrum composed of many multiplet sets is generated. Some electronic techniques are introduced in order to simplify the spectrum such as double resonance or spin decoupling and off-resonance decoupling. The two techniques involve the application of a strong magnetic field oscillating at a frequency close to the resonance frequency of a given nucleus which isolates it from interactions with other nuclei in the molecule. If sufficient radiative power is applied, the "isolated" nucleus will flip reversibly between  $\alpha$  and  $\beta$  states so rapidly that the states are indistinguishable to other nuclei and only an average orientation is observed. Coupling between this and other nuclei disappears and the signal collapses to give a

singlet line as seen in a  $^{13}\text{C}$  NMR fully proton decoupled spectrum. In the other case, that of off-resonance decoupling, a lower decoupling power is selected such that only the long range  $^{13}\text{C}$ - $^1\text{H}$  coupling are diminished. Then only the first neighbor  $^{13}\text{C}$ - $^1\text{H}$  couplings are observed. From the resultant simplified spectrum, assignments for each carbon can be made. Primary, secondary, tertiary and quaternary carbons in the structure are easily distinguishable as quartets, triplets, doublets and singlets, respectively.

### X-ray Diffraction

X-rays are rays of electromagnetic radiation of wavelength between .1-100  $\text{\AA}$ . The most useful region for analytical purpose is 0.7 to 2.0  $\text{\AA}$ . A single crystal X-ray diffraction method is based on the scattering of the X-rays by a single crystal of the compound being investigated. The scattering produces a diffraction pattern only when certain geometrical conditions as expressed by Bragg's law are satisfied. Diffraction patterns may be used to identify molecules and molecular structures.

To generate the X-rays, electrons from a hot cathode accelerated by a high voltage field (20-50 KV) attach to an anode target which is a plate made from a single element, of moderate atomic weight, in an evacuated chamber. When the beam of electrons impinges on the target, the electrons in general are decelerated by multiple interactions with

the electrons of the target. The energy lost is converted into a continuous or white radiation.

As the voltage is increased to a point where the energy is sufficient to knock an innermost (K shell) electron out of the target atom, then an electron from an outer shell falls back into the vacancy. A photon of X-ray is emitted with a wavelength determined by the energy difference between the levels, and hence is characteristic of the target element. When the vacancy is filled by an electron from the L shell or M shell, then the emitted X-radiations are designated as  $K_{\alpha}$  and  $K_{\beta}$  radiations, respectively as shown in Figure 7.

For the single crystal X-ray diffraction, a strong monochromatic beam is needed. This can be obtained by:

1. Choosing characteristic emission lines which are much stronger than the others, and isolating these with the aid of a filter.

2. Using a single crystal aligned at a particular angle and in a specific plane which acts as a monochromatic source. According to Bragg's law (36), the crystal will produce a specific wavelength of reflection.

When the beam of X-rays is passed into the single crystal, which is represented by layers of particles as shown in Figure 8, the X-rays are scattered by interaction with electrons in the atoms or ions of the crystal. The scattered waves would interfere constructively when the conditions for reflection of the beam satisfy the Bragg's




Figure 7. X-ray Spectra with Characteristic Peaks from a Molybdenum Anode  
X-ray Tube Operated at 35 kV




Figure 8. The Bragg Scattering Condition



law according to the following equation,

$$n\lambda = 2d \sin \theta$$

when  $n$  is an order of reflection. In single crystal X-ray diffraction,  $n$  is equal to 1.  $\lambda$  is the wavelength of the emitted x-radiation which strikes between parallel planes separated by a distance  $d$ , and at an angle  $\theta$  equal to the angle of reflection.

The intensities of the beams depend upon the type and the location of the atoms in the fundamental repetitive units of the crystal.

In practice the crystal is rotated around preselected axes such that various sets of planes assume the orientations for Bragg reflections. These planes are defined by three integers  $h$ ,  $k$  and  $l$  referred to as the Miller indices. In inverse space  $h$ ,  $k$  and  $l$  are defined by  $a/x'$ ,  $b/y'$  and  $c/z'$ , where  $a$ ,  $b$  and  $c$  are the unit cell dimensions and  $x'$ ,  $y'$  and  $z'$  are the fractions of  $a$ ,  $b$  and  $c$  cut by the plane, respectively. The intensities of the reflections from each plane are detected by a counter and electronically recorded by a dedicated computer.

From the relationship between the reflection intensities at each Miller index plane and Bragg law, the cell dimension ( $a$ ,  $b$ ,  $c$ ) and the angles of the unit cell ( $\alpha$ ,  $\beta$ ,  $\gamma$ ) are computed. Using this information, a crystal system can be identified. Each crystal system may belong to several different space groups, with different symmetries, depending upon the number of molecules in one unit cell

(Z). Z can be calculated knowing the molecular weight and the density (D) of the molecule from the relationship

$$D = 1.6602 * M.W. * Z/V$$

v is a volume of the unit cell. D is usually in the range 1.1-1.5 g/cc for most organic compounds and 1.6-2.1 g/cc for most organometallic compounds.

Space groups and the equivalent general position of any particular space group generated from symmetry elements which correspond to that space group are listed in the "International Table for X-ray Crystallography" (37).

Intensity data are reduced to obtain structural factors from which an electron density map can be calculated thereby locating the atomic positions in the molecule. Bond lengths and bond angles between atoms in the molecule are subsequently determined.

It is clear from the preceding remarks that X-ray crystallography is an unambiguous technique which provides a complete three dimensional picture of a molecule in the solid state. Structural information is frequently extrapolated from the solid to the solution state to complement information from other methods which are applicable to solutions, e.g. NMR and CD. These three methods have been used in this study of morphine derivatives and products from reactions of morphine with oxidizing agents.

## CHAPTER III

### EXPERIMENTAL

#### Instruments

Circular dichroism measurements were made on a JASCO (Japan Spectroscopic Company) model J-500A automatic recording spectropolarimeter. Data collection and manipulation was done either manually or with the data processor Model DP-500N accessory.

The instrument has a wavelength range from 180 to 800 nanometer (nm) and a sensitivity scale from .1 to 50 milli-degree per centimeter ( $m^{\circ}/cm$ ).

The light source is a 450 watt high pressure Xe arc lamp cooled by the flow of water at approximately 2L/min. The instrument is operated under 23 amp lamp current, and 105 volt line voltage. It is continuously purged with nitrogen gas boiled from a liquid nitrogen dewar at a flow rate of about 2L/min in order to prevent the production and accumulation of ozone from oxygen by UV radiation which would cause deterioration of the optics system.

The ellipticity scale is calibrated using a 0.05% (W/V) androsterone/dioxane solution in a 10 mm cell. The wavelength is scanned from 350 to 250 nm and the CD peak

height at 304 nm maximum is set equal to 96.2 mm of chart scale at sensitivity  $20 \text{ m}^0/\text{cm}$ . Baseline corrections were made for the solvent blank by spectral subtraction on the data processor. Sensitivity, scan rate, time constant and repeat time functions were selected which optimized the S/N ratio. The cells used for the analyses were made from quartz crystal. Path lengths of 1 or 10 millimeters were used, determined by the concentration of the solution under study.

UV-visible spectral measurements were made on a Perkin-Elmer UV-visible spectrophotometer Model 552, which is a double beam instrument. The light sources are tungsten-bromide and deuterium lamps with an automatic change made at 315 nm. This instrument has a wavelength range from 190 to 900 nm with  $\pm 0.5 \text{ nm}$  accuracy. The instrument is equipped with an automatic baseline adjustment and is connected to the Perkin-Elmer recorder Model 561 which has an accuracy less than  $\pm 0.4\%$  of span. The accessible absorbance range is from 0.000 to 3.000 absorbance unit (A). Measurements were made using a scan speed of 120 nm/min, a recording speed of 20 nm/cm and slit width 1 nm.

NMR spectra were run on a Varian XL-300 multinuclear NMR spectrometer with the fixed magnetic field strength 7.05 Telsa (1 Telsa =  $1\text{E}4$  Gauss) and the resonance frequencies of  $^1\text{H}$  and  $^{13}\text{C}$  nuclei at 299.944 and 75.45 MHz, respectively. Sample tubes of 5 mm o.d. were used for all

measurements. Chemical shifts were recorded from 0 to 10 ppm for  $^1\text{H}$  and from 0 to 200 ppm for  $^{13}\text{C}$  spectra, respectively.

X-ray structural determinations of the morphine derivatives and its oxidation product were done on a four-circle automated diffractometer Model Syntex P3. The instrument uses molybdenum radiation ( $\lambda = 0.71069 \text{ \AA}$ ) with a graphite monochromator fixed at an angle of  $6.1^\circ$  with respect to the beam as it emerges from the X-ray tube and impinges on the 110 plane.

Mass spectra were made on a CEC 21-110B high resolution spectrometer. It is equipped with a Data General NOVA 3112 DS-50S data system.

Weighings were made on two balances. A Sartorius balance (Model 2403) was used for sample weights in excess of 10 mg. A Cahn electrobalance (Model 2000RG) was used for the sample sizes in the 0.001 mg to 10.000 mg range.

### Chemicals

Morphine sulfate pentahydrate and morphine free base were the main starting reagents used in this study. Both were pure standard reagents, which were obtained from Mallinckrodt Inc., and from the National Institute for Drug Abuse (NIDA) via the Research Triangle Laboratories (RTI), respectively. Apomorphine hydrochloride was obtained from Applied Science Division and used without further purification. Distilled water and deionized water

were used as solvents for neutral and basic reactions while concentrated sulfuric acid, 98.0%, obtained from Fisher Scientific Company was used as the solvent in the acid reactions. Other reagents used in acid, neutral and basic reactions of morphine and its derivatives are listed in Tables I and II.

TABLE I  
REAGENTS USED IN ACID OXIDATION REACTION OF MORPHINE

Reagent	Supplier
Ammonium molybdate $(\text{NH}_4)_6\text{Mo}_7\text{O}_{24} \cdot 4\text{H}_2\text{O}$	Mallinckrodt
Ammonium vanadate $\text{NH}_4\text{VO}_3$	Fisher Scientific Company
Selenious acid (98%) $\text{H}_2\text{SeO}_3$	Aldrich Chemical Company, Inc.
Chloroform (spectra-analyzed grade) $\text{CHCl}_3$	Fisher Scientific Company

Methyl iodide and methanol were obtained from Fisher Scientific Company. Phenol, o,o'-biphenol (99%), anisole, 3,3'-bianisole (97%) and the NMR deuterated solvent, DMSO(d6) were obtained from Aldrich Chemical Company, Inc.

TLC-plates used in the separation of the mixtures were obtained from Analtech. All reagents are reagent grade unless specified and they were used without any further purification.

TABLE II  
REAGENTS USED IN NEUTRAL AND BASIC REACTIONS OF MORPHINE

Reagent	Supplier
Cadmium chloride	J. T. Baker Chemical Company
Copper sulfate pentahydrate	Fisher Scientific Company
Lead Chloride	J. T. Baker Chemical Company
Nickelous chloride	J. T. Baker Chemical Company
Potassium hydroxide	Mallinckrodt Inc.
Potassium ferricyanide	Fisher Scientific Company
Potassium permanganate	J. T. Baker Chemical Company
Silver nitrate	Fisher Scientific Company
Zinc chloride	Mallinckrodt Inc.

## Experimental Procedures

### Oxidation Reaction of Morphine in Basic Solution

Pseudomorphine was prepared by following the procedure reported by Bentley (11). The white water-insoluble product was recrystallized from a 1:1 mixture of concentrated ammonium hydroxide and distilled water. Fine crystals formed after two days. The recrystallized product was used for CD, UV and NMR studies. An attempt was made to grow single crystals of pseudomorphine and its derivatives such as dihydrochloride, dihydrobromide, monosulfate and monotartrate salts (38) for an X-ray structural study of pseudomorphine.

### Oxidation Reaction of Morphine by Metal Ions in Aqueous Solution

The metal salts used in this study were silver nitrate ( $\text{AgNO}_3$ ), cadmium chloride ( $\text{CdCl}_2$ ), copper sulfate pentahydrate ( $\text{CuSO}_4 \cdot 5\text{H}_2\text{O}$ ), lead chloride ( $\text{PbCl}_2$ ), nickelous chloride ( $\text{NiCl}_2$ ), zinc chloride ( $\text{ZnCl}_2$ ) and potassium ferricyanide ( $\text{K}_3\text{Fe}(\text{CN})_6$ ). The salts were dissolved in distilled water in which the mole ratio of the metal ion to morphine from base was equal to 10:1. For  $\text{Ag}^{+1}$  a systematic study was made in which the ratio was varied 1:1 to 10:1 in deionized water. The morphine free base concentrations used in these reactions were in the range



of  $1.5\text{E-}4$  to  $3.5\text{E-}5$  M.

Oxidation of Pseudomorphine. Potassium permanganate ( $\text{KMnO}_4$ ) in .2% KOH and the various metal ions as mentioned above were also used for experiments in which pseudomorphine was oxidized.

Oxidation of Morphine and Apomorphine  
in Strong Acid Solution

The acid used in this study was concentrated  $\text{H}_2\text{SO}_4$  which is a common solvent used in most of the color test reagents for screening opium alkaloids as shown in Table III. The color test reagents used were the Froehde, Mecke, and Mandelin reagents.

1. 200 mg portions of morphine sulfate were added to 50 ml concentrated  $\text{H}_2\text{SO}_4$  and to similar solutions containing the three color test reagents. CD and UV-visible spectra of the solutions were recorded after appropriate dilution with the solvent. The solutions were allowed to stand at room temperature for a period of time (2 days for color test reagents and 16 days for only concentrated sulfuric acid) before being diluted with distilled water and extracted with chloroform. Chloroform extracts were used for CD, UV-visible, mass spectrometry and thin-layer chromatography studies. The structure of one of the products (Compound I) was determined by X-ray crystallography.

TABLE III  
LIST OF COLOR TEST REAGENTS

Reagent	Composition	Color developed after a week
Froehde	700 mg of ammonium molybdate in 50 ml of concentrated sulfuric acid	dark purple to green
Mecke	250 mg of seleneous acid in 25 ml of concentrated sulfuric acid	dark brown to greenish gray
Mandelin	1 gm of ammonium vanadate in 100 ml concentrated sulfuric acid	dark brown purple

2. Apomorphine hydrochloride (20 mg) was added to 10 ml concentrated sulfuric acid and to 10 ml of Froehde reagent, respectively. The UV-visible and CD spectra of the solutions were recorded as a function of time. After equivalent times had lapsed, the solutions were diluted and extracted in the same manner as was described morphine sulfate in (1) and CD and UV-visible spectra of the solutions were recorded.

3. Thin-layer chromatography (TLC) was used as a technique to determine the number of components in the chloroform extracts. A mixture of chloroform, acetone,

methanol and dimethylamine in the ratio of 5:4:4:1 by volume was used as the eluant. A visual separation was confirmed and developed by UV light. The experiment was performed at  $5 \pm 2^{\circ}\text{C}$ .

### X-ray Structural Analysis

For X-ray structural determinations, a single crystal of morphine sulfate was obtained by the recrystallization of morphine sulfate from water. A methyliodide salt was also prepared by the reaction of excess methyl iodide with morphine free base, by refluxing in methanol for two hours. After cooling the solution in a refrigerator for 24 hours long needle-shaped crystal of morphine methyliodide formed. The crystals were recrystallized from water.

Single crystals of morphine sulfate, morphine methyliodide and Compound I were sealed in capillaries and mounted on the Syntex P3 automated diffractometer. Unit cell dimensions of the three crystals were determined by least squares refinement of the best angular position for fifteen independent reflections during normal alignment procedures. The data were collected at room temperature using a variable scan rate, a  $\theta$ - $2\theta$  scan mode and a scan width of  $1.2 \text{ \AA}^{\circ}$  below  $K_{\alpha 1}$  and above  $K_{\alpha 2}$  to a maximum  $2\theta$  value at  $116^{\circ}$ . Backgrounds were measured at each side of the scan for a combined time equal to the total scan time. The intensities of three standard reflections were measured after every 97 reflections and the intensities of these

reflections showed less than 8% variation. Corrections for decomposition were considered to be unnecessary. Data were corrected for Lorentz and polarization effects.

The structures of morphine sulfate and morphine methyliodide were solved from the Patterson synthesis (39) to locate the heavy atoms. Successive least squares, difference Fourier cycles allowed for the location of the remainder of the non-hydrogen atoms. The unknown structure of Compound I was solved by direct methods using MULTAN80 (40).

Refinement of the scale factor, and positional (x, y, z) and anisotropic thermal parameters ( $U_{11}$ ,  $U_{12}$ ,  $U_{13}$ ,  $U_{22}$ ,  $U_{33}$ ,  $U_{23}$ ) for all non-hydrogen atoms (41) was carried out to convergence. All hydrogen positional parameters of morphine methyliodide were determined from a difference Fourier synthesis. Positions of hydrogens on C6, C8 and C17 of morphine methyliodide, however, were calculated assuming normal geometry and a C-H distance of 0.97 Å. These hydrogen positional parameters and the associated isotropic thermal parameter (U) of Compound I were refined along with non-hydrogen parameters in the final cycles of refinement. The hydrogen positional parameters of morphine sulfate were refined along with non-hydrogen parameters in the final cycles of refinement but the associated isotropic thermal parameters were assigned to be 0.03 and were held invariant. In the case of the methyliodide salt, all parameters associated with hydrogen

atoms were held invariant with the isotropic thermal parameter of  $U = 0.03$ .

The final cycle of refinement (function minimized  $\Sigma(|F_{\text{obs}}| - |F_{\text{cal}}|)^2$ ) led to the ultimate agreement factor  $R$  defined as

$$R = (\Sigma ||F_{\text{obs}}| - |F_{\text{cal}}||) * 100 / \Sigma |F_{\text{obs}}|$$

Anomalous dispersion corrections were made for the iodide ion in the morphine methyliodide structure. Scattering factors for the elements were taken from Cromer and Mann (42). Unit weights were applied throughout the experiments.

#### Solid State Circular Dichroism Study of Morphine Methyliodide

The morphine methyliodide used in the previous section was used for a solid state circular dichroism (SSCD) study. The result was compared to the SSCD spectra for morphine free base and morphine sulfate reported earlier by Purdie and Bowen (31).

Samples for SSCD were prepared by pressing the morphine salt in KBr pellets formed in the usual way. Infra-red grade KBr (Aldrich Chemical Co.) and the compound were dried at  $110^{\circ}\text{C}$  at least 48 hours prior to being pressed. A typical specimen consisted of approximately 0.8 mg of the salt in 80 mg KBr. The KBr pellet must be thin (ca. 0.2 mm) for adequate UV transmission.

The spectrum was taken on the instrument in the

normal manner with the sample supported in a holder constructed to fit both the pellet press and the carriage in the sample compartment of the CD instrument in order to avoid any damage to the fragile specimen. Tests for the existence of linear dichroism components to the signal were made by rotating the sample by various intervals from  $0^{\circ}$  to  $180^{\circ}$  and repeating the spectrum. The sample was dissolved in distilled water and an isotropic solution CD spectrum was obtained. Experiments were repeated for the SSCD of morphine free base and morphine sulfate to obtain spectra under the same conditions on the same instrument for a legitimate comparison of the results.

## CHAPTER IV

### EXPERIMENTAL RESULTS

#### Oxidation of Morphine in Basic Solution

Pseudomorphine is believed to be the first product formed from the oxidation of morphine by dissolved oxygen in basic solution as mentioned earlier. Although a dimer has been confirmed, its structure has not been established with any certainty, hence the reason for an X-ray structure analysis. Every attempt to produce suitable single crystals of either pseudomorphine or its derivatives failed, although a variety of solvents and experimental conditions were explored. Therefore NMR spectroscopy was used in order to obtain some structural information.

#### NMR Study of Pseudomorphine

A comparison of the  $^1\text{H}$  NMR spectrum of pseudomorphine (Figure 9) with that for morphine sulfate (Figure 10) in DMSO( $d_6$ ) shows a significant change in the region of the aromatic protons ( $\delta = 6.2\text{-}6.6$  ppm). The quartet observed for morphine sulfate reduces to a singlet for pseudomorphine. Integrations of these two spectra show that the relative numbers of hydrogen(s) under the quartet and the




Figure 9.  $^1\text{H}$  NMR Spectrum of Pseudomorphine in  $\text{DMSO}-d_6$






Figure 10.  $^1\text{H}$  NMR Spectrum of Morphine Sulfate in DMSO ( $\text{d}_6$ )

singlet are two and one for morphine sulfate and pseudomorphine, respectively. The remaining peaks in both spectra are quite similar.

The fully decoupled and off-resonance decoupled  $^{13}\text{C}$  NMR spectra of pseudomorphine (Figures 11 and 12) on the other hand compare very well with those for morphine sulfate (Figures 13 and 14). The spectra show that both molecules are composed of seventeen different carbons, eight of which are unsaturated as expected. Chemical shifts for the unsaturated-carbon range from 119 to 150 ppm from the internal standard TMS. Assignments for the carbon chemical shifts for pseudomorphine (Table IV) were based on the equivalent assignments made for morphine sulfate. The latter assignments were made from data from the off-resonance experiments for morphine sulfate and codeine (Figure 15), as well as the fully decoupled spectrum for codeine (Figure 16). The carbon assignments for morphine sulfate were in accord with an earlier interpretation reported by Carroll, Moreland, Brine, and Kepler (43). Most of the carbon chemical shifts of pseudomorphine correspond with the analogous carbon chemical shifts for morphine sulfate to within  $\pm 1$  ppm, except for carbons C1, C2, C3, C11 and C12, which are the carbons directly affected by the dimerization of the two morphine molecules in producing the pseudomorphine structure. A reversal in the order of the C1 and C2 peaks for pseudomorphine compared to morphine




Figure 11.  $^{13}\text{C}$  NMR Spectra of Pseudomorphine in  $\text{DMSO-d}_6$ : a) Fully Decoupled, and b) Off-resonance




Figure 12. The Unsaturated Region Expansion of the Pseudomorphine  $^{13}\text{C}$  Off-resonance NMR Spectrum




Figure 13.  $^{13}\text{C}$  NMR (Fully Decoupled) Spectrum of Morphine Sulfate in DMSO ( $\text{d}_6$ )




Figure 14.  $^{13}\text{C}$  NMR (Off-resonance) Spectrum of Morphine Sulfate in DMSO ( $\text{d}_6$ )

TABLE IV  
 $^{13}\text{C}$  NMR CHEMICAL SHIFTS OF MORPHINE SULFATE AND  
PSEUDOMORPHINE IN DMSO ( $d_6$ )

Identification Carbon	Chemical Shifts (ppm) from TMS	
	Morphine Sulfate	Pseudomorphine
1	119.05	121.0
2	117.10	124.3
3	139.10	136.9
4	146.40	147.5
5	90.90	91.7
6	66.20	66.4
7	134.40	133.5
8	128.70	128.5
9	58.10	58.2
10	20.20	20.3
11	125.50	128.3
12	131.04	129.2
13	42.97	42.7
14	40.63	40.3
15	35.60	35.2
16	46.05	46.0
17	42.83	42.6




Figure 15.  $^{13}\text{C}$  NMR (Off-resonance) Spectrum of Codeine in DMSO ( $\text{d}_6$ )






Figure 16.  $^{13}\text{C}$  NMR (Fully Decoupled) Spectrum of Codeine in DMSO ( $\text{d}_6$ )

sulfate is suggested and supported by the comparison of the results with chemical shift data among monomers and dimers of two pairs of model compounds, phenol (Figures 17, 18), o,o'-biphenol (Figures 19, 20) and anisol (Figures 21, 22) and 3,3'-bianisole (Figures 23, 24). Chemical shift data for the model compounds are listed in Table V and Table VI. The chemical shifts assignments for phenol and anisole were based on the values reported in the Atlas of Carbon-13 NMR Data (43).

#### Oxidation of Pseudomorphine in Base

Pseudomorphine is not oxidized by air on standing in basic solutions (concentrated ammonium hydroxide or .2% KOH). However, pseudomorphine as the potassium salt is quickly oxidized to an unidentified higher oxidation product (HOP) on the addition of an excess of 1% solution of  $\text{KMnO}_4$  which is decolorized in the reaction. The dark brown precipitate of  $\text{MnO}_2$  is formed and the solution turns light brown in color. Evidence for a new compound is apparent from the changes which occur in the CD spectrum of pseudomorphine (Figure 25). The positive Cotton band whose maximum occurs at 320 nm for pseudomorphine is absent for HOP and the 275 nm negative maximum is shifted to the shorter wavelength at 263 nm. The CD spectral characteristics for morphine, pseudomorphine and HOP are presented in Table VII for easy comparison.




Figure 17.  $^{13}\text{C}$  NMR (Fully Decoupled) Spectrum of Phenol in DMSO ( $\text{d}_6$ )




Figure 18.  $^{13}\text{C}$  NMR (Off-resonance) Spectrum of Phenol in DMSO ( $\text{d}_6$ )




Figure 19.  $^{13}\text{C}$  NMR (Fully Decoupled) Spectrum of *o,o'*-Biphenol in  $\text{DMSO-}d_6$




Figure 20.  $^{13}\text{C}$  NMR (Off-resonance) Spectrum of *o,o'*-Biphenol in  $\text{DMSO-d}_6$




Figure 21.  $^{13}\text{C}$  NMR (Fully Decoupled) Spectrum of Anisole in DMSO ( $\text{d}_6$ )




Figure 22.  $^{13}\text{C}$  NMR (Off-resonance) Spectrum of Anisole in  $\text{DMSO-d}_6$






Figure 23.  $^{13}\text{C}$  NMR (Fully Decoupled) Spectrum of 3,3'-Bianisole in DMSO ( $\text{d}_6$ )




Figure 24.  $^{13}\text{C}$  NMR (Off-resonance) Spectrum of 3,3'-Bianisole in DMSO ( $\text{d}_6$ )

TABLE V  
 $^{13}\text{C}$  NMR CHEMICAL SHIFTS OF PHENOL AND  
o,o'-BIPHENOL IN DMSO (d6)

Identification Carbon	Chemical Shifts (ppm) from TMS	
	Phenol	o,o'-Biphenol
1	157.516	154.530
2	115.750	126.026
3	129.562	131.624*
4	119.140	118.986
5	129.562	128.179*
6	115.750	115.843

\*May be interchanged.

TABLE VI  
 $^{13}\text{C}$  NMR CHEMICAL SHIFTS OF ANISOLE AND  
3,3'-BIANISOLE IN DMSO (d6)

Identification Carbon	Chemical Shifts (ppm) from TMS	
	Anisole	3,3'-Bianisole
1	155.487	155.608
2	109.731	108.377
3	125.257	137.606
4	116.303	115.049
5	125.257	125.782
6	109.731	108.861
7	50.406	50.937




Figure 25. CD Spectra of a) Morphine in Distilled Water, b) Pseudomorphine in Distilled Water, and c) Oxidation Product of Pseudomorphine in Distilled Water


TABLE VII  
CD DATA OF MORPHINE, PSEUDOMORPHINE AND HOP  
IN DISTILLED WATER

Compound	$\lambda_{\text{max}}^{+}$ (nm)	$\lambda_{\text{max}}^{-}$ (nm)	$\lambda^0$ (nm)
Morphine	303, 247	285	296, 270, 236.5
Pseudomorphine	320, 243	274-5	303, 259, 229
HOP	—	263	247


## Oxidation of Morphine and Pseudomorphine by Metal Ions in Neutral Solution

A number of simple ions and two complex metal ions were used as potential oxidizing agents for morphine and pseudomorphine in neutral solutions. Evidence for new products was obtained from a CD study since both alkaloids and the possible products are CD active (Figure 26).


The oxidation of morphine to pseudomorphine is rapid when the mole-ratio of silver ion to morphine is 3:1 or higher. Within seven days the CD spectra (Figure 27) of the solutions correspond with the CD spectrum for HOP described for the basic  $\text{KMnO}_4$  reaction. Reduced silver is deposited on the wall of the containers. At lower mole-ratios (1:1 and 1:2), the oxidation appeared to terminate at pseudomorphine after 24 hours (Figure 28) and no evidence for the formation of HOP was found even though the solutions were allowed to stand at room temperature for one week. In excess  $\text{KMnO}_4$  and  $\text{K}_3\text{Fe}(\text{CN})_6$ , morphine was immediately oxidized to pseudomorphine and subsequently more slowly to HOP with time (Figure 29). For the metal ions, copper(II), lead(II), nickel(II), cadmium(II) and zinc(II), there were no significant changes in the CD spectrum of morphine even after the solutions were kept at room temperature for a week. When pseudomorphine was used as a reducing agent, instead of morphine, it was found that the pseudomorphine was oxidized to the HOP in



Morphine



Apomorphine



Pseudomorphine

Figure 26. Structures of Morphine, Apomorphine and Pseudomorphine. (\* Indicates chiral carbon.)






Figure 27. CD Spectra of Solutions from  $\text{Ag}^+$ -morphine at Mole-ratio 3:1 at Different Times.  
a) After mixing; b) After 15 minutes;  
c) After 24 hours; d) After a week.




Figure 28. CD Spectra of Solutions from Ag<sup>+</sup>-Morphine Reaction at Different Times. a) After Mixing, and b) After 24 Hours




Figure 29. CD Spectra of Solutions from Morphine in Excess of  $K_3Fe(CN)_6$  in Distilled Water. a) After mixing; b) After 1.5 hours; c) After 3 days.

excess of silver(I) ion, copper(II) ion,  $\text{KMnO}_4$  and  $\text{K}_3\text{Fe}(\text{CN})_6$  (Figure 30) but not with lead(II), nickel(II), cadmium(II) and zinc(II) salts.

#### Oxidation of Morphine and Apomorphine in Strong Acid Solutions

##### Morphine Sulfate

When morphine sulfate was added to concentrated sulfuric acid and kept in air at room temperature, the color of the solution changed with time. At the outset the color of the solution of morphine sulfate in concentrated sulfuric acid is brownish-orange. The UV-visible spectra (Figures 31, 32) show four bands with wavelengths maxima at 210, 240, 270 and 450 nm. The relative intensities of the bands around 240 and 270 nm changed with time showing a decrease for the 240 nm band and an increase for the 270 nm band. Simultaneously the band maxima shifted to longer wavelengths. At the same time the 450 nm band was observed to diminish and a shoulder to the band at 400 nm appeared. The data as a function of time are summarized in Table VIII.

CD spectra of an identical solution (Figure 33) were recorded over the same time frame; that is, immediately after mixing and at intervals over the next sixteen days before a chloroform extraction was performed. The CD spectrum of morphine in concentrated sulfuric acid is




Figure 30. CD Spectra of Solutions Containing 1:10 Mole Ratio of Pseudomorphine to  
a)  $\text{Cu}^{+2}$  2 Days After Mixing; b)  $\text{Ag}^{+}$  2 Days After Mixing; c)  $\text{KMnO}_4$  After Mixing; d)  $\text{K}_3\text{Fe}(\text{CN})_6$  After Mixing.




Figure 31. UV Spectra of Morphine Sulfate in Concentrated Sulfuric Acid.  
a) After Mixing, b) After Standing at Room Temperature for 2 Days;  
c) After Standing at Room Temperature for 16 Days.




Figure 32. Visible Spectra of Morphine Sulfate in Concentrated Sulfuric Acid.  
a) After Mixing, b) After Standing at Room Temperature for 2 Days,  
c) After Standing at Room Temperature for 16 Days.

TABLE VIII  
UV-VISIBLE DATA OF MORPHINE SULFATE IN  
CONCENTRATED SULFURIC ACID

Time Interval	Color Observed	$\lambda_{\text{max}}$ (nm)
0	brownish-orange	210, 240, 270, 450
2 days	green	210, 242, 278, 400 <sup>sh</sup> 635
16 days	brownish-gray	210, 245, 285, 400 <sup>sh</sup> 520
1 month	brownish-purple	—






Figure 33. CD Spectra of Morphine Sulfate in Concentrated Sulfuric Acid. a) After Mixing; b) After Standing at Room Temperature for 16 Days

inverted from the spectrum in aqueous solution. Maxima occur at 270 nm ( $\lambda_{\text{max}}^+$ ) and at 236 nm ( $\lambda_{\text{max}}^+$ ). Dilution with water did not restore the spectrum to that in aqueous acid which confirms that an irreversible reaction has occurred on solution (45). The  $\lambda_{\text{max}}^+$  of the acid solution was observed to shift from 270 nm to 282 nm and  $\lambda^0$  from 255 nm to 257 nm with time.

When the color test reagents of Froehde, Mandelin and Mecke which contain the Mn(VI), V(V) and Se(IV) oxidation states in complex oxyanions dissolved in concentrated sulfuric acid are used instead of just the acid, the color changed to be either dark blue-green or dark blue within one hour. UV-visible and CD spectra for these solutions were reported previously by Ketkeaw (45). Some contribution to the color of the final product could arise from the reduced forms of the oxyanions.

The chloroform extracts of these solutions in contrast were purple or purple-blue in color. Presumably inorganic salts were not extracted. The aqueous layer after extraction was green. The UV-visible spectra (Figure 34, 35) for all of the chloroform extracts are similar with  $\lambda_{\text{max}}$  at 570, 310 and 238 nm which is indicative of a common product or mixture of products. No CD spectrum was observed for the chloroform extract solutions.

A mass spectrum (Figure 36) of the morphine-Froehde extract shows that the maximum molecular ion has a mass equal to 325.25 g/mole.




Figure 34. UV-visible Spectra of the Chloroform Extracts of the Solutions Containing Morphine Sulfate with a) Froehde Reagent, b) Mandelin Reagent, and c) Mecke Reagent.




Figure 35. UV-visible Spectrum of the Chloroform Extract of the Solution Containing Morphine Sulfate and Concentrated Sulfuric Acid




Figure 36. Mass Spectrum of the Chloroform Extracted Solution from Morphine Sulfate-Froehde Reaction

Apomorphine Hydrochloride

Like morphine sulfate, a solution of apomorphine hydrochloride in concentrated sulfuric acid undergoes a series of color changes with time (Table IX) with accompanying changes in the UV-visible and CD spectra.

TABLE IX  
COLOR CHANGE DATA OF APOMORPHINE IN  
CONCENTRATED SULFURIC ACID

Time Interval	Color Observed
0	clear
2 days	purple-gray
15 days	brownish-pink

The UV spectrum (Figure 37) of a freshly prepared solution in concentrated sulfuric acid is quite different from that of the solution in distilled water (Figure 38). In the acid solution, the principal bands with maxima around 240 and 280 nm are shifted to longer wavelengths and the relative intensities again changed with time showing a decrease at 240 nm and an increase at 280 nm (Figures 39, 40). Unlike morphine sulphate, the CD




Figure 37. UV Spectrum of Apomorphine Hydrochloride in Concentrated Sulfuric Acid (After Mixing)




Figure 38. UV Spectrum of Apomorphine Hydrochloride in Distilled Water






Figure 39. UV Spectrum of Apomorphine Hydrochloride in Concentrated Sulfuric Acid After 2 Weeks




Figure 40. Visible Spectrum of Apomorphine Hydrochloride in Concentrated Sulfuric Acid After 2 Weeks

spectrum (Figure 41) of apomorphine hydrochloride in concentrated sulfuric acid is not inverted when compared to the spectrum in aqueous solution and has characteristic wavelengths  $\lambda_{\text{max}}^+$ ,  $\lambda_{\text{max}}^-$  and  $\lambda^0$  only slightly removed from those for apomorphine hydrochloride in distilled water (Figure 42). On dilution of the acid solution with water the spectrum changes to become essentially that for the aqueous solution. This indicates the presence of a reversible reaction. The air oxidation of apomorphine hydrochloride in concentrated sulfuric acid occurs gradually with time. The general appearance of the CD spectrum of apomorphine resembles that for morphine. Characteristic wavelengths differ by only a few nanometers:  $\lambda_{\text{max}}^+ = 285 \text{ nm}$ ,  $\lambda^0 = 257 \text{ nm}$  and  $\lambda_{\text{max}}^- = 232 \text{ nm}$ . Using the Froehde reagent the reaction occurred much more quickly. A dark purple color developed within an hour. After the solution was allowed to stand in air at room temperature for two days, it was extracted with chloroform. A purple color developed in the chloroform layer. The UV-visible spectrum (Figure 43) of this solution is identical to that for the corresponding morphine oxidation product extract in chloroform.

The apomorphine hydrochloride reaction with concentrated sulfuric acid alone was extracted in the same manner after it stood at room temperature, in air, for two weeks. The chloroform extract was clear in color. At the highest instrument sensitivity, however, the




Figure 41. CD Spectra of Apomorphine Hydrochloride in Concentrated Sulfuric Acid. a) After Mixing, b) After Standing at Room Temperature for 2 Weeks




Figure 42. CD Spectra of Apomorphine Hydrochloride in Distilled Water




Figure 43. UV-visible Spectra of the Chloroform Extracts of a) Apomorphine Hydrochloride with Concentrated Sulfuric Acid, b) with Froehde Reagent

UV-visible spectrum (Figure 43) obtained is similar to that for the extract of the apomorphine oxidation product after reaction with Froehde reagent. Relative band intensities are not exactly comparable, which may be attributed to different distributions of the products of the reactions. That mixtures are obtained is evident from the TLC study described in the next section.

#### Thin-layer Chromatography

The purple chloroform extracts from the final solutions of the morphine sulfate-color test reagents and apomorphine hydrochloride-Froehde reactions were spotted for a TLC separation. Under the conditions used, the mixture separated into three components with  $R_f$  values equal to .34, .72 and .90, respectively. The slowest advancing component was green in color. The other two components were only observable under UV light. TLC comparisons were run against morphine and apomorphine, under the same conditions. Only one spot was observed for morphine, a colorless, UV sensitive component with an  $R_f$  value equal to 0.39 and also only one green and UV active spot for apomorphine ( $R_f = 0.73$ ).

#### X-ray Structural Analysis

##### Structural Comparison of Morphine

##### Sulfate and Morphine Methyliodide

Recrystallized morphine sulfate and morphine methyliodide were obtained as transparent, colorless single

crystals. Photographic investigations indicated orthorhombic symmetry for both structures. Systematic absences are compatible with space groups  $P2_12_12_1$  and  $P2_12_12_1$ , respectively. The calculated unit cell dimensions are tabulated in Tables X and XI. Data were collected for a total of 4278 points for morphine sulfate and 3001 points for morphine methyliodide. After the removal of redundant and space group-forbidden data, 1754 and 2575 "legitimate" reflections were considered valid ( $I > 3.00\sigma(I)$ ) for morphine sulfate and morphine methyliodide, respectively (Appendix A and B).

The final atomic positional and thermal parameters; together with their standard deviations are presented in Tables XII, XIII, XIV, XV, XVI, XVII and XVIII. The least squares refinement produced the final R factor for morphine sulfate equal to 6.2% and 6.1% for morphine methyliodide. The bond angles and bond distances between atoms in these two molecules are listed in Tables XIX and XX.

The confirmations of ring D and ring E of the two morphine salts were compared with those for morphine hydrate. The torsional angles of ring D and ring E of morphine hydrate were calculated based on the structural information reported by Bye (46) and are listed in Table XXI.

Three dimensional ORTEP projections are presented in Figures 44 and 45, respectively. The packing diagram showing the hydrogen bonding network of morphine sulfate,



TABLE X  
CRYSTAL DATA FOR  $(C_{17}H_{20}NO_3^+)_2SO_4^{=}\cdot 3H_2O$

Formula	$C_{34}H_{38}N_2O_6\cdot H_2SO_4\cdot 3H_2O$
MWT	722.81
<u>a</u>	8.463(3) Å
<u>b</u>	6.825(2)
<u>c</u>	30.71(1)
$\alpha = \beta = \gamma$	90.0°
V	1774(1) Å <sup>3</sup>
F(000)	708
$\lambda_{MoK\alpha}$	0.71069 Å
$\mu_{MoK\alpha}$	1.6477
D <sub>calc</sub>	1.3529 g cm <sup>-3</sup>
Z	2
Space	P2 <sub>1</sub> <sup>2</sup> <sub>1</sub>
Systematic absences	00l, l ≠ 2n 0k0, k ≠ 2n
Obs. reflections	1754
Final R Factor	6.2%

TABLE XI  
CRYSTAL DATA FOR  $(C_{18}H_{22}NO_3)^+I^- \cdot H_2O$

Formula	$C_{18}H_{22}NO_3I \cdot H_2O$
MST	445.298
$\underline{a}$	9.574(1) Å
$\underline{b}$	10.614(2)
$\underline{c}$	17.549(4)
$\alpha = \beta = \gamma$	90.0°
V	1783.3(6) Å <sup>3</sup>
F(000)	856
$\lambda_{MoK_\alpha}$	0.7169 Å
$\mu_{MoK_\alpha}$	77.7342
$D_{calc}$	1.6582 g cm <sup>-3</sup>
Z	4
Space group	P2 <sub>1</sub> <sup>2</sup> <sub>1</sub> <sup>2</sup> <sub>1</sub>
Systematic absences	h00, h ≠ 2n 0k0, k ≠ 2n 00l, l ≠ 2n
Obs. reflections	2575
Final R factor	6.1%

TABLE XII  
 NONHYDROGEN POSITIONAL PARAMETERS FOR  $(C_{17}H_{20}NO_3^+)_2SO_4^{=}$

Atom	x( $\sigma(x)$ )	y( $\sigma(y)$ )	z( $\sigma(z)$ )
S1	0.6527(3)	0.5000	0.5000
O1	0.4294(7)	0.1037(11)	0.2791(2)
O2	0.3741(7)	0.6331(9)	0.3823(2)
O3	0.3315(6)	0.2453(9)	0.3623(2)
O111	0.7539(7)	0.3284(9)	0.4926(2)
O112	0.5537(8)	0.4730(11)	0.5380(2)
O113	0.5251(8)	0.5064(15)	0.3038(2)
O114	0.6911(8)	0.3011(10)	0.0108(2)
O115	0.9093(10)	0.0000(0)	0.0000(0)
N1	-.2470(7)	0.1135(9)	0.4185(2)
C1	-.0034(10)	0.0153(14)	0.2802(3)
C2	0.1561(10)	0.0129(14)	0.2669(3)
C3	0.2760(10)	0.0841(12)	0.2930(3)
C4	0.2333(9)	0.1507(11)	0.3332(3)
C5	0.2319(9)	0.3281(10)	0.3976(3)
C6	0.2207(11)	0.5536(11)	0.3900(3)
C7	0.1052(11)	0.6083(11)	0.3561(3)
C8	-.0372(10)	0.5260(10)	0.3567(3)
C9	-.2213(8)	0.2444(11)	0.3790(3)
C10	-.2102(8)	0.1260(13)	0.3366(3)
C11	-.0434(8)	0.0858(10)	0.3216(3)
C12	0.0791(8)	0.1382(10)	0.3484(3)
C13	0.0672(8)	0.2301(10)	0.3930(3)
C14	-.0746(8)	0.3714(9)	0.3900(3)
C15	0.0405(8)	0.0828(10)	0.4299(3)
C16	-.1114(9)	-.0283(10)	0.4252(3)
C17	-.4023(10)	0.0080(21)	0.4195(5)

TABLE XIII

HYDROGEN POSITIONAL PARAMETERS FOR  $(C_{17}H_{20}NO_3^+)_2SO_4^-$ 

Atom	x( $\sigma(x)$ )	y( $\sigma(y)$ )	z( $\sigma(z)$ )
HC1	-.079(10)	0.003(17)	0.256(28)
HC2	0.192(7)	-.016(11)	0.238(2)
HC5	0.307(12)	0.293(15)	0.422(3)
HC6	0.186(8)	0.610(10)	0.419(3)
HC7	0.155(6)	0.695(8)	0.338(2)
HC8	-.126(7)	0.563(8)	0.335(2)
HC9	-.322(9)	0.309(11)	0.376(2)
H1C10	-.299(13)	0.034(17)	0.331(3)
H2C10	-.274(12)	0.217(16)	0.316(3)
HC14	-.104(6)	0.424(7)	0.416(2)
H1C15	0.112(6)	0.004(9)	0.431(2)
H2C15	0.039(10)	0.153(13)	0.457(3)
H1C16	-.127(8)	-.086(9)	0.444(2)
H2C16	-.118(9)	-.107(12)	0.399(3)
H1C17	-.456(17)	0.118(21)	0.419(5)
H2C17	-.446(19)	-.077(22)	0.440(5)
H3C17	-.422(13)	-.098(17)	0.395(3)
HO1	0.489(15)	0.240(22)	0.294(4)
HO2	0.440(26)	0.706(34)	0.403(7)
HN1	-.250(7)	0.166(9)	0.450(2)

TABLE XIV  
NONHYDROGEN ATOM ANISOTROPIC THERMAL PARAMETERS  
FOR  $(C_{17}H_{20}NO_3^+)_2SO_4^-$

Atom	$10^3U_{11}$	$10^3U_{22}$	$10^3U_{33}$	$10^3U_{12}$	$10^3U_{13}$	$10^3U_{23}$
S1	31(1)	30(1)	47(2)	0(0)	0(0)	-2(1)
O1	39(3)	89(5)	63(5)	2(4)	15(3)	-15(4)
O2	75(4)	63(4)	42(5)	-39(4)	9(4)	3(3)
O3	35(3)	62(4)	21(4)	-9(3)	2(2)	-16(3)
O111	84(4)	62(3)	28(5)	36(4)	3(4)	-3(3)
O112	73(4)	80(5)	63(6)	-4(4)	-46(4)	11(4)
O113	71(4)	120(4)	67(6)	6(6)	25(4)	39(5)
O114	77(4)	78(4)	59(6)	-28(4)	3(4)	-13(4)
O115	48(4)	130(9)	53(8)	0(0)	0(0)	27(8)
N1	33(3)	34(3)	27(5)	-4(3)	3(3)	0(3)
C1	61(5)	45(4)	13(6)	-2(4)	-6(4)	-3(4)
C2	63(5)	46(4)	26(7)	0(5)	18(5)	-23(4)
C3	50(4)	43(4)	28(7)	9(4)	3(4)	-3(4)
C4	42(4)	40(4)	25(6)	3(3)	-5(4)	-5(4)
C5	38(4)	34(3)	40(7)	-6(3)	-2(4)	-2(4)
C6	63(5)	38(4)	18(8)	-17(4)	15(5)	-9(4)
C7	61(5)	28(3)	42(8)	-9(3)	17(5)	7(4)
C8	59(5)	30(5)	17(7)	11(3)	3(4)	3(3)
C9	35(3)	34(3)	36(7)	5(3)	0(4)	2(4)
C10	36(4)	48(4)	46(8)	-7(4)	-10(4)	3(5)
C11	43(4)	34(3)	8(6)	-3(3)	-7(4)	1(4)
C12	38(3)	24(3)	32(7)	0(3)	-3(4)	1(3)
C13	33(3)	27(3)	28(6)	2(3)	-3(3)	-2(3)
C14	39(3)	26(3)	35(6)	5(3)	8(4)	-9(3)
C15	32(3)	26(3)	49(7)	2(3)	-3(4)	8(4)
C16	46(4)	26(3)	40(7)	-1(3)	-1(4)	6(4)
C17	36(4)	66(6)	94(11)	-18(5)	7(5)	5(7)

TABLE XIV (Continued)

---


$$\exp[-2\pi^2(U_{11}h^2a^{*2} + U_{22}k^2b^{*2} + U_{33}l^2c^{*2} + 2U_{12}hka^{*}b^{*} \\ + 2U_{13}hla^{*}c^{*} + 2U_{23}klb^{*}c^{*})]$$

TABLE XV

HYDROGEN ISOTROPIC THERMAL PARAMETERS FOR  $(C_{17}H_{20}NO_3^+)_2SO_4^-$ 

Atom	$10^2U$	Atom	$10^2U$
HC1	10(3)	HC2	3(2)
HC5	7(3)	HC6	3(2)
HC7	0(1)	HC8	1(2)
HC9	5(2)	H1C10	11(4)
H2C10	8(4)	HC14	0(1)
H1C15	2(2)	H2C15	8(3)
H1C16	2(2)	H2C16	4(2)
H1C17	13(6)	H2C17	14(7)
H3C17	14(4)	H01	16(6)
HO2	47(13)	HN1	2(2)

TABLE XVI

NONHYDROGEN POSITIONAL PARAMETERS FOR  $(C_{18}H_{22}NO_3)^+I^-$ 

Atom	x( $\sigma(x)$ )	y( $\sigma(y)$ )	z( $\sigma(z)$ )
I1	-.0021(1)	0.0475(1)	-.1033(0)
N1	-.0120(10)	0.4205(6)	0.0414(4)
O1	-.1314(11)	0.3350(9)	0.4342(4)
O2	0.3102(11)	0.2273(9)	0.3486(4)
O3	0.1016(8)	0.3857(7)	0.3325(4)
O4	-.3310(12)	0.2172(9)	0.5114(5)
C1	-.2678(11)	0.2522(12)	0.2412(7)
C2	-.2637(11)	0.2782(11)	0.3211(6)
C3	-.1433(12)	0.3208(10)	0.3556(6)
C4	-.0291(9)	0.3492(8)	0.3106(6)
C5	0.1954(10)	0.3723(9)	0.2663(6)
C6	0.2831(12)	0.2557(12)	0.2695(8)
C7	0.2227(11)	0.1417(9)	0.2280(6)
C8	0.1537(11)	0.1540(9)	0.1634(6)
C9	0.0017(15)	0.2870(7)	0.0766(4)
C10	-.1358(11)	0.2423(10)	0.1145(6)
C11	-.1504(10)	0.2816(9)	0.1972(6)
C12	-.0382(9)	0.3343(8)	0.2322(5)
C13	0.0948(9)	0.3757(8)	0.1962(5)
C14	0.1275(10)	0.2845(8)	0.1305(6)
C15	0.0822(10)	0.5093(9)	0.1618(6)
C16	-.0314(10)	0.5140(8)	0.1042(6)
C17	-.1335(14)	0.4272(14)	-.0129(8)
C18	0.1151(12)	0.4487(14)	-.0061(7)



TABLE XVII  
NONHYDROGEN ATOM ANISOTROPIC THERMAL PARAMETERS FOR  
(C<sub>18</sub>H<sub>22</sub>NO<sub>3</sub>)<sup>+</sup>I<sup>-</sup>

Atom	10 <sup>3</sup> U <sub>11</sub>	10 <sup>3</sup> U <sub>22</sub>	10 <sup>3</sup> U <sub>33</sub>	10 <sup>3</sup> U <sub>12</sub>	10 <sup>3</sup> U <sub>13</sub>	10 <sup>3</sup> U <sub>23</sub>
I1	37(0)	40(0)	42(0)	5(0)	-1(0)	-3(0)
N1	30(4)	23(3)	28(4)	4(3)	2(4)	1(2)
O1	80(6)	49(4)	18(4)	-5(5)	11(4)	-4(3)
O2	88(7)	53(5)	26(5)	26(5)	-8(4)	9(4)
O3	38(4)	32(4)	20(4)	-1(3)	-1(3)	-4(3)
O4	85(7)	64(6)	34(5)	-16(6)	14(5)	6(4)
C1	25(5)	39(6)	39(8)	-2(4)	0(4)	-1(5)
C2	33(5)	38(6)	29(7)	3(4)	10(4)	12(5)
C3	40(6)	26(5)	28(6)	10(4)	3(4)	2(4)
C4	35(5)	19(4)	24(5)	1(3)	0(3)	0(3)
C5	28(4)	27(4)	27(6)	-6(3)	-3(4)	3(4)
C6	30(5)	38(6)	31(6)	6(4)	-3(4)	5(5)
C7	43(5)	22(5)	39(6)	11(4)	11(5)	11(4)
C8	43(5)	15(4)	36(6)	5(4)	3(4)	-2(4)
C9	39(4)	17(3)	23(4)	5(5)	5(5)	-5(3)
C10	39(5)	26(5)	32(6)	-7(4)	-1(4)	-1(4)
C11	29(4)	22(4)	27(6)	-1(3)	-1(4)	0(3)
C12	27(4)	18(4)	24(5)	1(3)	2(3)	-3(3)
C13	23(4)	17(4)	26(5)	-7(3)	-3(3)	-5(3)
C14	27(4)	10(4)	33(6)	0(3)	1(4)	0(3)
C15	34(4)	19(4)	35(6)	3(3)	-4(4)	9(4)
C16	45(6)	28(4)	19(5)	16(3)	0(4)	2(4)
C17	52(7)	54(9)	29(7)	3(6)	-15(6)	7(6)
C18	34(5)	37(6)	44(7)	2(5)	2(5)	-2(6)

$$\exp[-2\pi^2(U_{11}h^2a^{*2} + U_{22}k^2b^{*2} + U_{33}l^2c^{*2} + 2U_{12}hka^*b^* + 2U_{13}hla^*c^* + 2U_{23}klb^*c^*)]$$

TABLE XVIII  
HYDROGEN POSITIONAL PARAMETERS FOR  $(C_{18}H_{22}NO_3)^+I^-$

Atom	x	y	z
HC1	-.3511	0.2149	0.2171
HC2	-.3443	0.2633	0.3482
HC5	0.2511	0.4718	0.2722
HC6	0.3737	0.2741	0.2462
HC7	0.2528	0.0588	0.2473
HC8	0.1202	0.0788	0.1374
HC9	0.0000	0.2361	0.0286
H1C10	-.1400	0.1504	0.1115
H2C10	-.2055	0.3199	0.1031
HC14	0.2141	0.3046	0.1148
H1C15	0.0719	0.5718	0.2025
H2C15	0.1771	0.5134	0.1540
H1C16	-.0385	0.5969	0.0795
H2C16	-.1175	0.5142	0.1011
H1C17	0.1061	0.5365	-.0239
H2C17	0.1126	0.3947	0.0518
H3C17	0.1979	0.4429	0.0104
H1C18	-.1223	0.4950	-.0500
H2C18	-.1428	0.3474	-.0432
H3C18	-.2420	0.4418	0.0184
HO1	-.2154	0.3018	0.4819
HO2	0.3333	0.1420	0.3666

TABLE XIX

BOND LENGTHS(A°) AND ANGLES(°) FOR  $(C_{17}H_{20}NO_3^+)_2SO_4^{=}$ 

S1-0111	1.468(6)	S1-0112	1.448(7)
S1-0111'	1.468(6)	S1-0112'	1.448(7)
0111-S1-0111'	108.7(4)	0111'-S1-0112	109.3(4)
0111-S1-0112	111.1(4)	0111'-S1-0112'	111.1(4)
0111-S1-0112'	108.4(4)	0112-S1-0112'	109.3(4)
C1-C2	1.41(1)	C1-C11	1.40(1)
C2-C3	1.38(1)	C3-O1	1.37(1)
C3-C4	1.36(1)	C4-C12	1.39(1)
C4-O3	1.38(1)	C5-O3	1.48(1)
C5-C6	1.56(1)	C5-C13	1.55(1)
C6-O2	1.43(1)	C6-C7	1.48(1)
C7-C8	1.33(1)	C8-C14	1.50(1)
C9-C10	1.54(1)	C9-C14	1.55(1)
C10-C11	1.51(1)	C11-C12	1.37(1)
C12-C13	1.51(1)	C13-C14	1.54(1)
C13-C15	1.53(1)	C15-C16	1.50(1)
N1-C9	1.52(1)	N1-C16	1.52(1)
N1-C17	1.50(1)		
C1-C2-C3	121.9(8)	C2-C3-C4	116.6(8)
C2-C3-O1	123.1(8)	O1-C3-C4	120.1(8)
C3-C4-O3	125.6(7)	C3-C4-C12	122.2(7)
C12-C4-O3	112.2(7)	C11-C12-C4	121.7(8)
C13-C12-C4	110.0(6)	C11-C12-C13	127.0(6)
C10-C11-C12	118.4(7)	C1-C11-C12	116.8(7)
C1-C11-C10	124.5(7)	C11-C1-C2	120.0(8)
C4-O3-C5	108.0(6)	O3-C5-C6	107.6(7)
C13-C5-O3	106.3(6)	C13-C5-C6	110.9(6)
C5-C13-C12	101.6(6)	C6-C13-C14	115.8(6)
C5-C13-C15	110.3(6)	C12-C13-C14	104.9(6)
C12-C13-C15	114.1(6)	C14-C13-C15	109.9(6)
C5-C6-C7	113.3(7)	C5-C6-O2	110.2(6)

TABLE XIX (Continued)

O2-C6-C7	112.8(7)	C6-C7-C8	118.9(8)
C7-C8-C14	119.8(7)	C8-C14-C13	108.4(6)
C8-C14-C9	114.4(7)	C13-C14-C9	106.7(5)
C9-C10-C11	114.4(6)	C10-C9-C14	115.4(6)
N1-C9-C10	112.0(6)	N1-C9-C14	105.6(6)
C13-C15-C16	112.7(7)	C15-C16-N1	109.8(5)
C16-N1-C9	112.0(6)	C17-N1-C9	115.1(8)
C17-N1-C16	110.7(7)		

TABLE XX

BOND LENGTHS (Å) AND ANGLES (°) FOR  $(C_{18}H_{22}NO_3)^+I^-$ 

C1-C2	1.43(2)	C1-C11	1.40(2)
C2-C3	1.38(2)	C3-O1	1.40(2)
C3-C4	1.38(2)	C4-C12	1.39(1)
C4-O3	1.36(1)	C5-O3	1.48(1)
C5-C6	1.50(2)	C5-C13	1.56(1)
C6-O2	1.44(2)	C6-C7	1.53(2)
C7-C8	1.32(2)	C8-C14	1.52(1)
C9-C10	1.55(2)	C9-C14	1.53(2)
C10-C11	1.52(2)	C11-C12	1.36(1)
C12-C13	1.49(1)	C13-C14	1.54(1)
C13-C15	1.54(1)	C15-C16	1.49(1)
N1-C9	1.55(1)	N1-C16	1.50(1)
N1-C17	1.50(2)	N1-C18	1.50(2)
C1-C2-C3	121(1)	C2-C3-C4	119(1)
C2-C3-O1	123(1)	O1-C3-C4	119(1)
C3-C4-O3	129(1)	C3-C4-C12	120(1)
C12-C4-O3	112(1)	C11-C12-C4	123(1)
C13-C12-C4	109(1)	C11-C12-C13	127(1)
C10-C11-C12	118(1)	C1-C11-C12	118(1)
C1-C11-C10	123(1)	C11-C1-C2	118(1)
C4-O3-C5	108(1)	O3-C5-C6	113(1)
C13-C5-O3	104(1)	C13-C5-C6	113(1)
C5-C13-C12	101(1)	C5-C13-C14	117(1)
C5-C13-C15	112(1)	C12-C13-C14	108(1)
C12-C13-C15	111(1)	C14-C13-C15	107(1)
C5-C6-C7	115(1)	C5-C6-O2	108(1)
O2-C6-C7	111(1)	C6-C7-C8	121(1)
C7-C8-C14	120(1)	C8-C14-C13	109(1)
C8-C14-C9	112(1)	C13-C14-C9	107(1)
C9-C10-C11	114(1)	C10-C9-C14	113(1)
N1-C9-C10	112(1)	N1-C9-C14	109(1)

TABLE XX (Continued)

C13-C15-C16	111(1)	C15-C16-N1	113(1)
C16-N1-C9	109(1)	C17-N1-C9	110(1)
C17-N1-C16	112(1)	C17-N1-C18	105(1)
C18-N1-C9	111(1)	C18-N1-C16	110(1)

TABLE XXI  
TORSIONAL ANGLES OF MORPHINE AND MORPHINE SALTS

Dihedral Angle	Morphine methyl- iodide·H <sub>2</sub> O	Morphine sulfate·3H <sub>2</sub> O	Morphine·H <sub>2</sub> O
C5-C6-C7-C8	37.0(16)	46.3(11)	41.5(6)
C6-C7-C8-C14	-2.6(16)	-3.6(12)	-4.9(6)
C7-C8-C14-C13	-38.7(13)	-43.4(10)	-39.0(5)
C8-C14-C13-C5	48.0(10)	48.3(10)	46.9(5)
C14-C13-C5-C6	-17.1(12)	-9.9(11)	-13.1(5)
C13-C5-C6-C7	-24.2(13)	-36.4(10)	-29.2(5)
C15-C13-C14-C9	121.0(9)	122.2(8)	123.4(4)
C13-C14-C9-N1	64.5(9)	65.5(8)	65.0(4)
C14-C9-N1-C16	-59.8(11)	-64.8(7)	-62.5(5)
C9-N1-C16-C15	56.0(11)	57.4(8)	53.3(6)
N1-C16-C15-C13	-57.1(10)	-51.5(9)	-46.2(5)
C16-C15-C13-C14	60.1(10)	55.6(8)	54.1(5)




Figure 44. An ORTEP Projection of Morphinine Sulfate






Figure 45. An ORTEP Projection of Morphin Methyliodide

morphine methyliodide and also morphine hydrate are presented in Figures 46, 47 and 48, respectively

#### The Structure of Compound I

Compound I was crystallized from the chloroform extract of morphine sulfate-Froehde reaction as dark purple rectangular shaped crystals. Mass spectral data (Figure 49) show that the compound has a molecular weight of 279.09 corresponding to a molecular formula of  $C_{17}H_{13}NO_3$ . The X-ray crystallographic study indicated a monoclinic symmetry with a  $P2_1/n$  space group. Unit cell dimensions are 18.591(12), 7.774(4) and 8.762(5) Å and the  $\beta$  angle is 96.60(5)° (Table XXII). 1527 reflections from a total of 3476 points were considered valid (Appendix C). The atomic positional and thermal parameters are presented in Tables XXIII, XXIV, XXV and XXVI. The final bond lengths and bond angles between atoms in Compound I are listed in Table XXVII, with the R factor minimized to 5.8%.

An ORTEP diagram for Compound I is presented in Figure 50, and an ORTEP view of the hydrogen bonding between molecules of the Compound I is shown in Figure 51.

#### Solid State Circular Dichroism of Morphine Methyliodide

The SSCD of morphine hydrate not only shows a distinct difference from the spectrum for morphine sulfate as




Figure 46. An ORTEP Projection of Hydrogen Bonding Network of Morphine Sulfate




Figure 47. An ORTEP Projection of Hydrogen Bonding Network of Morphinethyliodide




Figure 48. An ORTEP Projection of Hydrogen Bonding Network of Morphine Hydrate




Figure 49. Mass Spectrum of Compound I

TABLE XXII  
CRYSTAL DATA FOR  $C_{17}H_{13}NO_3$

Formula	$C_{17}H_{13}NO_3$
MWT	279.09
<u>a</u>	18.591(12) Å <sup>o</sup>
<u>b</u>	7.774(4)
<u>c</u>	8.762(5)
$\beta$	96.60(5) <sup>o</sup>
V	1258(1) Å <sup>o3</sup>
F(000)	587
$\lambda_{Mok_{\alpha}}$	0.71069 Å
$\mu_{Mok_{\alpha}}$	4.86
$D_{calc}$	1.473 g cm <sup>-3</sup>
z	4
Space group	P2 <sub>1</sub> /n
Systematic absences	h0l, h + 1 $\neq$ 2n 0k0, k $\neq$ 2
Obs. reflections	1527
Final R factor	5.8%

TABLE XXIII  
NONHYDROGEN POSITIONAL PARAMETERS FOR  $C_{17}H_{13}NO_3$

Atom	x( $\sigma(x)$ )	y( $\sigma(y)$ )	z( $\sigma(z)$ )
O1	-.0827(2)	-.4205(4)	0.3100(4)
O2	-.0894(2)	0.1079(5)	0.0584(4)
O3	0.0394(2)	-.3324(4)	0.4674(4)
N1	0.1326(3)	0.4277(5)	0.2505(4)
C1	-.0591(3)	0.0064(6)	0.1530(5)
C2	-.0908(3)	-.1599(6)	0.1784(5)
C3	-.0556(3)	-.2655(6)	0.2834(5)
C4	0.0147(3)	-.2213(6)	0.3746(5)
C5	0.1607(3)	-.1101(6)	0.5335(5)
C6	0.2248(4)	-.0499(6)	0.6082(5)
C7	0.2491(3)	0.1166(6)	0.5837(5)
C8	0.2089(3)	0.2231(6)	0.4822(5)
C9	0.1050(3)	0.2720(5)	0.2844(5)
C10	0.0395(3)	0.2159(6)	0.2109(5)
C11	0.0119(3)	0.0555(6)	0.2427(5)
C12	0.0487(3)	0.0578(5)	0.3478(4)
C13	0.1175(3)	-.0044(5)	0.4268(4)
C14	0.1445(3)	0.1637(6)	0.3986(5)
C15	0.2303(4)	0.4074(6)	0.4632(6)
C16	0.2079(3)	0.4689(6)	0.3019(6)
C17	0.0963(5)	0.5334(7)	0.1262(6)



TABLE XXIV  
HYDROGEN POSITIONAL PARAMETERS FOR  $C_{17}H_{13}NO_3$

Atom	x( $\sigma(x)$ )	y( $\sigma(y)$ )	z( $\sigma(z)$ )
HO2	-.053(3)	-.476(8)	0.387(6)
HC2	-.144(2)	-.199(6)	0.109(5)
HC5	0.141(2)	-.228(6)	0.551(5)
HC6	0.257(3)	-.130(6)	0.679(5)
HC7	0.298(2)	0.163(6)	0.640(4)
HC10	0.008(2)	0.283(5)	0.139(4)
H1C15	0.207(3)	0.483(6)	0.530(5)
H2C15	0.284(3)	0.424(8)	0.494(6)
H1C16	0.245(4)	0.407(9)	0.225(8)
H2C16	0.214(2)	0.602(5)	0.293(5)
H1C17	0.104(4)	0.481(8)	0.022(8)
H2C17	0.117(3)	0.645(7)	0.131(5)
H3C17	0.036(5)	0.547(10)	0.136(9)

TABLE XXV  
 NONHYDROGEN ATOM ANISOTROPIC THERMAL PARAMETERS  
 FOR  $C_{17}H_{13}NO_3$

Atom	$10^3 U_{11}$	$10^3 U_{22}$	$10^3 U_{33}$	$10^3 U_{12}$	$10^3 U_{13}$	$10^3 U_{23}$
O1	49(4)	38(2)	58(2)	-10(2)	1(2)	7(2)
O2	45(4)	54(2)	66(2)	-1(2)	-16(2)	22(2)
O3	56(4)	36(2)	50(2)	-6(2)	-3(2)	16(2)
N1	34(4)	32(2)	42(2)	-3(2)	4(2)	10(2)
C1	40(6)	38(2)	38(2)	-2(3)	1(3)	3(2)
C2	29(5)	37(2)	42(2)	-4(3)	1(2)	5(2)
C3	36(6)	34(2)	41(2)	-7(3)	10(3)	0(2)
C4	36(5)	31(2)	32(2)	3(2)	8(2)	6(2)
C5	39(6)	34(2)	36(2)	2(3)	4(2)	10(2)
C6	39(6)	44(3)	34(2)	9(3)	1(2)	8(2)
C7	28(5)	40(2)	34(2)	1(3)	4(2)	-2(2)
C8	27(6)	35(2)	38(2)	2(2)	9(2)	2(2)
C9	34(6)	24(2)	33(2)	1(2)	9(2)	7(2)
C10	30(5)	34(2)	34(2)	7(2)	1(2)	10(2)
C11	31(5)	35(2)	28(2)	2(2)	2(2)	4(2)
C12	22(5)	28(2)	29(2)	4(2)	6(2)	5(1)
C13	35(5)	27(2)	26(2)	1(2)	10(2)	3(1)
C14	24(5)	31(2)	31(2)	3(2)	2(2)	3(2)
C15	51(6)	37(2)	48(2)	-7(3)	2(3)	2(2)
C16	40(6)	40(2)	53(3)	-13(3)	5(3)	16(2)
C17	68(8)	31(2)	45(3)	-2(3)	5(3)	13(2)

$$\exp[-2\pi^2(U_{11}h^2a^{*2} + U_{22}k^2b^{*2} + U_{33}l^2c^{*2} + hka^{*}b^{*} + 2U_{13}hla^{*}c^{*} + 2U_{23}klb^{*}c^{*})]$$

TABLE XXVI  
 HYDROGEN ISOTROPIC THERMAL PARAMETERS FOR  $C_{17}H_{13}NO_3$

Atom	$10^2U$	Atom	$10^2U$
HO1	9(2)	HC2	4(1)
HC5	4(1)	HC6	5(1)
HC7	4(1)	HC10	3(1)
H1C15	6(2)	H2C15	8(2)
H1C16	12(3)	H2C16	4(1)
H1C17	9(2)	H2C17	4(1)
H3C17	13(4)		

TABLE XXVII  
BOND LENGTHS (Å) AND ANGLES (°) FOR C<sub>17</sub>H<sub>13</sub>NO<sub>3</sub>

C1-O2	1.233(6)	C1-C1	1.448(7)
C1-C11	1.505(8)	C2-C3	1.345(7)
C3-O1	1.337(6)	C3-C4	1.493(8)
C4-C12	1.450(7)	C4-O3	1.238(8)
C5-C6	1.374(8)	C5-C13	1.421(6)
C6-C7	1.396(7)	C7-C8	1.372(7)
C8-C14	1.407(7)	C9-C10	1.391(8)
C9-C14	1.443(6)	C10-C11	1.389(7)
C11-C12	1.396(6)	C12-C13	1.443(8)
C13-C14	1.432(6)	C8-C15	1.502(7)
C15-C16	1.505(7)	N1-C9	1.360(6)
N1-C16	1.457(8)	N1-C17	1.466(7)
C1-C2-C3	118.4(5)	C2-C3-C4	123.6(5)
C2-C3-O1	120.4(5)	O1-C3-C4	116.0(4)
C3-C4-O3	115.2(4)	C3-C4-C12	119.1(4)
C12-C4-O3	125.7(5)	C11-C12-C4	118.0(5)
C13-C12-C4	123.5(4)	C11-C12-C13	118.4(4)
C10-C11-C12	122.3(5)	C1-C11-C12	121.2(4)
C1-C11-C10	116.4(4)	C11-C1-C2	119.6(4)
C13-C5-C6	120.7(4)	C5-C13-C12	123.5(4)
C5-C13-C14	117.2(5)	C12-C13-C14	119.3(4)
C5-C6-C7	121.4(5)	C6-C7-C8	119.8(5)
C7-C8-C14	120.4(4)	C8-C14-C13	120.3(4)
C7-C8-C15	121.0(5)	C8-C14-C9	120.3(4)
C13-C14-C9	119.4(4)	C9-C10-C11	121.1(4)
C10-C9-C14	119.4(4)	N1-C9-C10	120.9(4)
N1-C9-C14	119.8(4)	C14-C8-C15	118.4(4)
C8-C15-C16	111.0(4)	C15-C16-N1	111.9(5)
C9-N1-C17	120.6(4)	C9-N1-C16	120.0(4)
C16-N1-C17	116.8(4)		




Figure 50. An ORTEP Projection of Compound I




Figure 51. An ORTEP View of Hydrogen Bonding of Compound I

reported earlier but also from that for morphine methyl-iodide. The salts on the other hand have qualitatively similar SSCD spectra showing a positive band with  $\lambda_{\text{max}}^+$  around 255 nm and a negative band with  $\lambda_{\text{max}}^-$  around 295 to 300 nm. For morphine hydrate the SSCD is all positive with a shoulder around 297 nm and a band with  $\lambda_{\text{max}}^+$  of 260 nm.

The SSCD spectra of these three compounds (Figure 52) are different from the spectra of the compounds in the aqueous solutions. All the spectra of the compounds in distilled water possess the same CD patterns with slightly different values for  $\lambda^0$  (Figure 53). The characteristic wavelengths of these three morphine molecules are tabulated in Table XXVIII for easy comparison.




Figure 52. SS CD Spectra of Morphine (...), Morphine Sulfate (—) and Morphine Methyliodide (---)






Figure 53. CD Spectra of Morphine (···), Morphine Sulfate (—), and Morphine Methyl-iodide (---) in Distilled Water

TABLE XXVIII  
SSCD AND CD DATA OF MORPHINE AND ITS SALTS

Compound	$\lambda_{\text{max}}$ (nm)	
	SSCD	CD
Morphine	297 <sup>+</sup> (sh), 260 <sup>+</sup>	243 <sup>+</sup> , 285 <sup>-</sup> , 270 <sup>o</sup>
Morphine Sulfate	302 <sup>-</sup> , 254 <sup>+</sup> , 271 <sup>o</sup>	244 <sup>+</sup> , 284 <sup>-</sup> , 267 <sup>o</sup>
Morphine Methyl- iodide	294 <sup>-</sup> , 256 <sup>+</sup> , 278 <sup>o</sup>	244 <sup>+</sup> , 285 <sup>-</sup> , 264 <sup>o</sup>

+ indicates a positive band  
 - indicates a negative band  
 o indicates a crossing point  
 (sh) indicates a shoulder band

## CHAPTER V

### DISCUSSION

The purpose of this research was to study the oxidation-reduction reactions of morphine in three different media: basic, neutral and concentrated sulfuric acid solutions. Product identification is the critical theme of the work and is described for each system in turn.

#### Oxidation of Morphine in Basic Solution

The principal product of this reaction proposed previously (11) is the dimeric compound pseudomorphine. Working from this assumption a NMR study was done for structural confirmation.

Results from a  $^1\text{H}$  NMR of morphine sulfate and pseudomorphine in  $\text{DMSO}(d_6)$  show that the only significant changes in these two spectra occur in the region of the aromatic hydrogens. Since there are only two aromatic hydrogens in morphine the spectrum consists of a quartet of peaks. The pseudomorphine spectrum on the other hand consists of a singlet in this region. The integrations of those peaks confirm that the relative numbers of aromatic hydrogen atoms are two for morphine and one for pseudomorphine. The implication is that dimerization occurs via the aromatic

rings, by a combination yet undetermined but involving C1 or C2 and C1' or C2', a total of three possibilities.

Both compounds contain seventeen different carbons but only those in the region of the unsaturated carbons in the  $^{13}\text{C}$  NMR spectra show any significant difference in the positions of some of the chemical shifts and off-resonance peaks. This is further confirmation of dimerization via the aromatic nuclei, but it still does not reveal which carbons are involved in the link. The possibility of a C1 to C2' bond is, however, eliminated from consideration. The remaining two possibilities for dimerization are the symmetrical linkages between C1 and C1' or C2 and C2' (Figure 54).

To distinguish between these two possible structures, comparisons were made of the  $^{13}\text{C}$  NMR spectra for morphine, codeine, and simpler model compounds, namely phenol and o,o'-biphenol and anisole and m,m'-bianisole. Assignments of the chemical shifts for C1 and C2 in morphine sulfate were confirmed in a comparison with codeine. The methyl group on C3 of codeine produces a large ortho-effect on C2 causing it to move upfield relative to the C2 resonance in morphine sulfate. Accordingly the peaks at 119.05 ppm and 117.10 ppm in the fully decoupled  $^{13}\text{C}$  spectrum of morphine are assigned to the C1 and C2 resonances, respectively. The corresponding positions in codeine are 118.15 ppm and 112.91 ppm. The assignments are in accord with the results of Carroll et al. (43).




Figure 54. The Possible Structures of Pseudomorphine

The  $^{13}\text{C}$  NMR fully decoupled and off-resonance spectra of pseudomorphine identify the resonance at 121.0 ppm with a tertiary carbon atom and the peak at 124.3 ppm with a quaternary carbon atom, but does not distinguish which one is C1 and which is C2. In forming the quaternary carbon the aromatic proton has been substituted by a monomeric morphine. The chemical shift for the carbons forming the link would be expected to move to lower field compared to the corresponding carbon in morphine on substitution. This was confirmed from the C spectra of the four model compounds (Figure 55).

In 2 the C2 resonance peak was shifted downfield by 10.4 ppm on substitution relative to C2 on phenol. In addition the position of C1 on phenol was moved upfield and C6 was moved downfield on substitution to produce the dimer. For the anisole analog of the 1,1'-dimer, C3 shifted downfield by 12.35 ppm from the position on anisole on substitution. By analogy with these results it is concluded that a downfield shift of approximately 10 ppm for the bridging carbons might be anticipated when morphine is oxidatively dimerized to pseudomorphine. The largest downfield shift observed when the  $^{13}\text{C}$  spectra of morphine and pseudomorphine are compared is 7.2 ppm. This would require a reversal in the assignments of the C1 and C2 resonance peaks in the spectrum of pseudomorphine compared to morphine and is the most reasonable conclusion based upon the data for the model compounds. Bonding




Figure 55. Structures of Phenol, o,o'-Biphenol, Anisole and m,m'-Bianisole

between C2-C2' is consistent with this interpretation. The shifts observed for C3, C11 and C12 are also consistent with a new substituent at C2. In contrast, if the dimerization position was C1-C1', the ortho-effect of the hydroxyl group would produce an upfield shift for both C2 and C11 contrary to what is observed.



In conclusion, therefore, the evidence from the  $^{13}\text{C}$  NMR study supports the structural model in which the dimerization occurs at the carbon C2-C2' positions.

#### Oxidation of Morphine and Pseudomorphine by Inorganic Ions in Aqueous Solutions


In the previous section the oxidation of morphine in basic solution by either air or by ferricyanide ion was described. The only product is pseudomorphine if the oxidant is present in stoichiometric proportions.

The presence of pseudomorphine as a product of the oxidation reaction is easily confirmed from the CD spectrum of the solution. For each of the reactions in this section the product was again confirmed to be pseudomorphine from the measured CD spectrum, if stoichiometric quantities were used. The higher oxidation product (HOP), yet uncharacterized, is produced whenever the oxidant is present in large excess. HOP can also be distinguished by its CD spectrum. Accordingly the contributing half reactions may be formally described by the equations:





and



OA and RA are the oxidized and reduced forms of the inorganic reagent, and M and M<sub>2</sub> are morphine and pseudomorphine, respectively. Oxidative dimerization is a two electron transfer process of uncertain mechanism. It is represented here as a basic reaction involving OH<sup>-</sup> since the reaction readily proceeds in aqueous alkali, but not in aqueous acid. Other alternatives may be possible, such as a free radical reaction or a hydride transfer, but are considered to be less likely.

Eight potential oxidizing agents were selected whose standard reduction potentials (E<sup>0</sup>) ranged over 1.56 volt, and used to estimate reduction potentials for the morphine and pseudomorphine half reactions (Table XXIX). The morphine column in Table XXIX describes the results observed when 1:10 mole ratio of morphine to inorganic ions were used. For the last column, the starting material was pseudomorphine and the OA was present in large excess. Since the drug concentrations were on the order of 5\*1E-5 M, the E<sup>0</sup> data alone is not informative. Corrections to the E values in Table XXIX were made using the Nernst equation and the simple assumption that:

a) the oxidation products are the same for all the reactions.

TABLE XXIX  
REDUCTION POTENTIAL DATA OF METAL IONS IN  
THE NEUTRAL REDOX REACTIONS

Electrode Reaction	$E^{\circ}$ (volt)	Result	
		Morphine	Pseudo-morphine
$\text{Ag} + e^{-} \rightleftharpoons \text{Ag}$	0.799	R	R
$\text{MnO}_4^{-} + e^{-} \rightleftharpoons \text{MnO}_4^{-2}$	0.564	R	R
$\text{Fe}(\text{CN})_6^{-3} + e^{-}$ $\downarrow$ $\text{Fe}(\text{CN})_6^{-4}$	0.36	R	R
$\text{Cu}^{+2} + 2e^{-} \rightleftharpoons \text{Cu}$	0.337	NR	(R)
$\text{Pb}^{+2} + 2e^{-} \rightleftharpoons \text{Pb}$	-.126	NR	NR
$\text{Ni}^{+2} + 2e^{-} \rightleftharpoons \text{Ni}$	-.250	NR	NR
$\text{Cd}^{+2} + 2e^{-} \rightleftharpoons \text{Cd}$	-.403	NR	NR
$\text{Zn}^{+2} + 2e^{-} \rightleftharpoons \text{Zn}$	-.763	NR	NR

NR = non-spontaneous reaction  
R = spontaneous reaction  
(R) = reaction goes very slowly

b) at equilibrium, oxidation is effectively quantitative.

c) with OA in large excess,

$$[\text{OA}]_{\text{final}} \cong [\text{OA}]_{\text{initial}}$$

The reduction potential for the morphine half reaction is intermediate between that for ferricyanide and for Cu(II) ion; that is, between 0.42 volt and 0.24 volt after correction for concentration. For pseudomorphine there was some evidence for oxidation by Cu(II) ion, but the process was kinetically very slow, CD evidence for HOP developing over a period of several days. The lack of availability of an OA whose  $E^{\circ}$  value lies between 0.33 and -0.1 volt prevented a better estimate being obtained.

The results are consistent with the reported  $E^{\circ}$  values for oxidizing agents which function in basic media, such as  $\text{MnO}_4^-$  ( $E^{\circ} = 0.56$  volt) and dissolved oxygen ( $E^{\circ} = 0.41$  volt).

#### Oxidation of Morphine Sulfate and Apomorphine Hydrochloride in Concentrated Sulfuric Acid Solution

Concentrated sulfuric acid has excellent solvent properties, is a powerful dehydrating agent, a strong acid, a good oxidizing medium, and has the ability to stabilize free radicals (44). Even in the absence of the inorganic oxyanions used as oxidizing agents in the color

test reagents, the solvent alone, or in combination with dissolved oxygen, could produce a number of concurrent and/or consecutive reactions with morphine and apomorphine. Morphine has a pentacyclic structure and apomorphine is tetracyclic. There is a possibility that the latter is an intermediate in the reaction of morphine.

There is an immediate reaction on solution of morphine in concentrated sulfuric acid evidenced by inversion of the CD spectrum from that observed in dilute aqueous acid. Dilution of the concentrated sulfuric acid with water does not reverse the CD spectrum indicating that an intermediate product has been formed. There is no evidence for free radicals in the solution (45). Apomorphine on the other hand shows no change in the CD spectra between dilute and concentrated acid conditions.

As time progresses the morphine reaction sequence changes the color of the solution from brownish-orange to green to brownish-gray and finally to brownish-purple. For apomorphine the color sequence is from clear to purple gray and finally to brownish pink. The color of both solutions are very similar after two to four weeks when the reactions are considered to be complete. This does not mean that the compositions of both are identical but it strongly supports the idea that apomorphine is an intermediate in the morphine reaction.

The reactions which occur when the inorganic oxyanions are added appear to be similar to those already

described and the analyses of final solutions give comparable results. Reaction times are much shorter. The evidence supports the involvement of oxygen in the reactions when the compounds are reacted in only the solvent. For these reactions the ESR evidence supports the existence of free radical intermediates when the inorganic oxyanions are present, but only for morphine and not apomorphine (45).

Changes in the UV and CD spectra with time lend further support to the mechanism which includes apomorphine as an intermediate. The spectra of the final solutions under both reaction conditions are not exactly alike in quantitative terms but this might be explained as the same products being present in different amounts. For the color test reagents the reduced forms of the oxyanions may also contribute to the observed UV-visible spectra.

Chloroform extracts from the final morphine products in both acid alone and with color reagents added produce the same UV-visible spectrum. The same is true for the apomorphine reactions. Chloroform does not extract the inorganic ions. Optical activity is lost since none of the extracts give a CD spectrum.

Separation by TLC produces three spots in every case indicating the presence of at least three components. According to the mass spectrometric data, the heaviest component has a molecular weight equal to 325.25 g/mole.

Compound I has a molecular weight of 279.1 g/mole, and its structure was confirmed by X-ray diffraction. The  $R_f$  value of 0.34 is consistent with its being more polar than the other two components. Although there is some correspondence between the  $R_f$  values of the middle component and that for apomorphine alone, the product can not be unreacted apomorphine because the CD result is negative.

A possible mechanism for the formation of Compound I which includes apomorphine as an intermediate is given in Figure 56.

### X-ray Structural Analysis

#### Structural Comparison of Morphine

#### Sulfate and Morphine Methyliodide

Two reasons existed for initiating the solid state structural study of the morphine derivatives. First, the CD spectra of morphine and the sulfate salt are inverted in sign when the aqueous solution conditions are changed from acid to base. This implies a chirality change in the molecule which can not be interpreted without an x-ray structural analysis. Morphine sulfate is readily crystallized into a suitable form for diffraction studies. All attempts to prepare a metal morphinate salt so far have failed. The methyliodide derivative was prepared to explore the effects of N-substitution on the conformation of the four chiral center backbone of the morphine structure.




Figure 56. Postulated Mechanism of Formation of Compound I

The second reason for this study relates back to earlier work on solid state CD (31). In that study the CD spectra of morphine and morphine sulfate were seen to be of opposite signs. The rotational change must be either a consequence of a molecular conformational change or a change in the packing arrangement in the crystal, or perhaps both. A structural determination for morphine free base had already been reported (46); this is the first report for morphine sulfate and morphine methyl-iodide.

From the X-ray data all the bond distances and bond angles for morphine methyl iodide monohydrate (I) and morphine sulfate trihydrate (II) are similar. In (I) the average bond angles about the nitrogen atom bonded to four carbon atoms show no significant distortion from the  $109^\circ$  angle expected for  $sp^3$  hybridization. In (II) where the nitrogen is bound to one hydrogen and three carbons, the C-N-C angles show a larger average ( $112.6^\circ$ ) which may be expected considering the volumes of carbon and hydrogen and the expected compression of the C-N-H angles.

In both molecules ring C displays a distorted boat geometry and ring E is in the chair conformation. The conformation for ring D is the same for both salts. In (II) the methyl group, C17, occupies an equatorial position relative to ring E, the axial position being occupied by a hydrogen from  $H_2SO_4$ . Comparisons may be drawn between these two salts and the hydrochloride salt



(47). While structural details for the rest of the carbon skeleton are similar to those of the title structures, those of ring E of the hydrochloride more closely resemble those of the sulfate salt (II) than the methyliodide salt (I).

The parking arrangement of the molecules in the methyliodide analog involves a linear chain of hydrogen bonding which extends throughout the crystal parallel to the a axis of the unit cells by virtue of the phenolic hydrogen (H01) being bonded to the oxygen of the neighboring water (O4) molecules (H01-O4 distance equal to 1.517(11) Å). The water hydrogens bridge between the phenolic oxygen (O11i) and the hydroxyl oxygen (O2(i)) of another adjacent morphine molecule (O4-O1(i) = 3.08(2) Å, O4-O2(i) = 2.86(1) Å) which is symmetry related by  $x, y, z$  and  $(1/2) + x, (1/2) - y, 1 - z$ . The network then extends to a morphine at  $1 + x, y, z$  and is repeated throughout the crystal which traces out a clockwise turn on progression from the origin along the a direction. There is no hydrogen bonding involving the iodide ion.

On the other hand, the sulfate anion of morphine sulfate trihydrate is involved in hydrogen bonding. An oxygen (O111) of the sulfate group is hydrogen bonded to a nitrogen proton (HN1) of one morphine with a hydrogen bonded distance of 1.71(6) Å, and a second sulfate oxygen (O112) is hydrogen bonded to a hydroxyl hydrogen (HO2(i)) of an adjacent morphine molecule (O112-HO2(i) =

2.39(15) Å) with symmetry related by  $x, y, z$  to  $1+x, 1-y, 1-z$  which then in turn continues in a chain in the  $a$  direction. An oxygen from water (O113) is hydrogen bonded to a phenolic hydrogen (H01) at a hydrogen bonded distance of 1.87(15) Å and one hydrogen from the water molecule is bonded to a hydroxyl oxygen (O2) of the morphine. The hydrogen for the water molecule is not actually located but the observed O113-O2 distance is 2.862(9) Å. A clockwise turn similar to that in the methyliodide salt is again observed to occur along the  $a$  directional axis. The other two water molecules (O114, O115) in morphine sulfate trihydrate are within hydrogen bonding distance from each other (O114-O115 distance equal to 2.783(9) Å) but do not reach contact distance with the morphine molecules.

To complete the comparison of the X-ray structural information for the salts and morphine free base, the data from the earlier work by Bye on morphine monohydrate were used to obtain comparable projections. The conformational arrangements reported here for ring D and E are repeated in the morphine free base structure. The influence of the N1 substitution on the local geometry is negligible. Two types of hydrogen bonding networks exist in morphine monohydrate crystals. One network consists of strong hydrogen bonded interactions along the  $b$  axis and displays a clockwise turn as viewed down the  $b$  axis from the origin. The atoms involved are a nitrogen (N1) on one morphine molecule

and the hydrogen of a phenolic group (H01(i)) on another morphine molecule (N1-H01(i) distance = 1.57(3) Å).

These are symmetry related by  $x, y, z$  and  $-x, (1/2) + y, (1/2) - z$ . So the chain of morphine molecules directly hydrogen bonded to other morphine molecules extends along the  $y$  axis. The second network of hydrogen bonding involves a water molecule and is relatively weak when compared to the first network. The hydrogens (H104 and H204) of the water molecule bridge between the hydroxyl (O2) and phenolic oxygen (O1) atoms on the same morphine unit. The hydrogen bonding distances H104-O2 and H204-O1 are equal to 1.97(4) Å and 2.27(5) Å, respectively, and are typical of rather weak interactions. The oxygen of the same water molecule is also bonded to a hydroxyl (H02(i)) of another morphine molecule and is symmetry related by  $x, y, z$  and  $(-1/2) + x, (1/2) - y, 1 - z$ . The O4-H02(i) distance in this case is 1.91(3) Å. This network extends throughout the crystal along the  $a$  axis.

Considering these structural observations collectively, the reason for the inversion of the negative 286 nm Cotton band in the SSCD spectrum of morphine free base is difficult to interpret. One can eliminate from consideration the idea that the free base has a unique molecular conformation since all are essentially alike. Similarly, a conformation common to all three excludes from consideration any structural changes which might have accompanied the chirality changes around the N1 atom. Long range


crystal packing arrangements have analogous clockwise patterns established by hydrogen bonded networks in all three crystals.

Two differences exist but how these might affect the chirality in a SSCD measurement is not clear. First, the details of the primary hydrogen bonding networks are different in each case. For the salts either the sulfate ion or water is involved in bonding to the morphine moiety. For the free base the morphine molecules are directly bonded. The secondary network, which is much weaker involves both inter- and intramolecular hydrogen bonding. Secondly, the hydrogen bonding helical arrangement in the salts occurs along the short axis of the unit cell, while in the free base the helix is oriented along the long axis. Although the rotation of the helix is clockwise in every case, the pitch of the spiral is longer in the free base. It is known for example that the CD activity of cholesteric liquid crystals is affected by the pitch of the helical arrangement (51). It is also conceivable that the anisotropic crystal arrangements of the aromatic chromophores might orient the direction of polarization of the electronic transitions at different angles relative to the incident beam. The sign of a CD signal is determined by the difference in the molar absorbances of the left and right circularly polarized components. A change in sign might reflect a spatial effect on the relative rotatory strengths of the transitions. Anisotropy is lost


when the crystals are dissolved in polar and non-polar isotropic solvents.

### The Structure of Compound I

The numbering system adopted for the non-hydrogen atoms of Compound I is either based upon that for aporphine compounds (48) or 4H-dibenzo(de,g)quinoline (49) as shown in Figure 57. The molecule shows virtual coplanarity for all non-hydrogen atoms except C15, C16 and C17. The average deviation from the plane of C1 to C14, N1, O1, O2 and O3 is  $0.04061 \text{ \AA}^0$ . Atoms C15, C16 and C17 on the nitrogen ring deviate above the plane  $0.3143$ ,  $0.3597$  and  $0.1812 \text{ \AA}^0$ , respectively. Since this molecule is derived from morphine, the numbering and arbitrary naming of the rings as A, C, D and E is based on the morphine system. In ring A, C2-C3 ( $1.345 \text{ \AA}^0$ ) and C12-C11 ( $1.396 \text{ \AA}^0$ ) bond distances are typically aromatic (cf benzene) characteristic whereas C1-C11 ( $1.505 \text{ \AA}^0$ ) and C3-C4 ( $1.493 \text{ \AA}^0$ ) bond distances are characteristic of normal C-C single bond ( $1.53 \text{ \AA}^0$ ). The C-O bond distances are consistent with a C=O grouping on C1 and C4 and C-OH on C3 (C-O distances equal to  $1.233$ ,  $1.238$  and  $1.337 \text{ \AA}^0$ , respectively). The bond angles in the ring are in the range of  $118.0$  to  $123.6^\circ$ , close to the characteristic value for  $sp^2$  hybridization. Ring C and D are typically aromatic. Ring E, which by deduction is the only non-aromatic ring has a N1-C16 bond distance equal to  $1.457 \text{ \AA}^0$  which approximates to the



APORPHINE



4H-DIBENZO [de,g] QUINOLINE

Figure 57. The Numbering System of Compound I

normal C-N value ( $1.475 \text{ \AA}$ ) and the bond distance for C15-C16 ( $1.505 \text{ \AA}$ ) is typical of a C-C single bond. The bond angles C15-C16-N1 ( $111.9^\circ$ ) and C8-C15-C16 ( $111.0^\circ$ ) are also representative of  $sp^3$  hybridization as expected.

The molecules are interconnected to one other by hydrogen bonds between the hydroxyl hydrogen on C3 of one molecule and the carbonyl oxygen on C4 of another molecule. The O-H hydrogen bond distance is  $2.09 \text{ \AA}$ .

### Conclusion and Suggestions

#### for Further Work

From literature reports (3, 11) and the experimental results of this study, the reactions of morphine in both aqueous neutral and basic solutions with the oxidizing agents  $K_3Fe(CN)_6$ ,  $KMnO_4$ , and some metal ions produce the same products. The principal oxidation product is pseudomorphine which can undergo further oxidation when the conditions allow. The second oxidation product HOP can be differentiated by its CD spectrum from pseudomorphine. The details of the second oxidation and the identity of HOP are unknown. The uncertainty about the mode of dimerization of two morphine molecules to produce pseudomorphine has been tentatively resolved and is believed to occur by bridging between the C2-C2' carbon atoms of the aromatic rings. This conclusion is in agreement with the normal coupling of phenol compounds which usually takes place at positions which are either ortho or para to the -OH

positions (50). The stoichiometry of the reduction-oxidation reaction suggests that two protons are lost and two electrons are transferred in forming one molecule of pseudomorphine from two molecules of morphine.

On the other hand, the oxidation of morphine in concentrated sulfuric acid, a common solvent in color test reagents, produces totally different products. The evidence from CD, UV-visible and TLC data, strongly suggests that the products of the reactions in either the acid alone or in the presence of the color test reagents are the same, and that apomorphine is an intermediate. The reactions occur much more rapidly when the molybdate (Froehde), vanadate (Mandeline), and selenate (Mecke) reagents are used. One of the products is a quinoline type compound (Compound I) which contains the same numbers of C, N and O atoms as morphine. The structure of Compound I shows that the N-heterocyclic ring has been opened and reformed to produce a tetracyclic almost planar structure. The oxygen of the dihydrofuran ring becomes quinonoid in ring A.

The X-ray studies of morphine sulfate and morphine methyliodide to investigate the effect on the morphine C-C skeleton change on substitution at N1 indicated that the two molecules contain similar bond angles and distances. The only difference observed is in the C-N-C average bond angle which is greater in morphine sulfate than in morphine methyliodide as expected due to the compression of the



C-N-H angle in the sulfate salt. There is no significant conformational change between these two structures.

The hydrogen bonding networks in the packing diagram of morphine sulfate and morphine methyliodide show no significant differences but both of these are different from the main hydrogen bonding network of morphine hydrate. The structural data do not differ significantly and therefore cannot be used to interpret the unique SSCD spectrum of morphine free base.

In the future the first priority in the pseudomorphine study is to obtain a suitable single crystal for a definitive X-ray study of the structure of the dimer to confirm the present tentative structural assignment. For the color reactions the remaining two products must also be isolated for characterization to support the proposed mechanism. In addition a complete ESR study needs to be done to characterize the preliminary free radical intermediates. Numerous other empirical color reactions are worthy of study including other opiates, tetracyclines, barbiturates, marijuana and cocaine. The objectives would be to evaluate their specificity and to use the information to develop other better screening reagents.

The field of solid state CD is new and combined with X-ray diffraction, could be extremely important in understanding the crystal forces which are responsible for bulk chirality, and contribute to the knowledge of solid state physics in the interpretation and development of polarizing

materials. A related compound of immediate interest is heroin (3,6-diacetylmorphine) whose SSCD is entirely negative (31). Hydrogen bonding between the heroin molecules without involving the anion or water of crystallization would be consistent with the structural data in morphine free base and would confirm the importance of direct bonding on the SSCD spectra.

## BIBLIOGRAPHY

1. Allinger, N. L., Cava, M. P., Jongh, D. C., Johnson, C. R., Lebel, N. A. and Stevens, C. L., "Quinoline and Isoquinoline", in Organic Chemistry, Worth Publisher, Inc., New York, 1971, Chapter 28.
2. Serturner, Trommsdorff's Journal der Pharmazie, 13, 234 (1805).
3. Bentley, K. W., "The Chemistry of Alkaloids", The Clarendon Press, Oxford, London, 1954.
4. Manske, R. H. F., in "The Alkaloids", Vol. 2, Manske, R. H. F. and Holmes, H. L. (Eds.), Academic Press, New York, 1952.
5. Kerbosch, Arch. Pharm., 248, 536 (1910).
6. Broockmann and Polstorff, Ber., 13, 921 (1880).
7. Pelletier, J. Ann. Chem. U. Pharm., 16, 27 (1935).
8. Fulton, C. C., Am. Jour. Pharm., 1933, 511 (1933).
9. Bentley, K. W., "The Alkaloids", The Clarendon Press, Oxford, London, 1952.
10. Balls, A. K., J. Biol. Chem., 71, 537 (1926).
11. Bentley, K. W. and Dyke, S. F., J. Chem. Soc., 1959, 2574 (1959).
12. Small and Turnbull, J. Am. Chem. Soc., 59, 1541 (1937).
13. Danvin, W. D. and Cone, E. J., J. Pharm. Sci., 69(3), 253 (1980).
14. Fulton, C. C., Am. Jour. Pharm., 105, 503 (1933).
15. Siek, T. J., J. Forensic Sci., 21; 525 (1976).
16. Siek, T. J., J. Forensic Sci., 20, 18 (1975).

17. Bubley, C. and Leo, L., Bull. Narcotics (U.N. Dept. Social Affairs), 7(1), 1955.
18. Terui, T., Tori, K., Maedo, S. and Sawa, Y. K., Tetrahedron Lett., 33, 2853 (1975).
19. Okuda, S., Yamaguchi, S., Kawazoe, Y. and Tsuda, K., Chem. Pharm. Bull., 12(1), 104 (1964).
20. Carroll, F. I., Moreland, C. G., Brine, G. A. and John, A. K., J. Org. Chem., 41(6), 996 (1976).
21. Anders, M. W. and Mannung, G. T., Anal. Chem., 34, 730 (1962).
22. White, M. W., J. Chromatogr., 178, 229 (1979).
23. Clark, C. C., J. Forensic. Sci., 22, 418 (1977).
24. Crone, T. A., Purdie, N., Anal. Chem., 53(1), 17 (1981).
25. Fulton, C. C., Microchem. J., 6, 51 (1962).
26. Talor, J. F., in "Narcotic Drugs", D. H. Clouet (Ed.), Plenum Press, New York, 1971, p. 38.
27. Stewart, C. P. and Stolman, A., "Toxicology", V. 2, Academic Press, New York, 1961, p. 229.
28. Fulton, C. C. and Nebraska, O., J. Lab. and Clin. Med., 1928, 13 (1928).
29. Schiesen, D. W., J. Pharm. Sci., 53(8), 909 (1964).
30. Ahlers, G. and Auterhoff, H., Arch. Pharmaz., 308, 650 (1975).
31. Bowen, J. M. and Purdie, N., Anal. Chem., 52(3), 573 (1980).
32. Fresnel, A., Ann. Chim. Phys., 28, 147 (1825).
33. Condon, E. U., Altar, W., and Eyring, H., J. Chem. Phys., 5, 753 (1937).
34. Cotton, A., Ann. Chim. Phys., 8, 347 (1896).
35. Mitchell, S., "The Cotton Effect", G. Bell, London, 1933.
36. Bragg, W. H. and Bragg, W. L., "X-rays and Crystal Structure", George Bell & Sons, Ltd., London, 1924.

37. International Tables for X-ray Crystallography, The Kynock Press, Birmingham, England, 1965.
38. Small, L. and Lytz, R. E., J. A. C. S., 56, 1928 (1934).
39. Patterson, A. L., Z. Krist., A90, 517 (1935).
40. Main, P., Fiske, S. J., Hull, S. E., Lessinger, L., Germain, G., Declerg, J. P. and Woolfson, M. M., University of York, England, 1980.
41. Stewart, J. M., Ed., The X-ray System Version of 1980, Technical Report TR 446 of the Computer Center, University of Maryland, College Park, Maryland, 1980.
42. Cromer, D. T. and Mann, I. B., Acta Cryst., A24, 321 (1968).
43. Breitmaier, E., Haas, G. and Voelter, W. "Atlas of Carbon-13 NMR Data", London, 1979.
44. Liler, M., "Reaction Mechanisms in Sulfuric Acid", Academic Press, London, 1971.
45. Ketkeaw, N., Ph.D. Thesis, Oklahoma State University, Stillwater, Okla., 1982.
46. Bye, E., Acta. Chem. Scand. B, 30(6), 6 (1976).
47. Gylbert, L., Acta. Crystallogr. B, 29, 1630 (1973).
48. Gulland, Haworth, J. Chem. Soc., 1928, 582 (1928).
49. Patterson, A. M., Capell, L. T. and Walker, D. F., "The Ring Index", 2nd Ed., American Chemical Society, Washington, D. C. 1960, p. 704.
50. Finar, I. L., "Organic Chemistry", Vol. 1, 6th Ed., Longman Group Limited, London, 1973, p. 707.
51. Sackmann, E. and Mohwald, H., J. Chem. Phys., 58(2), 5407 (1973).

## APPENDICES

## APPENDIX A

### X-RAY DATA OF MORPHINE SULFATE

0.0.L			0.3.L			7	15	36	15	138	139
4	845	934	1	270	281	8	167	159	16	52	35
6	366	300	2	225	241	9	62	75	17	103	92
8	977	921	3	49	47	10	141	156	18	18	33
10	923	908	4	288	281	11	223	224	19	56	42
12	337	389	5	647	647	12	45	52	20	18	23
14	712	699	6	776	767	13	127	126			
16	76	120	7	480	502	14	205	223	0.3.L		
18	534	486	8	232	242	15	145	120			
20	45	58	9	341	302	16	62	38	3	212	208
			10	442	444	17	16	28	1	18	35
			11	98	129	18	154	130	2	18	10
			12	100	95	19	136	121	3	94	84
			13	97	57	20	71	84	4	40	58
			14	393	392				5	38	45
			15	76	86	0.6.L			6	60	46
			16	62	58	3	103	109	7	31	40
			17	103	99	1	25	43	8	120	117
			18	16	30	2	129	141	9	118	120
			19	192	168	3	62	73	10	45	13
			20	127	148	4	248	246	11	36	22
						5	54	40	12	18	20
			0.4.L			6	134	134	13	43	3
			0	114	117	7	116	109	14	94	85
			1	297	314	8	125	142	15	52	7
			2	205	200	9	25	38	16	36	4
			3	78	69	10	87	85	17	52	48
			4	384	393	11	192	189	18	29	43
			5	156	136	12	187	198	19	43	36
			6	14	20	13	16	2	20	38	11
			7	58	59	14	114	92	0.9.L		
			8	250	248	15	176	198			
			9	45	4	16	18	29	1	49	15
			10	129	122	17	71	71	2	18	13
			11	185	185	18	118	126	3	45	40
			12	103	107	19	38	20	4	34	46
			13	71	71	20	136	127	5	35	75
			14	172	179				6	80	69
			15	123	110	0.7.L			7	89	93
			16	232	197	1	141	137	8	43	28
			17	40	59	2	196	186	9	43	28
			18	94	95	3	54	107	10	58	64
			19	60	45	4	60	49	11	18	1
			20	65	67	5	25	33	12	20	3
			0.5.L			6	60	33	13	56	34
			1	183	183	7	65	73	14	20	36
			2	107	99	8	52	22	15	89	71
			3	69	73	9	167	174	16	74	67
			4	47	59	10	45	34	17	40	2
			5	319	298	11	114	115	18	52	37
			6	45	38	12	27	4	19	20	6
						13	137	133	20	63	31
						14	18	35			



0.10.L			6	47	11	15	38	8	1.2.L		
0	20	21	7	25	7	16	40	1	0	118	112
1	40	35	8	25	3	17	40	3	1	526	519
2	43	23	9	49	10	18	38	13	2	413	442
3	38	9	10	49	20	19	38	6	3	127	116
4	85	80	11	45	6	1.0.L			4	279	284
5	47	45	12	25	7	0	513	524	5	379	361
6	20	18	13	25	11	1	47	111	6	147	149
7	20	27	14	25	11	2	243	274	7	500	481
8	20	15	15	65	13	3	677	925	8	230	203
9	09	63	16	25	23	4	557	909	9	268	249
10	40	19	17	25	9	5	20	75	10	558	550
11	45	5	18	23	19	6	152	134	11	531	525
12	65	34	19	58	17	7	203	164	12	216	198
13	63	7	20	23	7	8	132	119	13	248	250
14	37	68	0.13.L			9	562	530	14	428	424
15	47	33	1	36	8	10	460	444	15	228	224
16	20	24	2	36	9	11	642	624	16	65	54
17	65	31	3	36	31	12	578	581	17	214	205
18	20	3	4	34	4	13	219	216	18	67	44
19	60	44	5	34	13	14	49	56	19	89	92
20	43	6	6	34	10	15	234	225	20	310	309
0.11.L			7	34	3	16	125	118	1.3.L		
1	23	24	8	34	17	17	143	123	0	65	68
2	78	45	9	34	3	18	207	192	1	170	192
3	20	8	10	34	6	19	125	137	2	216	218
4	20	1	11	56	17	20	36	36	3	393	410
5	20	10	12	31	14	1.1.L			4	147	138
6	20	18	13	31	14	0	105	82	5	502	507
7	38	55	14	31	10	1	406	417	6	241	231
8	74	26	15	31	7	2	524	519	7	417	409
9	20	8	16	31	14	3	366	359	8	174	173
10	20	8	17	45	1	4	1220	1213	9	259	242
11	43	17	18	31	16	5	531	527	10	207	201
12	40	18	19	31	3	6	660	678	11	125	100
13	56	36	20	31	1	7	506	469	12	127	134
14	56	5	0.14.L			8	283	289	13	87	78
15	20	3	0	45	2	9	319	293	14	125	133
16	43	34	1	45	9	10	80	68	15	49	34
17	20	14	2	45	1	11	288	288	16	112	117
18	23	0	3	43	2	12	475	467	17	98	69
19	58	28	4	63	10	13	361	366	18	141	157
20	23	5	5	43	2	14	279	272	19	94	62
0.12.L			6	43	11	15	56	43	20	52	80
0	29	15	7	43	3	16	335	318	1.4.L		
1	29	4	8	43	1	17	341	323	0	185	175
2	78	10	9	40	0	18	181	195	1	125	116
3	27	19	10	52	4	19	332	358	2	228	219
4	27	30	11	40	7	20	274	274	3	62	50
5	25	5	13	40	0				4	120	122
			14	40	0						

1.4.L			10	145	144	17	65	48	0	54	12
5	154	134	11	161	158	18	67	64	1	23	16
6	315	314	12	91	87	19	58	65	2	23	35
7	76	80	13	56	91	20	31	24	3	72	22
8	250	262	14	85	84	1.9.L			4	40	37
9	62	66	15	216	209	0	18	26	5	56	22
10	118	103	16	143	130	1	25	26	6	23	11
11	78	53	17	76	73	2	49	27	7	23	13
12	125	136	18	31	38	3	18	28	8	43	17
13	47	51	19	115	105	4	49	41	9	23	14
14	45	28	20	27	31	5	49	30	10	34	31
15	123	131	1.7.L			6	52	33	11	31	17
16	65	75	0	60	55	7	27	37	12	20	24
17	35	84	1	85	73	8	78	67	13	23	7
18	141	145	2	94	94	9	63	51	14	45	28
19	129	121	3	80	52	10	112	99	15	20	24
20	181	162	4	67	72	11	85	72	16	20	42
1.5.L			5	129	122	12	56	37	17	23	36
0	248	246	6	125	125	13	29	15	18	49	11
1	107	79	7	63	71	14	58	29	19	20	17
2	216	214	8	89	92	15	54	35	20	23	15
3	56	49	9	56	51	16	78	65	1.12.L		
4	127	133	10	87	70	17	72	48	0	29	12
5	143	145	11	71	82	18	87	88	1	29	3
6	281	272	12	83	76	19	23	39	2	29	15
7	89	81	13	78	84	20	69	57	3	29	19
8	223	219	14	18	20	1.10.L			4	29	8
9	252	249	15	52	28	0	58	15	5	29	17
10	56	50	16	49	28	1	36	36	6	60	16
11	89	83	17	43	23	2	56	31	7	27	33
12	190	191	18	89	99	3	20	33	8	40	27
13	89	84	19	49	43	4	20	12	9	27	17
14	89	51	20	63	64	5	56	41	10	45	4
15	89	101	1.8.L			6	49	29	11	27	1
16	132	158	0	18	29	7	20	26	12	63	15
17	132	139	1	103	101	8	60	36	13	40	9
18	62	72	2	67	57	9	45	52	14	47	11
19	127	115	3	105	56	10	20	32	15	27	13
20	16	49	4	58	64	11	23	45	16	25	10
1.6.L			5	36	44	12	40	19	17	25	15
0	114	102	6	38	11	13	20	25	18	38	11
1	123	113	7	18	3	14	20	3	19	25	11
2	125	146	8	18	35	15	47	33	20	25	24
3	127	124	9	112	107	16	20	31	1.13.L		
4	45	65	10	61	63	17	67	33	0	38	26
5	134	132	11	109	59	18	31	10	1	38	20
6	47	45	12	40	45	19	60	51	2	38	19
7	170	167	13	49	51	20	58	31	3	94	12
8	112	108	14	63	52	1.11.L			4	38	1
9	96	83	15	63	51	5	36	9	5	36	9
			16	67	67						



2.7.L			14	80	67	2.12.L			6	45	11
			15	21	28				8	45	5
9	98	92	16	20	31	0	29	12	9	45	6
10	69	66	17	38	11	1	52	23	10	45	7
11	58	60	18	65	64	2	29	7	11	45	5
12	67	72	19	20	20	3	72	19	12	43	4
13	60	42	20	34	37	4	25	11	13	43	12
14	58	57				5	29	11	15	43	7
15	107	100	2.10.L			6	29	16	16	43	3
16	58	74				7	54	6	18	40	5
17	85	77	0	20	19	8	29	30			
18	18	16	1	72	37	9	29	22	3.0.L		
19	18	31	2	52	43	10	27	9			
20	54	39	3	45	47	11	27	18	0	303	307
2.8.L			4	20	30	12	27	21	1	29	11
			5	54	26	13	27	12	2	415	394
0	13	18	6	20	10	14	27	4	3	43	15
1	56	35	7	38	23	15	58	9	4	11	18
2	31	18	8	20	16	16	27	12	5	47	57
3	65	57	9	54	51	17	27	19	6	366	377
4	31	17	10	20	8	18	43	20	7	395	384
5	103	58	11	20	16	19	25	8	8	495	501
6	65	53	12	20	21	20	25	21	9	292	306
7	58	48	13	56	51				10	167	177
8	45	45	14	65	34	2.13.L			11	270	264
9	18	34	15	20	39				12	43	46
10	76	54	16	20	18	0	38	9	13	502	475
11	38	51	17	36	25	1	85	18	14	453	452
12	45	59	18	52	25	2	38	17	15	192	194
13	127	115	19	20	31	3	38	9	16	265	273
14	45	43	20	38	24	4	38	16	17	120	115
15	60	70				5	38	17	18	76	67
16	60	45	2.11.L			6	38	22	19	29	7
17	31	14				7	38	11	20	47	48
18	40	50	0	56	47	8	72	10			
19	09	68	1	23	21	9	83	10	3.1.L		
20	80	77	2	47	14	10	36	14			
2.9.L			3	23	5	11	36	6	0	299	308
			4	23	14	12	36	4	1	170	153
0	18	21	5	23	18	13	36	8	2	154	170
1	58	59	6	23	38	14	67	5	3	225	234
2	36	47	7	36	4	15	34	10	4	172	162
3	18	18	8	23	38	16	34	12	5	308	305
4	58	48	9	69	11	17	34	3	6	254	238
5	52	61	10	23	32	18	34	13	7	370	368
6	83	82	11	23	17	19	65	9	8	49	54
7	31	8	12	47	25	20	31	5	9	134	143
8	69	54	13	23	9				10	181	165
9	18	23	14	23	11	2.14.L			11	62	33
10	58	42	15	20	9				12	288	280
11	09	48	16	23	23	0	45	7	13	326	326
12	71	50	17	23	27	2	47	3	14	196	197
13	20	9	18	47	28	3	45	11	15	199	212
			19	20	10	4	98	10			
			20	23	47	5	45	3			

3.1.L			3.4.L			3.7.L			3.10.L		
16	270	270	0	161	160	5	185	182	12	34	14
17	239	233	1	257	252	6	125	117	13	18	32
18	71	63	2	178	175	7	52	46	14	18	28
19	112	125	3	56	38	8	165	159	15	83	63
20	223	223	4	71	60	9	16	35	16	18	13
			5	83	67	10	107	105	17	18	22
			6	265	267	11	78	63	18	18	33
			7	143	140	12	56	53	19	18	11
			8	196	201	13	65	54	20	58	61
			9	129	143	14	76	70			
			10	94	96	15	65	66	3.9.L		
			11	109	122	16	38	37	0	65	33
			12	214	224	17	18	24	1	65	67
			13	330	328	18	120	114	2	52	46
			14	263	273	19	18	16	3	69	42
			15	207	200	20	67	61	4	43	23
			16	165	163				5	56	48
			17	221	222	3.7.L			6	18	13
			18	152	158	0	38	33	7	63	33
			19	120	110	1	127	132	8	20	19
			20	117	129	2	67	59	9	38	41
						3	18	17	10	18	23
						4	43	33	11	58	47
						5	18	43	12	34	29
						6	109	97	13	40	34
						7	85	78	14	18	11
						8	49	34	15	20	26
						9	74	77	16	20	16
						10	52	24	17	63	45
						11	71	63	18	20	32
						12	18	27	19	60	45
						13	76	66	20	20	16
						14	18	22			
						15	47	32	3.10.L		
						16	63	52	0	54	41
						17	98	117	1	20	4
						18	52	16	2	87	70
						19	31	33	3	20	6
						20	91	85	4	67	32
									5	40	27
						3.8.L			6	20	30
						0	34	33	7	52	28
						1	34	40	8	45	15
						2	89	62	9	20	19
						3	114	98	10	31	45
						4	31	27	11	20	13
						5	65	56	12	45	31
						6	71	42	13	56	20
						7	18	27	14	20	9
						8	74	71	15	47	21
						9	18	35	16	40	33
						10	36	59			
						11	40	25			

3.10.L			3.13.L			12	51	71	19	54	30
17	20	7	0	79	10	13	58	42	20	31	36
18	20	39	1	38	5	14	154	162	4.3.L		
19	20	8	2	38	7	15	250	246			
20	20	27	3	38	14	16	192	213	0	221	218
3.11.L			4	38	6	17	263	257	1	154	173
			5	60	18	18	45	3	2	236	223
0	23	34	6	38	9	19	80	59	3	337	354
1	23	4	7	38	15	20	16	41	4	34	6
2	31	17	8	38	7	4.1.L			5	310	315
3	23	31	9	36	10				6	65	68
4	23	35	10	39	7	0	109	89	7	239	233
5	23	23	11	36	4	1	145	167	8	65	28
6	23	15	12	36	6	2	522	524	9	274	282
7	56	25	13	36	7	3	263	258	10	123	127
8	40	23	14	63	11	4	448	435	11	181	161
9	23	40	15	36	8	5	301	305	12	116	119
10	23	6	16	36	5	6	62	44	13	116	115
11	45	28	17	36	11	7	230	224	14	174	167
12	23	14	18	34	4	8	219	218	15	181	177
13	38	15	19	34	11	9	74	82	16	172	176
14	58	22	20	34	5	10	165	151	17	138	136
15	23	19	3.14.L			11	170	152	18	16	24
16	36	16				12	205	219	19	234	224
17	54	16	1	45	7	13	207	215	20	56	18
18	43	17	2	79	8	14	165	183	4.4.L		
19	23	18	3	89	7	15	107	93			
20	23	24	4	45	11	16	80	76	0	132	124
3.12.L			6	45	2	17	149	140	1	158	149
			7	45	3	18	78	86	2	203	203
0	29	17	8	45	8	19	156	190	3	181	179
1	29	4	9	85	6	20	205	219	4	190	178
2	67	13	10	45	6	4.2.L			5	187	186
3	29	11	11	74	1				6	170	170
4	40	19	13	43	4	0	158	182	7	125	123
5	72	20	14	43	6	1	261	257	8	145	144
6	29	16	15	43	1	2	448	463	9	161	165
7	29	11	16	43	6	3	216	197	10	118	119
8	63	9	4.0.L			4	236	238	11	98	95
9	47	12				5	163	144	12	381	376
10	29	16	0	263	256	6	156	183	13	312	300
11	27	4	1	460	461	7	393	406	14	154	136
12	65	20	2	214	208	8	187	193	15	143	121
13	27	13	3	125	117	9	451	466	16	149	150
14	67	14	4	58	74	10	136	143	17	71	62
15	27	6	5	29	40	11	212	198	18	47	18
16	27	17	6	129	107	12	381	391	19	69	71
17	27	15	7	243	247	13	141	143	20	18	7
18	25	14	8	386	376	14	118	110	4.5.L		
19	25	4	9	274	257	15	96	109			
20	25	3	10	277	279	16	203	201	0	52	47
			11	96	78	17	156	149			
						18	203	224			

4.5.L			6	105	104	13	20	26	20	23	10
1	172	187	7	54	88	14	20	1			
2	118	124	8	40	22	15	20	40	4.12.L		
3	221	235	9	85	71	16	29	20	0	58	8
4	80	70	10	18	23	17	20	12	1	29	4
5	185	177	11	56	42	18	67	40	2	29	17
6	174	177	12	56	64	19	52	29	3	29	4
7	243	243	13	18	5	20	36	20	4	29	14
8	134	138	14	56	40				5	76	28
9	60	64	15	58	27	4.10.L			6	29	17
10	37	95	16	18	22	0	20	28	7	29	19
11	116	111	17	65	63	1	47	21	8	29	10
12	100	117	18	45	7	2	56	38	9	29	12
13	56	48	19	89	65	3	63	29	10	29	9
14	76	86	20	76	75	4	20	8	11	60	10
15	36	45		4.8.L			5	36	12	29	13
16	65	56				6	52	47	13	29	15
17	94	54	0	100	94	7	36	47	14	27	9
18	16	28	1	18	17	8	20	38	15	27	13
19	136	129	2	49	49	9	38	24	16	65	4
20	36	46	3	69	82	10	20	5	17	27	7
			4	18	12	11	49	3	18	43	8
4.6.L			5	74	54	12	45	30	19	52	5
0	60	75	6	20	29	13	20	15	20	27	12
1	109	105	7	36	31	14	20	6			
2	105	103	8	94	70	15	38	24	4.13.L		
3	40	34	9	36	54	16	63	38	0	38	14
4	78	83	10	56	49	17	38	21	2	38	18
5	67	59	11	29	37	18	23	25	3	38	9
6	89	86	12	36	30	19	49	4	4	65	6
7	56	45	13	20	31	20	36	25	5	38	8
8	78	69	14	20	26				6	38	6
9	91	82	15	27	39	4.11.L			7	38	6
10	85	72	16	43	54	0	47	12	8	38	6
11	34	47	17	43	39	1	23	18	10	38	4
12	96	67	18	43	38	2	23	24	11	38	11
13	112	106	19	20	8	3	40	7	12	36	6
14	76	63	20	20	54	4	47	41	13	36	8
15	65	64		4.9.L			5	60	14	36	10
16	109	115				6	23	2	15	36	6
17	25	14	0	67	48	7	40	17	16	36	2
18	52	42	1	74	77	8	23	14	17	36	11
19	34	40	2	69	37	9	78	9	18	56	2
20	52	52	3	18	16	10	23	11	19	34	7
			4	20	16	11	23	15	20	34	11
4.7.L			5	63	57	12	23	20			
0	34	3	6	20	19	13	56	24	4.14.L		
1	43	37	7	49	43	14	56	38	0	98	5
2	60	53	8	20	35	15	23	24	1	47	3
3	38	40	9	80	70	16	47	30	2	45	5
4	49	18	10	47	34	17	23	8	3	47	4
5	18	21	11	20	34	18	23	14			
			12	20	14	19	23	9			

4.14.L			16	129	129	5.4.L			5	16	25
			17	212	217				6	16	6
4	45	10	18	167	159	0	49	60	7	127	123
5	47	4	19	16	41	1	125	136	8	69	59
6	67	4	20	47	46	2	78	73	9	49	29
7	45	6				3	245	240	10	40	27
8	45	14	5.2.L			4	78	68	11	100	101
9	45	7				5	49	51	12	76	66
10	45	7	0	152	137	6	87	87	13	43	65
11	45	6	1	196	202	7	143	131	14	45	22
12	72	3	2	71	65	8	149	149	15	83	108
13	45	6	3	199	210	9	185	176	16	69	82
5.0.L			4	127	117	10	69	86	17	89	75
0	14	39	5	96	84	11	323	322	18	45	21
1	154	155	6	134	131	12	114	110	19	34	43
2	161	182	7	216	242	13	31	35	20	47	41
3	156	167	8	96	101	14	85	100	5.7.L		
4	103	106	9	205	207	15	60	60	0	18	29
5	62	74	10	118	118	16	83	77	1	49	40
6	286	294	11	118	106	17	165	157	2	71	70
7	165	162	12	71	61	18	56	61	3	56	53
8	89	69	13	203	215	19	103	97	4	118	116
9	118	129	14	67	76	20	100	93	5	80	65
10	60	69	15	315	338	5.5.L			6	43	36
11	185	194	16	103	94	0	96	94	7	18	16
12	47	48	17	143	144	1	91	81	8	76	69
13	181	180	18	83	74	2	52	32	9	76	59
14	203	208	19	245	241	3	45	9	10	89	72
15	40	33	20	221	230	4	103	113	11	18	10
16	103	107	5.3.L			5	116	113	12	71	60
17	62	61	0	183	173	6	16	32	13	43	37
18	261	255	1	132	117	7	62	60	14	89	89
19	393	392	2	60	49	8	143	145	15	40	32
20	232	237	3	112	100	9	83	97	16	54	23
5.1.L			4	156	159	10	60	74	17	18	16
0	120	133	5	118	126	11	16	13	18	18	32
1	76	79	6	167	166	12	16	26	19	38	44
2	116	126	7	67	72	13	38	40	20	20	47
3	109	110	8	277	282	14	103	107	5.8.L		
4	127	136	9	252	257	15	58	57	0	31	18
5	230	238	10	176	185	16	65	64	1	18	29
6	138	149	11	85	83	17	36	14	2	43	22
7	80	79	12	214	214	18	138	136	3	80	75
8	187	196	13	274	281	19	52	59	4	60	63
9	85	66	14	138	147	20	67	47	5	105	107
10	167	187	15	103	113	5.6.L			6	20	5
11	69	61	16	120	122	0	34	16	7	87	76
12	344	344	17	83	81	1	58	42	8	52	51
13	165	158	18	299	295	2	65	49	9	18	18
14	145	146	19	76	58	3	16	21	10	52	27
15	56	40	20	91	98	4	47	45	11	69	49



5.8.L			17	31	20	0	38	3	20	23	31
12	29	26	18	45	30	1	85	1			
13	60	54	19	23	22	2	38	19	6.1.L		
14	60	19	20	54	22	3	38	7			
15	20	32				4	38	12	0	127	130
16	20	4	5.11.L			5	58	0	1	145	148
17	47	38	0	25	30	6	38	18	2	268	275
18	20	23	1	23	20	7	38	10	3	62	51
19	40	49	2	43	11	8	38	12	4	40	50
20	58	45	3	43	17	9	38	5	5	74	81
			4	23	3	10	38	4	6	27	25
5.9.L			5	23	19	11	81	5	7	67	75
0	20	26	6	52	19	12	38	2	8	161	168
1	36	25	7	36	14	13	36	6	9	80	65
2	20	27	8	60	23	14	60	3	10	123	106
3	20	34	9	43	26	15	36	3	11	152	152
4	31	55	10	56	25	16	36	10	12	234	225
5	49	49	11	23	7	17	36	4	13	190	181
6	67	55	12	47	24	18	36	5	14	286	296
7	45	42	13	43	24	19	65	7	15	167	171
8	29	22	14	23	18	20	34	7	16	138	129
9	20	32	15	23	7	5.14.L			17	100	107
10	20	13	16	23	15				18	38	18
11	38	11	17	23	10	0	45	1	19	47	72
12	45	54	18	23	23	1	45	9	20	196	193
13	31	23	19	23	10	2	45	4	6.2.L		
14	56	59	20	23	18	3	45	2			
15	20	11				4	45	6	0	60	52
16	20	29	5.12.L			5	65	4	1	145	156
17	18	20	0	29	15	6	45	7	2	54	32
18	47	28	1	29	18	7	45	9	3	74	46
19	34	14	2	29	25	6.0.L			4	230	240
20	20	20	3	29	15				5	62	62
			4	81	10	0	350	357	6	103	108
5.10.L			5	83	8	1	49	7	7	112	115
0	20	12	6	29	11	2	85	90	8	120	116
1	20	18	7	29	6	3	69	52	9	69	68
2	20	4	8	29	4	4	277	281	10	167	172
3	20	42	9	29	12	5	299	312	11	129	138
4	20	13	10	29	8	6	214	241	12	100	82
5	20	19	11	29	16	7	207	217	13	125	134
6	20	32	12	65	5	8	125	116	14	56	55
7	52	18	13	29	9	9	100	96	15	158	159
8	20	19	14	29	14	10	91	81	16	71	49
9	23	30	15	29	14	11	40	13	17	65	52
10	43	12	16	72	4	12	98	101	18	216	200
11	31	8	17	27	6	13	248	236	19	112	94
12	76	45	18	56	10	14	129	143	20	78	81
13	43	6	19	69	9	15	127	126	6.3.L		
14	23	1	20	27	9	16	181	177			
15	52	38	5.13.L			17	31	19	0	74	65
16	69	54				18	96	106	1	156	168
						19	29	30			

6.3.L			7	107	105	14	18	21	6.10.L		
2	91	103	8	25	34	15	60	51	0	58	30
3	214	206	9	47	17	16	31	46	1	20	24
4	80	60	10	47	41	17	54	22	2	20	11
5	38	40	11	71	72	18	58	24	3	20	32
6	54	52	12	18	37	19	18	16	4	36	18
7	129	124	13	25	45	20	20	26	5	31	17
8	100	107	14	76	71	6.8.L			6	20	23
9	149	162	15	87	80	0	29	21	7	20	13
10	78	59	16	18	16	1	56	64	8	20	23
11	239	245	17	58	55	2	34	18	9	20	10
12	112	130	18	49	52	3	20	21	10	40	26
13	132	130	19	71	40	4	118	113	11	23	17
14	112	113	20	38	16	5	40	30	12	69	45
15	165	190	6.6.L			6	20	39	13	20	14
16	145	142	0	54	53	7	56	51	14	20	3
17	100	93	1	67	60	8	74	33	15	45	13
18	89	70	2	18	2	9	20	36	16	20	29
19	67	78	3	120	105	10	76	42	17	34	7
20	58	64	4	118	93	11	54	36	18	47	43
6.4.L			5	38	17	12	31	28	19	45	18
0	277	294	6	36	22	13	63	31	20	45	16
1	58	62	7	52	36	14	60	32	6.11.L		
2	85	73	8	56	61	15	20	15	0	23	1
3	67	62	9	47	25	16	58	33	1	38	3
4	112	110	10	105	91	17	31	20	2	47	16
5	69	80	11	112	108	18	49	28	3	23	24
6	125	119	12	38	37	19	20	17	4	23	12
7	105	103	13	100	93	20	40	16	5	49	23
8	120	120	14	52	40	6.9.L			6	23	5
9	89	88	15	125	112	0	45	27	7	23	2
10	16	7	16	49	40	1	20	45	8	23	9
11	67	84	17	56	40	2	38	32	9	25	27
12	67	68	18	76	72	3	69	65	10	23	13
13	49	52	19	34	36	4	49	20	11	38	19
14	174	172	20	43	19	5	20	32	12	23	3
15	63	74	6.7.L			6	47	34	13	43	14
16	78	82	0	18	32	7	38	23	14	23	14
17	105	84	1	45	58	8	47	42	15	54	12
18	18	38	2	63	18	9	69	36	16	47	17
19	18	5	3	31	39	10	20	22	17	38	17
20	18	38	4	63	56	11	65	65	18	23	9
6.5.L			5	63	72	12	38	12	19	23	9
0	40	50	6	100	89	13	20	21	20	23	8
1	45	51	7	85	66	14	20	9	6.12.L		
2	136	134	8	96	85	15	20	13	0	94	19
3	58	32	9	65	52	16	34	16	1	29	10
4	25	15	10	34	42	17	56	26	2	31	23
5	47	36	11	58	64	18	20	23			
6	89	75	12	20	28	19	23	35			
			13	98	81	20	31	43			

6.12.L			8	259	256	15	91	91	7.5.L		
3	29	1	9	40	1	16	116	116	0	94	75
4	07	3	10	134	134	17	18	20	1	69	75
5	29	12	11	27	8	18	18	26	2	18	23
6	29	9	12	145	143	19	65	53	3	69	62
7	29	6	13	181	181	20	85	88	4	65	58
8	29	10	14	103	88	7.3.L			5	105	95
9	29	10	15	134	134	0	134	148	6	63	56
10	29	6	16	152	138	1	109	96	7	69	57
11	29	5	17	65	33	2	76	63	8	18	23
12	29	8	18	178	151	3	78	88	9	67	63
13	27	3	19	152	139	4	78	64	10	80	75
14	67	11	20	152	130	5	145	170	11	78	73
15	29	3	7.1.L			6	98	94	12	18	23
16	27	8	0	317	320	7	116	57	13	143	129
17	27	10	1	167	164	8	118	137	14	85	88
18	29	6	2	279	280	9	49	18	15	43	28
19	27	4	3	123	131	10	152	154	16	34	17
20	78	4	4	54	32	11	52	59	17	18	41
6.13.L			5	125	119	12	114	100	18	74	61
0	36	1	6	199	213	13	199	188	19	18	26
1	38	10	7	103	100	14	105	85	20	63	45
2	38	10	8	105	98	15	67	71	7.6.L		
3	38	5	9	377	393	16	18	28	0	38	44
4	38	8	10	67	63	17	18	43	1	56	48
5	38	5	11	83	57	18	18	13	2	18	38
6	38	4	12	85	80	19	67	48	3	54	29
7	38	15	13	141	135	20	58	30	4	87	81
8	38	4	14	161	160	7.4.L			5	31	50
9	38	13	15	109	112	0	27	36	6	114	107
10	36	2	16	78	72	1	56	68	7	83	75
11	36	2	17	141	143	2	58	85	8	80	40
12	36	10	18	114	102	3	149	145	9	18	18
13	36	6	19	129	110	4	85	91	10	89	59
14	36	5	20	74	60	5	67	65	11	109	110
15	36	2	7.2.L			6	34	41	12	60	50
16	36	2	0	87	100	7	18	34	13	18	19
17	36	10	1	114	122	8	89	87	14	74	68
18	36	9	2	167	161	9	71	57	15	18	47
19	52	15	3	96	87	10	116	120	16	65	70
20	34	3	4	223	236	11	107	68	17	18	20
7.0.L			5	225	239	12	83	66	18	58	24
0	292	299	6	183	191	13	60	57	19	43	37
1	71	71	7	252	267	14	74	75	20	89	88
2	178	173	8	145	154	15	58	59	7.7.L		
3	158	151	9	127	127	16	105	110	0	31	29
4	303	325	10	118	115	17	38	17	1	98	79
5	216	221	11	152	144	18	54	43	2	83	61
6	137	158	12	185	181	19	54	37	3	80	79
7	178	182	13	107	84	20	63	63			
			14	107	55						

7.7.L			9	49	25	16	52	4	3	143	145
			10	63	15	17	23	7	4	34	15
4	43	26	11	65	41	18	45	5	5	152	151
5	40	16	12	58	15	19	23	20	6	89	68
6	38	51	13	38	40	20	23	9	7	74	52
7	56	46	14	76	28				8	54	42
8	40	11	15	47	33	7.12.L			9	120	122
9	125	111	16	20	14				10	83	69
10	74	39	17	20	10	0	29	4	11	58	48
11	40	39	18	20	23	1	29	13	12	76	76
12	20	21	19	20	10	2	29	14	13	31	35
13	31	28	20	52	30	3	29	6	14	63	65
14	38	15				4	29	10	15	69	65
15	36	37	7.10.L			5	54	12	16	45	38
16	20	28				6	58	6	17	109	105
17	100	70	0	20	26	7	29	9	18	38	40
18	43	38	1	20	14	8	45	8	19	29	7
19	20	7	2	20	31	9	29	7	20	18	30
20	20	25	3	20	10	10	74	9			
7.8.L			4	23	37	11	29	14	8.1.L		
			5	49	22	12	29	12			
0	91	89	6	23	20	13	29	0	0	38	15
1	20	16	7	20	24	14	29	6	1	40	13
2	43	51	8	20	13	15	29	7	2	125	128
3	20	35	9	20	13	16	29	16	3	65	40
4	31	52	10	45	1	17	25	9	4	87	88
5	65	42	11	36	5	18	29	4	5	91	87
6	76	66	12	69	30	19	27	2	6	181	163
7	69	56	13	23	6	20	27	3	7	105	97
8	63	57	14	23	19				8	45	29
9	20	13	15	20	15	7.13.L			9	67	53
10	20	21	16	36	20				10	129	120
11	38	11	17	20	6	0	38	8	11	60	35
12	47	26	18	23	14	1	38	2	12	96	79
13	65	34	19	23	11	2	38	4	13	40	44
14	40	25	20	23	8	3	38	9	14	43	57
15	58	24				4	38	12	15	36	30
16	63	38	7.11.L			5	38	12	16	80	77
17	20	18				6	38	1	17	54	40
18	38	32	0	23	7	7	38	5	18	34	24
19	34	25	1	25	18	8	38	5	19	67	72
20	67	44	2	45	9	9	38	8	20	49	20
7.9.L			3	63	5	10	38	1			
			4	49	5	11	38	6	8.2.L		
0	20	52	5	43	24	12	38	4			
1	34	16	6	23	2	13	36	9	0	85	80
2	40	37	7	25	16	14	36	7	1	71	81
3	49	27	8	25	2	15	36	10	2	18	37
4	20	13	9	25	11	16	36	5	3	67	69
5	56	59	10	25	7				4	38	45
6	20	22	11	36	4	8.0.L			5	114	111
7	20	19	12	60	9				6	69	89
8	66	20	13	23	19	0	16	29	7	63	54
			14	49	16	1	69	53	8	80	84
			15	23	3	2	47	46			

8.2.L			14	18	34	8.7.L			5	20	27
			15	18	6				6	20	27
9	62	64	16	34	22	0	69	75	7	20	16
10	52	45	17	76	60	1	43	30	8	20	30
11	36	37	18	54	58	2	89	73	9	20	24
12	18	14	19	54	36	3	18	25	10	31	42
13	63	75	20	49	15	4	36	45	11	31	28
14	69	70	8.5.L			5	31	28	12	23	17
15	67	60				6	20	29	13	20	17
16	27	9				7	54	33	14	20	7
17	18	39	0	65	68	8	34	34	15	49	34
18	18	17	1	52	30	9	69	30	16	23	21
19	18	16	2	49	58	10	98	70	17	20	10
20	49	18	3	67	64	11	52	24	18	20	18
8.3.L			4	67	63	12	67	43	19	54	11
			5	43	34	13	31	11	20	36	18
0	40	53	6	49	28	14	63	53	8.10.L		
1	100	56	7	36	34	15	20	28			
2	83	83	8	27	19	16	20	38	0	34	4
3	31	24	9	69	78	17	20	1	1	49	20
4	74	78	10	69	63	18	49	32	2	47	26
5	18	23	11	78	82	19	36	20	3	20	10
6	18	35	12	45	40	20	43	49	4	23	13
7	18	23	13	34	26	8.8.L			5	20	10
8	18	19	14	100	96				6	23	9
9	31	46	15	60	13	0	20	13	7	20	18
10	116	105	16	38	28	1	47	34	8	20	9
11	52	44	17	20	36	2	20	40	9	23	12
12	45	40	18	63	24	3	20	28	10	45	2
13	69	59	19	40	44	4	38	26	11	20	16
14	18	27	20	54	33	5	38	20	12	23	16
15	96	79	8.6.L			6	25	17	13	43	16
16	29	35				7	20	30	14	23	11
17	18	9	0	83	80	8	56	21	15	23	13
18	29	39	1	49	21	9	76	63	16	38	24
19	56	26	2	45	15	10	20	7	17	63	9
20	63	55	3	118	100	11	56	37	18	20	13
8.4.L			4	27	21	12	49	29	19	52	20
			5	96	86	13	60	25	20	23	11
0	18	2	6	18	12	14	20	19	8.11.L		
1	34	38	7	83	87	15	23	36			
2	47	45	8	20	27	16	36	33	0	49	1
3	47	44	9	47	43	17	87	75	1	56	6
4	67	60	10	63	45	18	20	15	2	23	10
5	107	91	11	20	38	19	23	35	3	23	12
6	38	37	12	49	19	20	43	46	4	25	11
7	103	92	13	36	34	8.9.L			5	47	15
8	18	35	14	18	19				6	23	14
9	56	41	15	20	21	0	20	26	7	23	7
10	45	18	16	18	21	1	20	7	8	23	10
11	69	60	17	65	43	2	20	37	9	23	11
12	45	39	18	59	57	3	45	12			
13	18	20	19	52	9	4	20	15			
			20	20	21						

8.11.L			2	69	57	9	67	63	16	94	81
10	58	10	3	69	49	10	78	76	17	103	97
11	23	13	4	38	5	11	60	41	18	60	30
12	23	9	5	147	148	12	60	56	19	89	66
13	23	6	6	112	95	13	18	7	20	20	16
14	23	6	7	87	90	14	49	25			
15	23	7	8	103	111	15	40	42	9.5.L		
16	23	16	9	78	83	16	18	11	0	65	64
17	23	9	10	18	10	17	67	53	1	29	32
18	23	4	11	18	43	18	87	83	2	38	36
19	23	16	12	34	13	19	47	25	3	20	44
20	23	10	13	94	97	20	76	63	4	67	40
			14	87	73				5	18	19
			15	112	87	9.3.L			6	43	42
8.12.L			16	31	5	0	149	153	7	43	24
0	29	1	17	34	36	1	60	56	8	63	36
1	29	10	18	18	46	2	89	97	9	38	67
2	29	5	19	18	16	3	45	57	10	56	57
3	54	7	20	38	5	4	25	4	11	78	52
4	29	3				5	63	69	12	49	53
5	29	8	9.1.L			6	49	24	13	20	29
6	29	8	0	63	60	7	91	58	14	76	63
7	29	8	1	94	84	8	80	87	15	83	61
8	29	7	2	85	85	9	18	18	16	52	54
9	29	6	3	83	71	10	51	75	17	83	63
10	29	5	4	96	82	11	43	47	18	20	24
11	45	4	5	54	37	12	31	39	19	63	47
12	29	8	6	116	108	13	18	33	20	31	45
13	29	13	7	43	42	14	65	66			
14	76	4	8	65	41	15	47	33	9.6.L		
15	72	16	9	18	41	16	38	43	0	60	51
16	29	11	10	45	42	17	78	62	1	20	22
17	29	2	11	65	64	18	54	37	2	67	41
18	27	8	12	85	84	19	29	26	3	45	59
19	45	5	13	18	23	20	58	51	4	54	53
20	65	13	14	78	68				5	40	57
			15	31	40	9.4.L			6	20	23
8.13.L			16	69	63	0	52	47	7	45	39
			17	43	40	1	52	48	8	29	22
0	38	8	18	18	12	2	49	44	9	20	16
1	38	6	19	47	22	3	36	33	10	72	63
2	38	5	20	54	37	4	40	47	11	20	38
3	38	3				5	58	31	12	58	48
4	38	2	9.2.L			6	18	16	13	58	56
5	38	7	0	31	1	7	100	88	14	38	26
6	38	12	1	54	32	8	67	53	15	96	72
7	38	7	2	98	132	9	78	62	16	20	36
8	38	6	3	58	25	10	58	38	17	20	12
9	69	5	4	112	110	11	58	63	18	65	60
			5	98	91	12	72	45	19	40	21
9.J.L			6	56	36	13	67	47	20	45	15
0	149	159	7	45	54	14	18	29			
1	109	94	8	13	53	15	20	30			

9.7.L			5	36	24	12	23	4	20	36	22
0	20	3	6	20	19	13	25	9			
1	34	20	7	20	10	14	45	3	10.1.L		
2	42	13	8	43	16	15	23	7			
3	47	56	9	20	7	16	52	5	0	54	36
4	83	49	10	23	6	17	23	6	1	78	53
5	52	21	11	23	11	18	23	16	2	74	60
6	69	41	12	20	8	19	25	2	3	63	61
7	38	32	13	45	24	20	36	6	4	71	65
8	38	13	14	23	8				5	67	50
9	20	33	15	34	11	9.12.L			6	18	20
10	20	31	16	40	3	0	29	4	7	69	58
11	20	9	17	20	4	1	29	1	8	49	31
12	58	32	18	23	11	2	29	5	9	34	39
13	20	21	19	47	5	3	29	5	10	18	37
14	20	10	20	80	20	4	29	1	11	47	34
15	40	37	9.10.L			5	29	6	12	18	21
16	45	18				6	29	5	13	18	39
17	20	15	0	78	18	7	65	7	14	58	39
18	52	34	1	23	12	8	54	15	15	20	23
19	20	26	2	54	18	9	29	11	16	31	35
20	20	7	3	52	10	10	29	8	17	18	3
			4	23	20	11	56	9	18	18	29
9.8.L			5	20	7	12	29	3	19	63	52
			6	23	10	13	29	9	20	20	18
0	20	16	7	23	6	14	29	4	10.2.L		
1	40	25	8	49	3	15	29	3			
2	20	8	9	23	15	16	27	6	0	18	15
3	45	12	10	23	2	17	27	12	1	18	27
4	20	26	11	43	18	18	27	8	2	67	74
5	45	1	12	34	12	19	27	7	3	52	35
6	47	39	13	45	2				4	83	77
7	20	14	14	23	12	10.0.L			5	56	54
8	54	22	15	31	5				6	83	91
9	20	15	16	23	13	0	103	97	7	31	16
10	23	9	17	23	11	1	18	19	8	18	38
11	20	12	18	23	7	2	58	55	9	49	58
12	63	31	19	23	5	3	18	1	10	31	17
13	20	36	20	23	16	4	43	37	11	63	21
14	65	37				5	18	7	12	56	47
15	38	25	9.11.L			6	18	31	13	63	61
16	45	31				7	18	27	14	18	11
17	43	33	0	25	21	8	27	29	15	63	30
18	20	14	1	25	11	9	43	17	16	47	42
19	23	17	2	25	14	10	40	51	17	38	49
20	20	9	3	23	2	11	27	12	18	36	44
			4	25	6	12	18	44	19	43	20
9.9.L			5	25	11	13	18	13	20	20	31
			6	25	4	14	18	1			
0	23	5	7	23	2	15	40	22	10.3.L		
1	20	4	9	23	13	16	47	42			
2	20	18	9	56	12	17	18	2	0	56	22
3	54	22	10	25	9	18	38	17	1	80	79
4	20	24	11	34	8	19	18	23			

10.3.L			7	83	61	14	23	11	10.10.L		
2	54	39	8	63	4	15	38	23	0	23	14
3	18	44	9	63	44	16	23	19	1	23	5
4	40	49	10	20	12	17	20	16	2	23	23
5	34	27	11	103	90	18	23	10	3	31	9
6	43	13	12	67	5	19	63	15	4	23	7
7	58	48	13	80	73	20	45	11	5	23	3
8	34	45	14	58	9	10.9.L			6	23	16
9	18	17	15	85	82	0	20	1	7	23	15
10	58	43	16	20	8	1	20	24	8	23	5
11	47	44	17	56	46	2	20	19	9	23	2
12	18	12	18	45	17	3	20	5	10	23	12
13	56	45	19	20	35	4	20	21	11	23	12
14	18	7	20	40	22	5	45	16	12	23	13
15	58	65	10.6.L			6	56	16	13	69	4
16	47	22	0	63	15	7	20	9	14	23	7
17	38	31	1	38	20	8	56	12	15	58	6
18	31	38	2	20	27	9	43	20	16	23	11
19	76	62	3	20	15	10	20	18	17	23	4
20	72	58	4	20	28	11	47	16	18	23	7
10.4.L			5	20	26	12	45	6	19	23	4
0	20	15	6	20	26	13	80	21	20	23	9
1	18	30	7	78	60	14	43	7	10.11.L		
2	38	28	8	20	14	15	20	12	0	23	9
3	18	16	9	94	70	16	23	6	1	25	8
4	20	39	10	23	34	17	45	11	2	34	6
5	52	47	11	40	20	18	20	12	3	25	8
6	94	65	12	56	30	19	36	9	4	25	10
7	43	31	13	20	4	20	38	15	5	23	5
8	107	89	14	52	45	10.9.L			6	54	5
9	87	66	15	20	8	0	20	9	7	25	11
10	20	15	16	49	29	1	23	4	8	23	7
11	54	32	17	49	19	2	20	18	9	23	18
12	20	14	18	31	13	3	20	13	10	40	2
13	20	25	19	49	33	4	23	7	11	23	9
14	20	9	20	36	15	5	23	25	12	25	3
15	56	33	10.7.L			6	23	5	13	23	2
16	49	20	0	56	42	7	20	10	14	23	4
17	60	58	1	20	7	8	20	4	15	25	7
18	67	33	2	69	52	9	20	3	16	31	2
19	20	23	3	49	24	10	20	3	17	23	4
20	85	75	4	45	13	11	47	3	18	23	9
10.5.L			5	20	13	12	52	13	19	31	10
0	29	4	6	38	10	13	23	2	20	69	6
1	100	101	7	26	13	14	23	10	10.12.L		
2	54	14	8	47	14	15	43	4	0	29	11
3	20	45	9	34	15	16	23	6	1	29	2
4	18	4	10	34	14	17	23	5	2	29	7
5	20	7	11	20	23	18	23	6			
6	36	38	12	52	22	19	20	15			
			13	20	5	20	23	12			



1J.12.L			18	94	75	1	20	31	8	20	12
3	29	2	19	47	37	2	54	38	9	20	34
4	29	4	20	20	15	3	89	59	10	20	13
5	29	4	11.2.L			4	52	66	11	36	13
6	29	10				5	100	98	12	20	7
7	29	4	0	20	35	6	34	23	13	38	5
8	29	2	1	74	54	7	65	52	14	20	14
9	29	9	2	52	15	8	20	25	15	20	14
10	29	2	3	20	19	9	45	24	16	20	34
11.J.L			4	74	66	10	45	23	17	49	20
0	67	66	5	20	11	11	43	22	18	20	11
1	31	10	6	20	42	12	67	47	19	49	24
2	20	2	7	76	75	13	20	41	20	36	29
3	136	121	8	49	57	14	20	6	11.7.L		
4	96	88	9	63	47	15	20	1			
5	16	16	10	49	36	16	40	12	0	20	34
6	18	43	11	33	44	17	36	28	1	20	15
7	56	24	12	20	21	18	20	36	2	20	22
8	18	23	13	20	21	19	52	42	3	20	12
9	49	32	14	20	22	20	36	13	4	20	14
10	20	25	15	20	26	11.5.L			5	20	19
11	80	62	16	20	42	0	54	25	6	52	20
12	47	34	17	54	38	1	34	5	7	20	12
13	18	4	18	20	26	2	56	47	8	20	14
14	52	28	19	87	59	3	49	32	9	36	22
15	72	50	20	20	34	4	49	38	10	67	7
16	40	12	11.3.L			5	20	23	11	67	3
17	40	30	0	20	15	6	31	9	12	58	30
18	20	7	1	60	20	7	20	34	13	20	12
19	20	21	2	63	53	8	56	4	14	20	22
20	52	39	3	89	78	9	20	30	15	20	7
11.1.L			4	65	56	10	20	27	16	40	9
0	118	56	5	43	8	11	49	46	17	43	8
1	94	78	6	80	41	12	58	45	18	20	22
2	76	63	7	49	41	13	20	3	19	20	12
3	52	33	8	36	9	14	31	42	20	43	8
4	18	26	9	20	27	15	20	23	11.8.L		
5	71	52	10	79	69	16	20	31	0	54	12
6	71	50	11	20	34	17	20	1	1	40	10
7	18	12	12	49	37	18	20	22	2	23	9
8	60	52	13	54	18	19	20	15	3	23	18
9	60	49	14	65	53	20	20	18	4	23	10
10	20	21	15	49	12	11.6.L			5	20	6
11	63	52	16	54	31	0	34	27	6	20	4
12	54	26	17	69	32	1	20	14	7	20	6
13	20	38	18	43	29	2	47	13	8	20	11
14	43	26	19	20	9	3	56	26	9	23	19
15	34	11	20	20	7	4	20	6	10	20	9
16	20	20	11.4.L			5	36	25	11	23	18
17	56	52	0	20	2	6	20	5	12	23	12
						7	45	37	13	40	2

11.8.L			19	23	3	5	20	25	12	20	27
14	23	3	20	23	8	6	20	29	13	38	21
15	23	9	11.11.L			7	20	17	14	36	15
16	20	5				8	20	18	15	20	25
17	31	8	0	34	4	9	20	17	16	52	17
18	60	12	1	25	7	10	20	41	17	20	21
19	40	11	2	34	7	11	72	41	18	54	22
20	23	6	3	25	9	12	20	13	19	20	3
			4	25	14	13	36	25	20	20	17
			5	63	4	14	40	30			
11.9.L			6	25	10	15	31	27	12.4.L		
0	49	21	7	49	5	16	20	14	0	20	15
1	23	14	8	23	3	17	65	43	1	65	29
2	38	7	9	25	9	18	63	22	2	20	2
3	20	7	10	23	3	19	20	31	3	20	20
4	58	8	11	23	6	20	20	22	4	45	16
5	23	9	12	60	5	12.2.L			5	43	28
6	23	6	13	58	11				6	47	29
7	23	7	14	23	11	0	20	13	7	20	20
8	46	20	15	25	4	1	20	12	8	45	7
9	23	15	16	23	7	2	20	14	9	40	36
10	23	13	17	23	0	3	43	31	10	20	20
11	23	8				4	65	47	11	20	13
12	23	5	12.0.L			5	20	29	12	52	11
13	23	4				6	65	15	13	36	19
14	23	5	0	47	38	7	20	36	14	20	11
15	49	4	1	63	33	8	29	26	15	20	16
16	23	8	2	20	13	9	20	16	16	20	20
17	23	6	3	18	31	10	20	34	17	20	15
18	20	14	4	20	12	11	43	16	18	20	21
19	23	6	5	18	8	12	20	18	19	75	25
20	31	8	6	20	5	13	60	14	20	23	9
			7	56	45	14	20	10			
11.10.L			8	76	56	15	52	45	12.5.L		
0	56	0	9	40	13	16	47	26	0	20	18
1	23	15	10	29	24	17	34	34	1	20	16
2	23	8	11	69	42	18	20	32	2	67	27
3	23	3	12	20	15	19	29	24	3	20	29
4	23	7	13	20	17	20	20	5	4	20	6
5	23	4	14	36	41	12.3.L			5	43	9
6	23	13	15	20	22				6	29	21
7	23	16	16	31	22	0	31	34	7	20	13
8	52	8	17	20	30	1	20	22	8	40	23
9	23	5	18	34	9	2	52	34	9	31	23
10	23	8	19	69	31	3	54	23	10	20	8
11	23	10	20	36	30	4	20	18	11	20	14
12	23	3	12.1.L			5	47	15	12	38	22
13	45	9				6	20	32	13	20	18
14	23	2	0	38	40	7	20	19	14	20	10
15	23	12	1	20	18	8	63	33	15	20	9
16	38	12	2	58	42	9	45	21	16	56	16
17	23	4	3	20	29	10	31	25	17	20	4
18	47	11	4	20	27	11	20	24			

12.5.L			12.8.L			5	23	11	7	76	49
18	23	10	0	47	17	6	45	5	8	20	7
19	43	10	1	23	8	7	23	7	9	54	29
20	36	21	2	23	6	8	36	9	10	52	18
12.6.L			3	23	11	9	23	11	11	38	4
0	23	32	4	38	11	10	23	9	12	31	19
1	20	21	5	36	18	11	54	8	13	58	34
2	63	11	6	49	7	12	23	9	14	20	21
3	20	6	7	23	11	13	23	3	15	20	19
4	49	10	8	23	10	14	23	2	16	20	21
5	52	20	9	23	8	15	23	10	17	20	5
6	20	10	10	49	4	16	49	8	18	45	23
7	23	12	11	54	14	17	23	11	19	40	23
8	20	4	12	23	14	18	36	8	20	54	16
9	65	15	13	23	8	19	23	3	13.2.L		
10	20	16	14	23	11	20	23	7	0	56	55
11	20	15	15	23	8	12.11.L			1	20	20
12	23	15	16	23	9	0	47	10	2	20	28
13	49	13	17	38	8	1	38	8	3	20	21
14	20	4	18	23	4	13.0.L			4	49	22
15	20	19	19	23	8	0	45	7	5	49	23
16	23	16	20	38	5	1	43	5	6	23	15
17	34	11	12.9.L			2	34	16	7	36	14
18	34	24	0	23	2	3	40	17	8	36	38
19	23	25	1	31	2	4	43	54	9	20	17
20	23	6	2	72	20	5	69	22	10	43	28
12.7.L			3	52	6	6	63	19	11	60	20
0	23	1	4	23	14	7	78	60	12	43	23
1	23	7	5	23	10	8	23	18	13	20	8
2	34	18	6	54	11	9	20	19	14	23	16
3	38	10	7	23	8	10	49	1	15	49	22
4	45	13	8	40	3	11	43	24	16	20	25
5	20	3	9	23	2	12	20	20	17	20	15
6	20	10	10	67	7	13	23	3	18	31	27
7	20	16	11	23	2	14	60	60	19	54	26
8	56	9	12	23	11	15	23	8	20	34	16
9	40	17	13	23	2	16	38	35	13.3.L		
10	23	19	14	34	13	17	23	20	0	38	12
11	23	13	15	49	1	18	45	12	1	20	21
12	52	13	16	23	12	19	20	8	2	43	17
13	20	11	17	60	7	20	23	9	3	36	23
14	23	2	18	40	4	13.1.L			4	20	4
15	23	9	19	23	10	0	34	23	5	20	13
16	23	6	20	40	6	1	20	12	6	20	18
17	38	12	12.10.L			2	40	32	7	20	11
18	40	9	0	23	1	3	63	21	8	23	24
19	49	15	1	23	6	4	20	7	9	49	41
20	36	13	2	56	9	5	54	42	10	47	6
			3	52	13	6	56	44	11	58	11
			4	23	5				12	40	18
									13	23	25

13.3.L			19	20	12	2	43	10	9	23	5
14	34	8	20	23	6	3	23	8			
15	20	20	13.6.L			4	43	17	14.0.L		
16	38	15				5	40	4			
17	60	31	0	23	11	6	23	11	0	23	13
18	20	14	1	43	4	7	63	19	1	20	19
19	23	19	2	52	19	8	23	9	2	38	13
20	23	12	3	23	6	9	49	1	3	20	6
			4	23	14	10	52	5	4	20	9
13.4.L			5	36	11	11	34	14	5	40	25
0	20	18	6	23	11	12	45	9	6	20	21
1	40	16	7	20	13	13	56	5	7	38	13
2	63	23	8	23	27	14	58	7	8	23	42
3	20	4	9	23	4	15	47	4	9	23	1
4	20	2	10	23	14	16	80	3	10	23	19
5	23	6	11	23	15	17	23	4	11	47	35
6	20	30	12	31	3	18	23	4	12	23	18
7	54	17	13	43	5	19	58	6	13	23	23
8	23	20	14	38	4	20	23	9	14	47	18
9	23	16	15	23	15	13.9.L			15	49	11
10	34	6	16	23	19				16	23	10
11	23	21	17	23	6	0	54	15	17	56	13
12	36	26	18	23	7	1	36	11	18	69	52
13	20	8	19	23	14	2	23	15	19	63	24
14	52	19	20	23	17	3	38	7	20	23	25
15	43	23	13.7.L			4	23	7	14.1.L		
16	34	11				5	23	4			
17	23	4	0	23	2	6	23	7	0	20	10
18	23	7	1	23	12	7	23	2	1	23	27
19	54	18	2	38	5	8	23	4	2	23	22
20	23	29	3	23	3	9	23	11	3	45	27
			4	20	6	10	23	8	4	31	19
13.5.L			5	40	22	11	23	9	5	43	16
0	20	4	6	52	10	12	38	5	6	34	13
1	23	26	7	40	7	13	31	7	7	47	16
2	23	18	8	40	4	14	23	2	8	45	20
3	31	11	9	23	15	15	47	4	9	45	9
4	43	7	10	36	6	16	23	2	10	23	29
5	54	22	11	23	3	17	23	3	11	65	34
6	23	5	12	40	8	18	23	4	12	23	9
7	20	8	13	74	24	19	38	5	13	52	17
8	54	6	14	23	8	20	23	8	14	23	17
9	36	23	15	29	15	13.10.L			15	52	19
10	23	12	16	23	2				16	43	23
11	23	12	17	23	1	0	23	13	17	23	24
12	40	6	18	23	9	1	23	6	18	40	14
13	23	16	19	23	4	2	52	5	19	20	8
14	23	16	20	23	3	3	23	5	20	56	27
15	26	13	13.8.L			4	23	4	14.2.L		
16	23	3				5	23	1			
17	47	35	0	23	9	6	54	7	0	36	12
18	23	13	1	52	6	7	23	4	1	49	24
						8	34	8			

14.2.L			7	23	7	14	38	10	14.9.L		
2	34	4	8	23	8	15	23	9	0	54	6
3	22	20	9	23	14	16	23	6	1	23	2
4	23	33	10	20	18	17	29	3	2	74	13
5	23	11	11	22	15	18	23	3	3	20	8
6	23	24	12	69	11	19	47	3	4	56	4
7	65	19	13	52	16	20	23	5	5	23	4
8	20	24	14	20	3	14.7.L			6	38	5
9	43	24	15	23	16	0	23	5	7	54	4
10	23	8	16	36	19	1	67	9	8	23	6
11	20	8	17	23	9	2	23	6	9	23	4
12	23	0	18	23	6	3	23	12	10	67	3
13	49	13	19	23	3	4	20	13	11	23	9
14	20	15	20	49	2	5	23	4	15.0.L		
15	23	19	14.5.L			6	23	11	0	20	0
16	23	9	0	20	13	7	23	8	1	56	11
17	56	23	1	20	12	8	47	6	2	23	2
18	20	8	2	23	4	9	52	15	3	20	10
19	23	2	3	23	7	10	23	14	4	47	7
20	23	18	4	23	9	11	40	9	5	34	5
14.3.L			5	20	20	12	23	8	6	23	7
0	45	8	6	31	20	13	47	4	7	60	19
1	23	19	7	56	19	14	31	6	8	20	13
2	20	5	8	23	19	15	23	3	9	23	15
3	20	5	9	23	3	16	23	15	10	20	3
4	23	28	10	23	10	17	43	4	11	40	6
5	36	20	11	23	2	18	23	5	12	23	3
6	34	31	12	23	8	19	23	6	13	52	2
7	23	11	13	43	12	20	23	14	14	36	2
8	20	15	14	36	4	14.8.L			15	23	15
9	20	10	15	23	7	0	23	12	16	31	0
10	23	15	16	36	7	1	23	5	17	49	16
11	20	6	17	23	5	2	23	10	18	52	5
12	31	18	18	38	14	3	23	5	19	23	14
13	63	19	19	49	9	4	63	5	20	23	11
14	23	14	20	60	10	5	23	7	15.1.L		
15	34	10	14.6.L			6	63	9	0	23	20
16	23	11	0	23	6	7	45	4	1	58	3
17	23	6	1	34	22	8	38	9	2	43	9
18	23	23	2	56	11	9	23	3	3	54	7
19	23	11	3	40	16	10	23	5	4	54	10
20	23	9	4	23	9	11	49	9	5	47	8
14.4.L			5	23	8	12	23	2	6	23	15
0	43	2	6	23	11	13	23	8	7	52	6
1	20	9	7	29	16	14	23	2	8	31	12
2	23	18	8	58	5	15	23	4	9	23	4
3	23	31	9	45	16	16	23	4	10	40	14
4	52	11	10	23	3	17	54	7	11	36	5
5	45	11	11	23	9	18	63	13			
6	20	8	12	56	12	19	23	5			
			13	69	13	20	60	5			

15.1.L			19	23	13	4	23	8	13	23	3
12	20	11	20	23	4	5	40	7	14	36	9
13	23	4	15.3.L			6	23	5	15	34	7
14	40	12	0	20	5	7	31	10	16	20	8
15	36	9	1	60	8	8	23	1	17	40	6
16	23	10	2	23	14	9	23	12	18	45	8
17	20	5	3	23	2	10	23	5	19	23	2
18	34	7	4	36	17	11	23	9	20	31	7
19	23	8	5	23	7	12	23	8	15.6.L		
20	23	5	6	23	12	13	23	10	0	43	7
15.2.L			7	31	7	14	23	8	1	23	5
0	36	5	8	63	7	15	23	7	2	23	7
1	23	7	9	20	6	16	23	8	3	60	8
2	36	6	10	23	9	17	58	0	4	45	9
3	20	16	11	20	5	18	58	3	5	40	12
4	34	10	12	52	11	19	31	3	6	23	4
5	47	14	13	38	3	20	23	2	7	23	6
6	67	3	14	23	7	15.5.L			8	63	11
7	49	8	15	23	1	0	23	2	9	23	7
8	43	3	16	20	12	1	23	11	10	58	6
9	20	7	17	23	2	2	23	11	11	23	13
10	23	8	18	23	13	3	23	12	12	23	7
11	20	13	19	43	3	4	45	9	13	23	3
12	23	6	20	23	6	5	31	13	14	49	10
13	23	13	15.4.L			6	20	7	15	23	2
14	23	4	0	34	16	7	23	13	16	36	7
15	34	10	1	23	13	8	23	6	17	23	5
16	23	1	2	23	6	9	45	5			
17	23	3	3	23	9	10	58	4			
18	23	4	3	23	9	11	23	4			
						12	38	6			

APPENDIX B

X-RAY DATA OF MORPHINE METHYLIODIDE

0,0,L			7	260	233	0	301	322	14	138	130
2	218	219	8	464	423	1	665	631	15	34	36
4	1333	1216	9	821	776	2	240	220	16	46	60
6	629	572	10	73	82	3	1278	1147	17	150	152
8	875	851	11	437	420	4	73	49	0,9,L		
10	901	914	12	476	477	5	65	51	1	139	138
12	401	410	13	248	218	6	21	30	2	320	310
14	646	658	14	309	312	7	1000	943	3	146	132
16	263	262	15	405	418	8	316	284	4	340	306
0,1,L			16	428	434	9	348	312	5	296	272
1	481	497	17	22	38	10	229	210	6	331	297
3	68	35	0,4,L			11	559	546	7	58	20
4	1220	1093	0	869	870	12	110	97	8	491	464
5	875	786	1	343	328	13	355	365	9	248	243
6	1384	1254	2	442	373	14	207	189	10	207	212
7	507	462	3	1207	1084	15	90	84	11	150	150
8	1379	1301	4	830	695	16	41	8	12	372	371
9	845	800	5	294	269	17	253	260	13	22	37
10	252	250	6	379	352	0,7,L			14	185	172
11	38	67	7	1071	998	1	576	581	15	95	115
12	979	970	8	646	602	2	782	712	16	399	390
13	396	412	9	277	277	3	321	313	17	65	18
14	22	86	10	459	439	4	219	212	0,10,L		
15	160	167	11	585	548	5	597	565	0	624	675
16	491	497	12	75	80	6	204	182	1	126	125
17	24	11	13	561	540	7	32	44	2	21	17
0,2,L			14	211	208	8	342	325	3	21	38
0	68	107	15	265	246	9	586	560	4	270	261
3	1051	916	16	102	123	10	21	12	5	38	19
4	1243	1079	17	678	678	11	382	346	6	473	438
5	189	131	0,5,L			12	70	58	7	22	9
6	787	758	1	1508	1391	13	48	26	8	311	301
7	471	432	2	518	456	14	22	33	9	168	172
8	418	363	3	133	148	15	473	474	10	401	384
9	467	444	4	216	188	16	143	150	11	55	87
10	1101	1040	5	1462	1289	17	48	70	12	22	30
11	321	318	6	338	291	0,8,L			13	107	118
12	19	16	7	167	158	0	316	327	14	269	272
13	297	305	8	197	184	1	141	155	15	24	14
14	709	706	9	841	793	2	609	597	16	70	78
15	233	249	10	70	68	3	479	438	17	109	105
16	401	414	11	614	567	4	474	410	0,11,L		
17	144	163	12	53	24	5	60	25	1	124	128
0,3,L			13	476	458	6	335	308	2	439	427
3	212	197	14	235	226	7	592	560	3	21	31
4	360	320	15	277	274	8	265	247	4	263	248
5	473	448	16	85	75	9	195	177	5	22	6
6	602	516	17	41	59	10	398	380	6	381	359
			0,6,L			11	223	218			
						12	55	40			
						13	163	146			



0,11,L			1,1,L			14 277 266			6 721 667		
7	151	143	0	49	68	15	311	314	7	309	296
8	335	300	1	838	855	16	314	321	8	258	259
9	68	81	2	221	220	17	53	53	9	83	77
10	48	8	3	1149	1063	1,4,L			10	736	716
11	73	2	4	364	315	0	1707	1775	11	139	136
12	252	248	6	304	277	1	677	694	12	173	186
13	46	48	7	219	213	2	274	211	13	190	188
14	24	33	8	180	162	3	377	336	14	632	638
15	284	273	9	940	866	4	1501	1333	15	136	142
16	68	26	10	121	112	5	253	196	16	269	261
17	102	120	11	641	627	6	1171	1081	17	126	107
0,12,L			12	333	326	7	194	187	1,7,L		
0	272	305	13	357	358	8	469	432	0	82	9
1	122	119	14	102	118	9	308	294	1	185	211
2	92	80	15	580	602	10	894	852	2	1015	989
3	155	148	16	105	118	11	153	123	3	219	215
4	282	269	17	95	111	12	95	78	4	508	489
5	73	23	1,2,L			13	199	199	5	529	502
6	190	154	0	428	412	14	422	422	6	568	538
7	133	123	1	986	1025	15	105	96	7	117	126
8	55	57	2	306	274	16	214	213	8	624	590
9	148	143	4	1288	1147	17	151	161	9	428	423
10	352	348	5	530	465	1,5,L			10	122	116
11	170	166	6	874	791	0	253	240	11	206	197
12	87	18	7	1203	1109	1	275	254	12	546	534
13	88	70	8	605	553	2	1028	978	13	153	149
14	92	85	9	527	495	3	534	505	14	83	77
15	104	98	10	709	681	4	578	533	15	303	303
16	61	32	11	712	699	5	85	64	16	165	179
17	68	14	12	73	28	6	826	760	17	36	25
1,0,L			13	481	470	7	17	24	1,8,L		
1	585	637	14	440	431	8	1062	998	0	685	747
2	207	200	15	240	256	9	92	83	1	445	483
3	1319	1196	16	267	271	10	316	299	2	296	286
4	58	29	17	201	201	11	19	30	3	549	534
5	99	67	1,3,L			12	756	754	4	615	569
6	66	54	0	495	555	13	160	152	5	223	222
7	1776	1644	1	510	516	14	342	332	6	309	277
8	195	208	3	357	310	15	139	150	7	464	447
9	309	296	4	508	479	16	367	390	8	177	160
10	347	351	5	1090	972	17	109	96	9	280	266
11	847	819	6	1054	946	1,6,L			10	379	362
12	17	11	7	503	450	0	988	1054	11	394	385
13	677	662	8	957	881	1	447	463	12	105	128
14	199	204	9	886	830	2	399	383	13	382	380
15	104	97	10	233	219	3	355	355	14	277	286
16	95	107	11	231	241	4	994	945	15	155	156
17	466	485	12	666	657	5	180	139	16	218	215
			13	350	337				17	228	230

1,9,L			11 180 169			6 884 842			2,4,L		
0	156	205	12	76	77	7	269	251	0	122	160
1	549	599	13	131	141	8	1246	1210	1	563	637
2	357	350	14	155	136	9	280	257	2	445	475
3	428	421	15	321	322	10	369	368	3	1256	1249
4	241	232	16	148	141	11	58	80	4	651	619
5	661	660	17	26	51	12	762	749	5	126	116
6	126	106	1,12,L			13	150	125	6	607	582
7	165	166	0	80	59	14	331	336	7	884	849
8	243	225	1	184	207	15	202	209	8	224	217
9	325	316	2	53	9	16	530	555	9	326	326
10	109	106	3	311	320	17	21	5	10	394	390
11	320	297	4	129	143	2,2,L			11	694	685
12	197	192	5	48	39	1	918	998	12	158	155
13	272	264	6	102	96	2	61	80	13	578	586
14	48	61	7	262	270	3	792	799	14	167	180
15	126	140	8	129	122	4	1234	1177	15	289	297
16	121	128	9	100	73	5	399	384	16	199	200
17	104	87	10	97	83	6	687	644	17	391	402
1,10,L			11	136	148	7	539	519	2,5,L		
0	112	116	12	24	22	8	629	629	0	180	154
1	311	348	13	189	187	9	325	306	1	1125	1209
2	99	112	14	82	93	10	739	722	2	410	418
3	578	593	15	26	69	11	456	429	3	561	570
4	46	24	16	73	53	12	116	102	4	139	122
5	119	106	17	202	213	13	320	322	5	947	938
6	165	156	2,0,L			14	707	720	6	194	166
7	461	461	1	131	176	15	141	147	7	165	151
8	85	83	2	1095	1203	16	367	372	8	139	134
9	148	135	3	597	567	17	296	301	9	668	631
10	102	102	5	90	47	2,3,L			10	156	141
11	240	236	6	1168	1167	0	420	457	11	513	502
12	24	26	7	92	94	1	648	716	12	105	107
13	284	282	8	87	81	2	1095	1149	13	343	355
14	97	90	9	437	400	3	161	163	14	199	194
15	24	45	10	1234	1212	4	525	500	15	457	462
16	72	77	11	257	274	5	1205	1160	16	124	129
17	199	202	12	105	90	6	850	806	17	22	13
1,11,L			13	195	204	7	350	325	2,6,L		
0	58	41	14	751	772	8	529	519	0	51	17
1	343	382	15	192	221	9	814	780	1	510	555
2	100	89	16	313	329	10	53	61	2	372	372
3	155	174	17	105	98	11	396	396	3	1052	1061
4	22	39	2,1,L			12	714	693	4	255	242
5	398	393	0	133	162	13	495	510	5	87	81
6	56	42	1	836	915	14	71	47	6	224	209
7	95	92	3	267	263	15	427	432	7	1034	976
8	139	140	4	981	955	16	235	250	8	68	65
9	284	284	5	433	413	17	21	15	9	342	342
10	41	26							10	109	104

2,6,L			1	274	285	14	95	108	7	231	224
11	364	343	2	464	495	15	102	82	8	337	325
12	206	203	3	97	123	16	133	151	9	1056	1027
13	379	368	4	297	293	17	26	36	10	38	53
14	119	108	5	287	302	2,12,L			11	724	742
15	43	72	6	396	393	0	245	261	12	391	390
16	144	139	7	109	109	1	161	169	13	496	507
17	318	319	8	362	352	2	112	106	14	104	123
2,7,L			9	240	255	3	168	176	15	347	361
0	214	234	10	177	173	4	282	294	16	170	189
1	517	568	11	180	192	5	55	46	17	34	41
2	433	464	12	253	240	6	236	251	3,2,L		
3	267	288	13	22	36	7	136	143	0	1045	1139
4	354	357	14	121	106	8	107	101	1	700	760
5	777	753	15	70	91	9	87	118	2	427	453
6	185	175	16	219	230	10	240	239	3	654	715
7	143	139	17	72	29	11	82	82	4	699	713
8	241	231	2,10,L			12	56	67	5	127	114
9	558	546	0	364	398	13	126	102	6	399	398
10	122	109	1	66	55	14	78	111	7	1112	1089
11	313	312	2	182	198	15	26	43	8	471	478
12	389	379	3	141	133	16	99	59	9	469	469
13	126	143	4	444	449	17	117	83	10	588	575
14	49	64	5	133	155	3,0,L			11	226	206
15	284	303	6	219	215	1	1008	1111	12	214	191
16	61	88	7	85	77	2	218	214	13	474	482
17	22	43	8	163	179	3	1390	1467	14	190	191
2,8,L			9	109	117	4	31	55	15	178	175
0	637	701	10	447	435	5	547	548	16	114	97
1	425	440	11	68	67	6	87	50	17	374	395
2	32	22	12	36	49	7	1280	1286	3,3,L		
3	216	233	13	55	43	8	15	12	0	83	74
4	243	253	14	267	256	9	444	464	1	1032	1097
5	155	139	15	105	98	10	54	71	2	677	735
6	564	569	16	92	81	11	692	712	3	304	315
7	413	414	17	60	28	12	178	184	4	688	680
8	326	306	2,11,L			13	501	516	5	746	719
9	224	205	0	148	160	14	29	35	6	656	655
10	301	289	1	61	61	15	170	172	7	107	119
11	282	293	2	360	380	16	150	168	8	891	860
12	48	30	3	73	60	17	371	399	9	484	480
13	263	276	4	190	197	3,1,L			10	68	91
14	270	257	5	55	88	0	252	289	11	204	217
15	105	100	6	228	229	1	896	1003	12	620	623
16	134	131	7	53	65	2	717	794	13	184	190
17	112	162	8	384	392	3	321	325	14	194	194
2,9,L			9	122	112	4	530	505	15	362	370
0	75	38	10	75	64	5	1134	1142	16	226	228
			11	46	59	6	722	734	17	63	84
			12	241	235						
			13	97	89						

3,4,L			11	318	325	3	85	101	16	90	46
0	1011	1068	12	85	81	4	189	206	17	39	11
1	83	67	13	124	129	5	493	510	3,12,L		
2	158	178	14	381	396	6	199	207	0	58	55
3	495	505	15	177	165	7	60	20	1	250	276
4	830	854	16	178	196	8	189	173	2	75	80
5	61	59	17	65	61	9	422	428	3	146	155
6	847	828	3,7,L			10	88	103	4	76	74
7	267	266	0	66	87	11	338	341	5	41	64
8	490	467	1	325	341	12	192	195	6	58	73
9	194	182	2	576	620	13	136	134	7	235	245
10	949	916	3	330	346	14	85	67	8	127	131
11	520	518	4	372	383	15	267	253	9	55	38
12	311	321	5	326	342	16	80	95	10	133	129
13	153	157	6	388	391	17	38	26	11	180	215
14	583	586	7	185	187	3,10,L			12	53	57
15	141	148	8	503	494	0	221	215	13	158	149
16	255	264	9	430	438	1	246	281	14	46	69
17	22	45	10	21	48	2	87	43	15	26	75
3,5,L			11	185	182	3	539	571	16	26	5
0	15	13	12	549	545	4	21	14	17	100	126
1	158	158	13	107	102	5	44	35	4,0,L		
2	1198	1246	14	158	157	6	71	52	0	918	1017
3	136	153	15	214	208	7	435	443	1	65	87
4	505	527	16	274	294	8	75	42	2	54	23
5	240	215	17	112	127	9	280	288	3	270	280
6	864	849	3,8,L			10	75	77	4	602	684
7	238	221	0	343	373	11	221	223	5	382	366
8	952	916	1	301	327	12	46	37	6	1518	1551
9	39	38	2	337	350	13	265	256	7	15	22
10	255	247	3	503	524	14	112	111	8	386	430
11	93	81	4	445	453	15	51	69	9	48	54
12	442	441	5	187	195	16	63	69	10	765	792
13	178	177	6	563	576	17	236	242	11	229	231
14	156	148	7	462	462	3,11,L			12	88	26
15	87	94	8	136	120	0	138	154	13	194	197
16	369	385	9	156	153	1	371	405	14	755	792
17	55	60	10	320	314	2	78	85	15	216	232
3,6,L			11	240	238	3	182	184	16	357	376
0	1139	1222	12	44	52	4	63	45	17	21	14
1	75	74	13	194	187	5	330	333	4,1,L		
2	335	367	14	202	220	6	51	69	0	71	43
3	116	110	15	51	66	7	55	67	1	362	385
4	790	789	16	141	125	8	66	64	2	1413	1572
5	248	231	17	173	162	9	252	264	3	272	264
6	345	337	3,9,L			10	78	61	4	1076	1115
7	376	378	0	267	255	11	173	181	5	571	583
8	512	520	1	437	480	12	24	51	6	702	728
9	201	195	2	447	468	13	182	214			
10	775	769				14	61	38			
						15	129	141			

4,1,L			4,4,L			11	335	337	3	83	98
7	297	305	0	541	555	12	46	41	4	143	154
8	648	658	1	626	679	13	461	477	5	192	201
9	100	111	2	165	167	14	49	57	6	263	265
10	204	204	3	1088	1137	15	95	82	7	65	57
11	138	145	4	384	382	16	53	78	8	294	288
12	687	713	5	112	112	17	297	314	9	185	195
13	172	186	6	335	319	4,7,L			10	136	140
14	294	300	7	746	749	0	187	181	11	95	94
15	194	208	8	138	107	1	532	563	12	275	266
16	362	398	9	212	218	2	338	372	13	22	38
17	21	4	10	204	217	3	308	321	14	38	39
4,2,L			11	564	574	4	282	297	15	122	117
0	1044	1111	12	235	248	5	682	703	16	136	141
1	309	356	13	318	314	6	144	135	17	24	13
2	432	463	14	182	201	7	265	272	4,10,L		
3	682	737	15	76	81	8	257	265	0	469	496
4	1248	1289	16	126	124	9	484	487	1	60	36
5	151	138	17	270	272	10	124	142	2	97	88
6	671	708	4,5,L			11	214	230	3	102	107
7	505	529	0	219	166	12	253	239	4	287	297
8	413	422	1	617	677	13	199	205	5	55	43
9	187	200	2	241	259	14	90	90	6	314	334
10	709	730	3	284	315	15	233	231	7	109	108
11	493	499	4	243	235	16	83	80	8	209	208
12	214	226	5	874	883	17	44	21	9	75	82
13	474	498	6	124	102	4,8,L			10	272	274
14	437	458	7	180	176	0	452	491	11	105	107
15	126	129	8	297	294	1	311	337	12	24	42
16	211	214	9	479	484	2	129	148	13	100	109
17	236	245	10	63	68	3	262	278	14	258	249
4,3,L			11	466	471	4	488	537	15	24	52
0	151	136	12	44	53	5	218	220	16	90	97
1	733	789	13	379	385	6	231	239	17	39	11
2	539	595	14	61	54	7	320	326	4,11,L		
3	731	746	15	416	437	8	51	59	0	60	56
4	221	240	16	65	63	9	131	141	1	41	45
5	882	921	17	104	92	10	331	311	2	258	275
6	365	382	4,6,L			11	228	244	3	22	52
7	150	180	0	250	253	12	44	38	4	263	271
8	678	673	1	568	600	13	212	224	5	92	103
9	498	512	2	80	83	14	299	306	6	294	311
10	277	274	3	736	788	15	76	56	7	22	39
11	216	239	4	158	175	16	148	136	8	263	273
12	668	687	5	65	65	17	131	144	9	117	101
13	447	454	6	178	180	4,9,L			10	22	10
14	21	32	7	614	620	0	107	102	11	105	118
15	467	478	8	151	164	1	99	107	12	224	231
16	229	254	9	410	407	2	552	599	13	87	85
17	76	82	10	158	153						

4,11,L			5	734	794	5,4,L			11	136	145
			6	328	319				12	21	44
14	133	128	7	172	181	0	1015	1055	13	121	111
15	24	16	8	330	350	1	435	437	14	362	381
16	114	120	9	864	923	2	320	350	15	109	95
17	73	18	10	189	191	3	127	135	16	112	116
4,12,L			11	467	489	4	714	734	17	75	56
			12	367	381	5	258	253	5,7,L		
0	212	217	13	284	287	6	522	519			
1	133	146	14	85	108	7	432	454	0	82	64
2	168	177	15	371	378	8	435	449	1	257	260
3	163	166	16	73	71	9	131	131	2	576	603
4	340	344	17	22	22	10	726	743	3	127	129
5	109	103	5,2,L			11	161	179	4	291	310
6	199	202				12	184	187	5	381	384
7	105	105	0	741	752	13	184	182	6	386	400
8	80	90	1	388	416	14	352	369	7	189	199
9	72	66	2	80	98	15	22	46	8	362	373
10	178	173	3	772	842	16	177	185	9	275	290
11	122	77	4	603	623	17	163	169	10	122	131
12	80	82	5	156	178	5,5,L			11	122	145
13	129	116	6	252	259				12	355	370
14	73	86	7	663	705	0	17	9	13	92	92
15	70	57	8	156	156	1	165	158	14	131	138
16	124	109	9	461	467	2	789	846	15	102	104
17	66	67	10	474	489	3	61	62	16	209	223
5,0,L			11	444	476	4	513	521	17	24	27
			12	258	260	5	75	70	5,8,L		
1	933	1017	13	444	453	6	675	703			
2	14	5	14	231	253	7	243	247	0	478	480
3	894	1012	15	44	58	8	653	659	1	122	146
4	267	251	16	87	108	9	185	183	2	270	296
5	704	703	17	270	289	10	146	160	3	371	382
6	87	56	5,3,L			11	114	107	4	320	330
7	858	909				12	452	471	5	158	172
8	82	89	0	117	87	13	43	52	6	340	348
9	377	406	1	590	626	14	127	135	7	498	500
10	126	124	2	801	837	15	102	116	8	158	170
11	345	354	3	199	180	16	321	338	9	138	136
12	311	325	4	304	317	17	116	98	10	258	270
13	415	433	5	692	718	5,6,L			11	228	235
14	109	90	6	518	534				12	51	42
15	204	209	7	99	103	0	624	651	13	168	173
16	21	12	8	801	807	1	284	294	14	126	122
17	379	417	9	269	269	2	250	259	15	58	88
5,1,L			10	53	75	3	19	47	16	129	146
			11	345	354	4	537	560	17	170	162
0	43	21	12	369	386	5	95	97	5,9,L		
1	838	943	13	105	111	6	573	583			
2	104	108	14	129	140	7	316	310	0	82	74
3	447	469	15	253	268	8	449	445			
4	199	212	16	184	210	9	70	94			
			17	22	36	10	554	560			

5,9,L			12	73	46	4	592	628	17	46	42
1	347	365	13	114	124	5	286	313	6,4,L		
2	275	295	14	49	50	6	544	567	0	369	358
3	143	156	15	141	139	7	121	129	1	432	450
4	207	222	16	66	51	8	627	663	2	192	195
5	413	418	17	26	19	9	66	75	3	688	729
6	75	59	5,12,L			10	182	173	4	333	332
7	131	150	0	82	99	11	204	238	5	104	117
8	201	195	1	143	148	12	510	560	6	182	192
9	291	312	2	76	65	13	178	185	7	578	594
10	22	34	3	240	237	14	143	148	8	92	75
11	343	341	4	38	54	15	158	172	9	170	175
12	201	200	5	126	133	16	325	367	10	337	353
13	148	147	6	92	96	17	70	76	11	399	398
14	49	16	7	114	121	6,2,L			12	156	166
15	201	211	8	100	106	0	569	577	13	382	394
16	49	69	9	143	144	1	318	346	14	63	79
17	60	69	10	90	79	2	173	186	15	160	176
5,10,L			11	187	193	3	512	543	16	68	75
0	117	118	12	24	18	4	529	560	17	355	375
1	207	214	13	121	116	5	184	186	6,5,L		
2	53	51	14	72	65	6	697	715	0	56	56
3	405	432	15	87	88	7	529	549	1	636	677
4	82	71	16	26	34	8	331	337	2	134	126
5	22	52	17	92	110	9	206	225	3	325	338
6	68	82	6,0,L			10	508	531	4	51	19
7	416	421	0	1197	1252	11	289	290	5	658	663
8	22	50	1	185	189	12	155	161	6	136	118
9	139	143	2	573	610	13	167	172	7	365	386
10	105	106	3	110	90	14	452	492	8	260	263
11	199	216	4	845	923	15	109	118	9	544	552
12	88	98	5	224	218	16	148	156	10	88	99
13	148	160	6	622	690	17	143	141	11	308	327
14	73	106	7	46	28	6,3,L			12	63	83
15	48	61	8	376	408	0	102	102	13	318	319
16	55	27	9	19	31	1	620	645	14	139	136
17	124	138	10	643	706	2	541	560	15	255	270
5,11,L			11	48	40	3	277	281	16	83	68
0	71	48	12	100	37	4	320	324	17	38	41
1	325	341	13	21	13	5	619	636	6,6,L		
2	68	48	14	309	346	6	342	334	0	155	160
3	127	140	15	21	13	7	151	130	1	493	525
4	22	12	16	211	227	8	433	452	2	233	233
5	223	235	17	22	9	9	371	385	3	685	693
6	104	107	6,1,L			10	182	178	4	122	131
7	82	65	0	34	31	11	211	214	5	134	126
8	49	34	1	311	342	12	355	374	6	177	192
9	279	285	2	881	939	13	138	147	7	430	446
10	110	112	3	209	220	14	143	135			
11	104	98				15	357	378			
						16	110	130			

6,6,L			6,9,L			11 24 34			4 78 61		
8	163	166	0	38	32	12	163	160	5	739	793
9	284	310	1	117	105	13	112	101	6	187	192
10	93	106	2	280	295	14	109	92	7	78	91
11	323	339	3	161	163	15	44	28	8	212	194
12	82	95	4	223	233	16	82	92	9	423	446
13	299	306	5	240	236	17	26	21	10	71	59
14	76	67	6	250	254	6,12,L			11	367	398
15	78	81	7	68	58	0	212	209	12	99	94
16	58	62	8	301	300	1	87	83	13	243	248
17	190	196	9	114	109	2	68	45	14	82	81
6,7,L			10	182	197	3	95	94	15	292	312
0	55	65	11	104	117	4	163	164	16	53	58
1	398	403	12	245	238	5	75	51	17	43	23
2	479	494	13	76	54	6	218	222	7,2,L		
3	223	223	14	85	55	7	102	121	0	306	320
4	53	69	15	110	104	8	124	111	1	716	756
5	440	446	16	158	156	9	93	98	2	126	115
6	104	121	17	44	8	10	141	156	3	610	647
7	129	116	6,10,L			11	110	118	4	425	449
8	206	195	0	376	386	12	39	60	5	158	161
9	357	361	1	32	9	13	58	64	6	245	275
10	117	120	2	61	85	14	90	77	7	537	567
11	204	218	3	46	56	15	26	48	8	204	209
12	209	202	4	406	417	16	26	56	9	202	228
13	127	128	5	22	47	17	65	40	10	299	320
14	93	76	6	134	152	7,0,L			11	296	309
15	199	223	7	99	94	1	289	307	12	56	49
16	112	107	8	107	100	2	41	73	13	289	315
17	24	27	9	92	82	3	524	591	14	269	285
6,8,L			10	252	240	4	38	15	15	85	78
0	520	546	11	24	48	5	49	68	16	144	137
1	245	249	12	24	45	6	126	114	17	153	160
2	51	27	13	93	72	7	653	700	7,3,L		
3	241	234	14	211	229	8	48	10	0	114	132
4	311	302	15	72	54	9	260	273	1	372	383
5	95	98	16	99	110	10	107	140	2	661	683
6	197	195	17	26	25	11	433	454	3	207	209
7	277	283	6,11,L			12	90	92	4	352	349
8	151	174	0	22	17	13	445	456	5	411	435
9	117	115	1	51	54	14	39	16	6	433	435
10	206	214	2	330	350	15	63	54	7	207	210
11	121	137	3	46	43	16	34	12	8	512	529
12	73	71	4	155	165	17	328	361	9	338	362
13	206	216	5	68	71	7,1,L			10	109	120
14	182	194	6	216	227	0	141	174	11	161	168
15	95	70	7	24	42	1	583	614	12	440	446
16	109	96	8	207	211	2	194	183	13	116	120
17	117	138	9	90	89	3	432	474	14	66	56
			10	34	11				15	246	257



7,3,L			6	445	440	7,9,L			11	168	160
16	124	103	7	121	106	0	22	37	12	56	36
17	22	12	8	270	283	1	335	328	13	146	158
7,4,L			9	87	103	2	153	153	14	90	54
0	767	777	10	379	373	3	226	223	15	127	131
1	122	127	11	65	53	4	184	188	16	26	47
2	122	108	12	85	95	5	374	371	17	26	42
3	371	384	13	60	36	6	39	5	7,12,L		
4	644	655	14	313	316	7	117	119	0	24	42
5	133	139	15	44	55	8	202	188	1	180	162
6	508	508	16	211	194	9	229	232	2	24	36
7	82	72	17	55	53	10	24	31	3	136	141
8	275	272	7,7,L			11	105	105	4	146	155
9	126	113	0	80	53	12	102	99	5	60	53
10	515	519	1	76	54	13	90	85	6	24	42
11	53	37	2	498	512	14	34	2	7	190	183
12	43	60	3	161	161	15	109	113	8	70	66
13	126	132	4	286	285	16	51	39	9	24	42
14	246	245	5	265	270	17	65	56	10	24	22
15	70	60	6	347	342	7,10,L			11	95	110
16	85	94	7	75	79	0	83	50	12	38	14
17	148	155	8	343	340	1	156	153	13	93	98
7,5,L			9	311	301	2	83	100	14	70	68
0	44	56	10	102	116	3	301	302	15	41	57
1	214	217	11	105	81	4	60	31	16	26	50
2	532	531	12	265	264	5	88	79	17	119	113
3	63	42	13	87	92	6	68	37	8,0,L		
4	275	281	14	78	83	7	277	280	0	1039	1059
5	21	27	15	160	174	8	43	36	1	168	180
6	440	437	16	150	163	9	95	103	2	313	332
7	36	54	17	60	94	10	24	12	3	76	73
8	425	430	7,8,L			11	173	177	4	603	627
9	53	61	0	415	420	12	66	76	5	29	36
10	92	101	1	185	172	13	102	97	6	435	443
11	63	35	2	112	99	14	24	34	7	136	142
12	391	408	3	287	286	15	88	92	8	250	250
13	131	125	4	240	245	16	26	38	9	277	318
14	219	211	5	41	52	17	116	92	10	469	495
15	122	127	6	192	208	7,11,L			11	21	50
16	241	248	7	364	359	0	66	21	12	31	16
17	38	39	8	78	91	1	173	181	13	48	1
7,6,L			9	124	117	2	39	50	14	301	338
0	403	382	10	253	249	3	24	29	15	68	57
1	38	38	11	236	242	4	48	10	16	56	59
2	185	200	12	75	24	5	236	233	17	49	29
3	190	180	13	241	239	6	51	72	8,1,L		
4	478	482	14	156	156	7	83	79	0	287	296
5	170	171	15	105	94	8	58	60	1	398	415
			16	104	91	9	151	159			
			17	77	95	10	24	28			

8,1,L			13	178	168	5	110	95	8,9,L		
2	588	614	14	61	65	6	68	62	0	83	83
3	116	103	15	241	248	7	440	425	1	121	93
4	376	397	16	151	165	8	119	125	2	287	274
5	224	234	17	82	99	9	153	156	3	105	102
6	384	389	8,4,L			10	110	101	4	104	100
7	49	34	0	282	296	11	170	172	5	153	132
8	517	548	1	326	312	12	22	13	6	216	214
9	245	268	2	56	48	13	189	187	7	95	65
10	153	143	3	462	475	14	146	136	8	156	156
11	21	18	4	156	157	15	78	93	9	134	130
12	362	378	5	63	46	16	72	49	10	55	73
13	76	100	6	270	270	17	156	165	11	87	67
14	75	77	7	493	488	8,7,L			12	255	240
15	99	96	8	156	150	0	41	16	13	68	17
16	238	249	9	168	154	1	282	290	14	63	41
17	46	69	10	202	207	2	287	301	15	99	103
8,2,L			11	323	313	3	104	94	16	180	171
0	471	464	12	21	10	4	165	159	17	26	19
1	297	293	13	262	266	5	320	307	8,10,L		
2	260	276	14	104	112	6	170	167	0	325	318
3	473	483	15	148	159	7	22	32	1	73	48
4	321	324	16	22	22	8	178	174	2	58	69
5	163	173	17	296	290	9	214	226	3	48	51
6	274	301	8,5,L			10	22	46	4	199	194
7	241	249	0	92	62	11	253	251	5	56	59
8	250	248	1	643	620	12	110	97	6	173	158
9	192	179	2	136	153	13	116	131	7	41	35
10	393	402	3	201	190	14	49	43	8	139	146
11	177	179	4	126	115	15	144	141	9	110	111
12	82	38	5	530	518	16	73	68	10	223	221
13	126	124	6	68	47	17	24	13	11	63	53
14	226	243	7	99	95	8,8,L			12	80	69
15	97	97	8	48	67	0	365	362	13	55	70
16	177	189	9	394	393	1	82	80	14	177	172
17	121	142	10	21	35	2	124	118	15	61	45
8,3,L			11	214	211	3	250	233	16	104	95
0	141	147	12	58	68	4	224	210	17	39	41
1	415	417	13	160	150	5	38	26	8,11,L		
2	418	413	14	34	45	6	194	190	0	24	6
3	143	158	15	192	192	7	226	213	1	24	33
4	221	221	16	75	84	8	141	133	2	253	257
5	343	342	17	34	26	9	114	119	3	24	12
6	253	250	8,6,L			10	161	157	4	121	114
7	178	190	0	216	234	11	146	155	5	73	20
8	269	262	1	428	407	12	24	10	6	161	157
9	345	331	2	41	52	13	150	123	7	44	33
10	21	31	3	457	454	14	110	106	8	187	202
11	228	235	4	75	57	15	24	23			
12	201	216				16	100	87			
						17	90	94			

8,11,L			0	197	208	13	155	147	5	121	113
9	88	98	1	467	460	14	122	125	6	187	182
10	75	55	2	151	139	15	107	110	7	151	141
11	46	46	3	309	307	16	127	137	8	99	84
12	105	112	4	116	85	17	63	30	9	60	41
13	43	61	5	551	539	9,4,L			10	304	297
14	78	53	6	153	158	0	394	381	11	92	78
15	73	67	7	105	119	1	246	250	12	22	40
16	41	62	8	90	91	2	138	124	13	70	57
17	58	44	9	379	374	3	73	49	14	282	270
8,12,L			10	82	59	4	513	505	15	39	29
0	185	184	11	195	214	5	73	58	16	99	94
1	73	66	12	104	108	6	364	356	17	24	17
2	66	51	13	136	134	7	143	136	9,7,L		
3	117	98	14	38	49	8	163	158	0	110	110
4	189	196	15	272	284	9	76	73	1	148	144
5	24	26	16	92	93	10	291	300	2	326	317
6	85	80	17	24	40	11	112	102	3	129	119
7	60	61	9,2,L			12	22	23	4	182	169
8	70	62	0	236	228	13	133	132	5	195	168
9	48	41	1	219	200	14	199	202	6	160	144
10	163	149	2	53	56	15	104	91	7	75	62
11	90	85	3	575	567	16	99	103	8	280	266
12	65	42	4	219	223	17	51	69	9	131	116
13	83	61	5	114	115	9,5,L			10	109	92
14	90	75	6	286	285	0	21	46	11	114	116
15	100	78	7	416	405	1	80	87	12	216	202
16	26	32	8	229	237	2	399	392	13	139	140
17	70	47	9	173	176	3	71	43	14	63	63
9,0,L			10	253	248	4	224	204	15	95	72
1	396	389	11	229	240	5	21	16	16	119	112
2	133	106	12	75	72	6	292	277	17	77	56
3	389	389	13	172	189	7	73	72	9,8,L		
4	19	8	14	139	133	8	372	362	0	248	232
5	19	19	15	32	28	9	75	62	1	185	176
6	75	77	16	58	79	10	105	93	2	112	100
7	503	517	17	129	148	11	22	24	3	212	198
8	51	66	9,3,L			12	330	313	4	195	192
9	129	145	0	216	201	13	68	67	5	126	101
10	44	48	1	289	284	14	53	61	6	218	219
11	304	309	2	483	489	15	43	34	7	163	156
12	38	39	3	165	164	16	102	114	8	83	87
13	342	339	4	209	198	17	24	33	9	95	82
14	76	79	5	357	352	9,6,L			10	112	109
15	48	67	6	314	309	0	449	433	11	167	152
16	66	26	7	87	60	1	92	115	12	56	34
17	175	185	8	280	272	2	83	74	13	144	133
9,1,L			9	180	178	3	39	60	14	66	94
0	197	208	10	70	53	4	359	348	15	24	20
1	467	460	11	185	198				16	72	95
2	151	139	12	267	261						

9,8,L			7	99	97	10,1,L			11	119	101
17	90	92	8	95	76	0	48	44	12	139	125
9,9,L			9	141	131	1	212	194	13	216	204
0	24	6	10	41	32	2	500	460	14	100	87
1	192	183	11	26	49	3	60	50	15	122	111
2	150	146	12	26	8	4	241	236	16	148	142
3	144	133	13	70	68	5	190	170	17	46	51
4	105	107	14	26	7	6	279	263	10,4,L		
5	270	260	15	136	114	7	43	40	0	71	50
6	110	105	16	53	32	8	289	271	1	335	296
7	56	77	17	55	17	9	55	41	2	150	148
8	66	83	9,12,L			10	41	47	3	367	341
9	218	193	0	102	72	11	124	114	4	151	147
10	82	60	1	73	33	12	291	278	5	112	89
11	97	68	2	44	71	13	92	87	6	136	116
12	26	61	3	172	154	14	22	30	7	309	274
13	95	71	4	80	54	15	53	23	8	100	68
14	78	59	5	48	29	16	119	115	9	95	94
15	82	86	6	72	73	17	83	73	10	48	53
16	41	46	7	90	87	10,2,L			11	236	213
17	90	74	8	55	29	0	525	496	12	22	54
9,10,L			9	66	40	1	253	241	13	173	164
0	51	45	10	83	76	2	80	60	14	34	35
1	180	179	11	99	95	3	255	239	15	102	81
2	63	51	12	26	16	4	320	300	16	72	76
3	209	208	13	27	63	5	76	55	17	151	147
4	63	28	14	75	76	6	185	171	10,5,L		
5	75	49	15	60	41	7	178	165	0	60	41
6	24	35	16	87	76	8	151	141	1	335	311
7	211	193	17	63	64	9	104	107	2	112	92
8	24	5	10,0,L			10	260	246	3	121	95
9	72	97	0	508	491	11	124	117	4	22	17
10	68	57	1	71	66	12	22	18	5	369	328
11	92	95	2	173	164	13	143	152	6	88	96
12	41	32	3	21	39	14	228	212	7	66	51
13	139	125	4	355	335	15	107	84	8	22	36
14	43	35	5	99	92	16	82	80	9	240	222
15	109	111	6	323	317	17	119	112	10	80	70
16	68	32	7	38	7	10,3,L			11	168	147
17	97	106	8	127	124	0	88	80	12	82	61
9,11,L			9	70	43	1	301	283	13	80	79
0	70	54	10	282	292	2	219	203	14	95	91
1	143	134	11	109	98	3	194	174	15	258	234
2	75	84	12	41	24	4	117	115	16	44	65
3	95	93	13	104	86	5	328	303	17	44	58
4	77	83	14	212	219	6	138	120	10,6,L		
5	206	180	15	88	97	7	102	98	0	155	146
6	87	49	16	153	144	8	236	232	1	122	82
			17	24	32	9	297	289	2	22	8
						10	105	95			



11,3,L			1	100	90	14	78	88	6	26	10
11	124	120	2	36	16	15	68	58	7	26	25
12	95	93	3	138	113	16	80	44	8	95	70
13	85	88	4	235	207	17	63	50	9	73	58
14	124	117	5	138	108	11,9,L			10	26	24
15	78	52	6	160	130	0	26	26	11	27	43
16	131	117	7	24	42	1	124	122	12	26	16
17	85	37	8	151	118	2	134	117	13	41	49
11,4,L			9	109	97	3	76	33	14	38	31
0	311	268	10	201	186	4	68	57	15	41	54
1	88	48	11	78	65	5	172	153	16	26	13
2	134	113	12	24	34	6	63	35	17	26	24
3	199	169	13	24	8	7	87	81	11,12,L		
4	265	233	14	122	130	8	77	56	0	26	27
5	117	77	15	26	25	9	153	134	1	26	29
6	214	187	16	68	65	10	48	22	2	61	55
7	78	73	17	60	32	11	102	85	3	87	84
8	134	112	11,7,L			12	73	39	4	61	42
9	90	91	0	56	38	13	77	44	5	82	44
10	277	248	1	141	112	14	53	59	6	41	23
11	143	141	2	187	167	15	85	76	7	82	75
12	85	59	3	170	163	16	46	61	8	43	28
13	83	65	4	117	97	17	26	40	9	78	30
14	165	145	5	44	52	11,10,L			10	83	61
15	49	51	6	107	93	0	73	53	11	70	47
16	61	76	7	66	47	1	148	134	12	55	21
17	73	48	8	184	157	2	95	68	13	27	40
11,5,L			9	88	86	3	180	153	14	55	31
0	63	51	10	78	52	4	38	33	15	27	40
1	102	85	11	117	85	5	80	58	16	27	15
2	359	326	12	136	131	6	56	5	17	60	49
3	41	25	13	55	61	7	107	102	12,0,L		
4	129	118	14	65	48	8	26	7	0	284	247
5	71	27	15	26	16	9	85	73	1	22	42
6	246	212	16	68	68	10	55	22	2	121	96
7	24	53	17	63	25	11	61	63	3	22	21
8	236	206	11,8,L			12	26	33	4	304	270
9	80	57	0	219	201	13	143	103	5	88	79
10	43	28	1	199	186	14	26	21	6	153	125
11	60	16	2	24	36	15	26	17	7	22	19
12	280	253	3	90	70	16	26	15	8	22	17
13	85	84	4	131	108	17	66	72	9	55	27
14	24	39	5	39	28	11,11,L			10	204	189
15	51	23	6	107	95	0	99	88	11	92	95
16	85	66	7	156	126	1	167	145	12	24	41
17	24	11	8	68	19	2	66	16	13	77	87
11,6,L			9	48	33	3	58	68	14	177	155
0	364	336	10	88	88	4	70	15	15	148	134
			11	78	66	5	124	120	16	87	91
			12	82	48				17	24	19
			13	102	79						

12,1,L			11	114	76	3	257	208	16	26	34
0	22	14	12	92	90	4	51	26	17	68	62
1	83	61	13	148	131	5	46	33	12,9,L		
2	318	289	14	95	72	6	80	70	0	66	34
3	110	87	15	121	88	7	207	174	1	53	45
4	158	127	16	102	91	8	24	11	2	158	129
5	73	58	17	56	32	9	109	84	3	65	41
6	146	121	12,4,L			10	51	43	4	104	95
7	93	87	0	114	100	11	121	90	5	61	47
8	192	170	1	216	186	12	51	26	6	150	125
9	39	31	2	38	43	13	151	127	7	66	20
10	92	53	3	246	201	14	72	42	8	92	65
11	80	72	4	36	38	15	26	27	9	26	68
12	211	179	5	24	42	16	26	18	10	75	37
13	24	40	6	77	80	17	109	75	11	65	31
14	55	54	7	236	202	12,7,L			12	131	111
15	24	32	8	24	18	0	60	53	13	26	18
16	82	77	9	82	60	1	167	127	14	63	5
17	24	24	10	48	43	2	93	88	15	26	28
12,2,L			11	156	140	3	68	61	16	68	44
0	299	264	12	90	75	4	146	107	17	46	2
1	93	87	13	138	112	5	233	196	12,10,L		
2	55	13	14	61	31	6	68	44	0	95	74
3	209	186	15	83	50	7	24	22	1	63	33
4	265	220	16	63	40	8	129	92	2	39	51
5	63	49	17	92	88	9	105	87	3	70	6
6	194	161	12,5,L			10	26	48	4	167	134
7	109	89	0	24	31	11	102	92	5	26	13
8	138	115	1	190	163	12	88	77	6	121	100
9	66	74	2	24	8	13	87	74	7	85	35
10	131	112	3	105	64	14	95	40	8	107	65
11	155	135	4	24	34	15	65	65	9	27	51
12	75	52	5	280	227	16	80	49	10	97	86
13	150	122	6	24	34	17	26	12	11	56	15
14	150	134	7	112	97	12,8,L			12	26	23
15	46	19	8	24	10	0	134	123	13	26	34
16	82	64	9	201	176	1	24	17	14	102	81
17	75	61	10	70	63	2	24	37	15	27	30
12,3,L			11	95	55	3	131	109	16	53	44
0	22	15	12	26	12	4	143	105	17	27	4
1	236	207	13	92	64	5	70	28	12,11,L		
2	197	170	14	41	27	6	143	131	0	26	36
3	156	117	15	151	145	7	161	128	1	72	22
4	80	69	16	26	18	8	80	74	2	121	108
5	252	216	17	26	20	9	44	35	3	26	48
6	167	132	12,6,L			10	122	85	4	78	62
7	34	33	0	83	76	11	121	80	5	68	32
8	83	69	1	170	123	12	26	34	6	70	42
9	155	132	2	44	34	13	92	61			
10	70	37				14	78	84			
						15	61	43			

APPENDIX C

X-RAY DATA OF COMPOUND I



-10,0,L			2	70	66	-10,6,L			1	17*	11
0	39*	25	3	141	124	1	98	94	2	32*	39
2	250	246	4	126	125	2	39*	29	3	128	113
4	13*	22	5	211	201	3	109	111	4	49*	47
6	90	92	6	309	314	4	96	111	5	19*	29
8	28*	36	7	39*	47	5	15*	29	6	49*	53
10	17*	1	8	85	100	6	32*	19	7	19*	21
12	51*	21	9	17*	22	7	54*	40	8	19*	12
14	19*	6	10	145	144	8	17*	2	9	19*	3
16	45*	17	11	98	92	9	51*	58	10	19*	23
-10,1,L			12	19*	12	10	17*	13	11	22*	7
1	36*	26	13	19*	19	11	19*	11	12	54*	3
2	81	83	14	22*	20	12	54*	15	13	36*	16
3	24*	30	15	32*	16	13	32*	31	-10,10,L		
4	13*	17	16	22*	12	14	22*	2	1	92	80
5	75	77	-10,4,L			15	19*	5	2	43*	8
6	122	129	1	198	203	-10,7,L			3	139	125
7	62	62	2	130	131	1	105	93	4	64	62
8	43*	32	3	109	99	2	53	52	5	34*	7
9	17*	11	4	15*	3	3	39*	23	6	19*	1
10	17*	18	5	15*	6	4	60	41	7	19*	8
11	17*	26	6	70	67	5	70	35	8	19*	15
12	19*	11	7	47*	47	6	17*	1	9	30*	2
13	51*	3	8	81	79	7	17*	2	10	22*	13
14	45*	18	9	17*	24	8	94	109	11	22*	1
15	32*	4	10	66	45	9	30*	8	12	22*	3
16	43*	6	11	17*	16	10	19*	5	13	22*	14
17	19*	3	12	39*	14	11	19*	4	-9,0,L		
-10,2,L			13	19*	7	12	51*	11	1	495	487
1	32*	39	14	32*	4	13	22*	11	3	13*	2
2	24*	13	15	19*	19	14	22*	12	5	107	104
3	13*	1	16	22*	11	15	22*	6	7	68	70
4	132	125	-10,5,L			-10,8,L			9	15*	9
5	111	115	1	81	76	1	17*	34	11	47*	55
6	226	228	2	43*	45	2	17*	1	13	19*	9
7	45*	37	3	15*	8	3	45*	56	15	45*	44
8	41*	32	4	58	69	4	54*	53	17	41*	1
9	39*	0	5	15*	8	5	43*	31	-9,1,L		
10	119	128	6	15*	8	6	17*	6	1	75	87
11	51*	51	7	47*	24	7	36*	38	2	326	316
12	17*	6	8	51*	25	8	19*	16	3	132	130
13	19*	2	9	34*	10	9	19*	20	4	87	85
14	19*	1	10	17*	10	10	19*	5	5	28*	18
15	22*	4	11	19*	0	11	19*	1	6	81	82
16	49*	1	12	19*	4	12	56*	9	7	15*	6
17	22*	6	13	47*	14	13	51*	44	8	39*	17
-10,3,L			14	19*	8	14	22*	4	9	15*	17
1	158	149	15	19*	2	-10,9,L			10	55	57
			16	22*	9				11	68	68



-8,2,L			7	49*	52	6	45*	30	-7,1,L		
13	19*	18	8	17*	8	7	17*	10	1	215	212
14	19*	12	9	36*	14	8	45*	42	2	55	61
15	19*	9	10	58*	28	9	19*	25	3	241	244
16	22*	4	11	32*	62	10	19*	20	4	196	189
17	22*	2	12	45*	18	11	43*	26	5	39*	30
-8,3,L			13	54*	20	12	43*	4	6	137	126
1	81	79	14	32*	12	13	39*	28	7	39*	43
2	205	194	15	41*	14	14	22*	30	8	73	71
3	43*	49	16	19*	11	-8,9,L			9	75	66
4	24*	14	-8,6,L			1	60	43	10	41*	34
5	24*	3	1	41*	30	2	49*	34	11	55	19
6	222	220	2	39*	39	3	17*	18	12	19*	19
7	26*	30	3	55	69	4	70	44	13	45*	18
8	258	265	4	177	186	5	58*	62	14	47*	24
9	143	146	5	36*	20	6	19*	9	15	19*	7
10	17*	6	6	68	85	7	19*	9	16	54*	9
11	19*	9	7	17*	24	8	47*	6	17	22*	1
12	19*	28	8	28*	46	9	34*	2	-7,2,L		
13	19*	1	9	54*	11	10	47*	17	1	11*	3
14	28*	15	10	19*	42	11	22*	26	2	90	87
15	19*	17	11	32*	27	12	22*	12	3	503	489
16	39*	12	12	19*	3	13	49*	47	4	469	465
17	19*	4	13	19*	1	14	45*	4	5	64	70
-8,4,L			14	19*	9	-8,10,L			6	15*	35
1	39*	24	15	36*	1	1	117	91	7	15*	40
2	43*	36	16	47*	3	2	94	74	8	45*	25
3	109	107	-8,7,L			3	79	62	9	17*	4
4	132	126	1	147	154	4	41*	40	10	43*	12
5	49*	39	2	77	75	5	43*	12	11	26*	16
6	75	64	3	34*	9	6	30*	3	12	19*	25
7	43*	8	4	39*	13	7	19*	1	13	43*	28
8	98	97	5	36*	10	8	19*	4	14	68	4
9	83	98	6	100	106	9	19*	3	15	19*	2
10	17*	21	7	17*	17	10	19*	9	16	45*	7
11	34*	17	8	79	76	11	22*	7	17	22*	8
12	30*	18	9	36*	20	12	22*	2	-7,3,L		
13	19*	38	10	19*	3	13	22*	8	1	13*	13
14	19*	18	11	47*	10	-7,0,L			2	100	89
15	30*	23	12	43*	5	1	166	172	3	245	241
16	36*	20	13	36*	1	3	109	102	4	218	211
-8,5,L			14	22*	6	5	232	235	5	15*	17
1	62	63	15	49*	10	7	15*	2	6	15*	34
2	53	57	-8,8,L			9	34*	34	7	43*	45
3	66	72	1	264	262	11	17*	1	8	190	193
4	68	66	2	154	150	13	32*	37	9	171	173
5	62	64	3	51*	33	15	19*	15	10	30*	14
6	15*	12	4	17*	34	17	62*	1	11	19*	19
			5	17*	8				12	47*	2

-7,3,L			8	39*	9	9	19*	5	-6,2,L		
13	19*	9	9	17*	3	10	58*	35	1	11*	21
14	19*	5	10	17*	1	11	54*	31	2	403	390
15	43*	5	11	19*	17	12	22*	13	3	288	271
16	49*	21	12	19*	5	13	49*	12	4	800	786
17	39*	15	13	39*	13	14	22*	3	5	205	206
-7,4,L			14	19*	9	-7,10,L			6	126	124
1	109	102	15	22*	2	1	54*	38	7	49*	51
2	68	64	16	22*	0	2	39*	34	8	15*	16
3	28*	42	-7,7,L			3	28*	30	9	17*	28
4	15*	22	1	66	64	4	56*	31	10	17*	30
5	15*	26	2	166	181	5	51*	12	11	17*	4
6	26*	24	3	39*	58	6	19*	24	12	41*	46
7	17*	32	4	58	55	7	19*	10	13	19*	27
8	122	133	5	28*	9	8	19*	23	14	45*	1
9	113	126	6	17*	12	9	19*	8	15	19*	13
10	17*	17	7	81	83	10	36*	10	16	22*	4
11	75	58	8	28*	19	11	58*	28	17	22*	1
12	34*	10	9	43*	18	12	34*	5	-6,3,L		
13	19*	7	10	19*	16	13	22*	6	1	66	62
14	45*	32	11	19*	2	-6,0,L			2	388	383
15	39*	27	12	19*	16	2	19*	17	3	51	48
16	22*	26	13	22*	7	4	130	126	4	34*	44
-7,5,L			14	19*	3	6	34*	45	5	26*	21
1	83	90	15	45*	5	8	158	154	6	70	65
2	15*	7	-7,8,L			10	81	84	7	43*	51
3	49	16	1	145	131	12	45*	21	8	28*	28
4	30*	16	2	147	138	14	19*	38	9	128	116
5	36*	6	3	34*	22	16	22*	7	10	17*	12
6	126	124	4	83	68	-6,1,L			11	17*	12
7	75	70	5	17*	3	1	166	162	12	17*	23
8	34*	15	6	85	81	2	11*	12	13	19*	4
9	17*	23	7	34*	32	3	352	343	14	19*	6
10	32*	29	8	17*	2	4	403	398	15	22*	18
11	19*	31	9	34*	2	5	232	237	16	22*	18
12	19*	6	10	19*	4	6	102	104	17	22*	2
13	47*	25	11	19*	7	7	53	57	-6,4,L		
14	19*	16	12	41*	24	8	73	68	1	352	356
15	19*	3	13	34*	7	9	169	166	2	55	53
16	30*	16	14	22*	14	10	17*	19	3	94	87
-7,6,L			-7,9,L			11	32*	20	4	15*	4
1	126	126	1	17*	18	12	17*	0	5	26*	28
2	66	72	2	54*	23	13	34*	35	6	156	151
3	45*	28	3	32*	29	14	19*	34	7	85	87
4	15*	26	4	17*	1	15	19*	1	8	143	144
5	70	49	5	60	51	16	22*	2	9	96	92
6	64	72	6	36*	7	17	22*	3	10	17*	12
7	68	71	7	17*	9	-7,10,L			11	17*	19
			8	19*	7	1	54*	38	12	28*	17
						2	39*	34	13	19*	8
						3	28*	30			
						4	56*	31			
						5	51*	12			
						6	19*	24			
						7	19*	10			
						8	19*	23			
						9	19*	8			
						10	36*	10			
						11	58*	28			
						12	34*	5			
						13	22*	6			

-6,4,L			10	19*	10	13	22*	10	-5,3,L		
14	19*	1	11	28*	1				1	572	560
15	22*	7	12	19*	8		-5,0,L		2	262	248
16	51*	13	13	19*	0	1	1143	1155	3	45	40
-6,5,L			15	22*	2	3	151	146	4	53	58
1	100	91	-6,8,L			5	132	131	5	58	47
2	15*	9				7	45*	24	6	15*	3
3	58	68	1	81	59	9	169	175	7	92	105
4	15*	23	2	43*	36	11	58	45	8	15*	25
5	75	72	3	47*	24	13	39*	35	9	26*	16
6	87	82	4	34*	6	15	19*	3	10	17*	26
7	15*	16	5	17*	28	17	22*	3	11	19*	32
8	17*	14	6	124	120	-5,1,L			12	19*	10
9	60	39	7	56*	19	1	70	69	13	51*	36
10	43*	10	8	19*	4	2	386	379	14	58*	45
11	54*	58	9	19*	27	3	85	88	15	22*	5
12	41*	16	10	19*	2	4	352	344	16	22*	1
13	36*	40	11	43*	32	5	13*	14	17	22*	4
14	30*	9	12	39*	44	6	15*	9	-5,4,L		
15	43*	0	13	22*	11	7	51	48	1	79	76
16	51*	2	14	22*	3	8	15*	25	2	22*	23
-6,6,L			-6,9,L			9	17*	31	3	87	81
1	45*	57	1	36*	28	10	41*	12	4	92	80
2	15*	35	2	92	58	11	17*	1	5	15*	31
3	15*	16	3	45*	5	12	34*	22	6	354	351
4	124	125	4	56*	45	13	54*	48	7	256	269
5	17*	36	5	19*	23	14	39*	30	8	81	67
6	55	61	6	54*	4	15	30*	1	9	100	92
7	58	41	7	49*	25	16	70	2	10	17*	15
8	17*	15	8	19*	11	17	32*	2	11	41*	13
9	41*	20	9	19*	20	-5,2,L			12	19*	24
10	19*	5	10	19*	3	1	113	117	13	19*	41
11	19*	20	11	134	115	2	79	92	14	19*	30
12	19*	9	12	22*	14	3	58	60	15	22*	6
13	19*	5	13	22*	2	4	595	589	16	22*	5
14	32*	8	14	22*	6	5	64	59	-5,5,L		
15	32*	2	-6,10,L			6	96	97	1	162	153
16	22*	6				7	47*	51	2	102	104
-6,7,L			1	39*	11	8	15*	13	3	111	111
1	179	188	2	19*	12	9	43*	51	4	28*	36
2	171	196	3	28*	27	10	17*	1	5	122	121
3	39*	19	4	28*	3	11	17*	24	6	70	66
4	41*	27	5	19*	15	12	19*	3	7	75	71
5	53	33	6	36*	33	13	19*	22	8	51*	63
6	62	52	7	19*	33	14	32*	0	9	17*	8
7	58	49	8	45*	16	15	22*	7	10	17*	18
8	17*	24	9	19*	19	16	22*	10	11	54*	32
9	19*	5	10	22*	17	17	22*	1	12	19*	36
			11	70	11				13	19*	11
			12	22*	25				14	19*	2

-5,5,L			12	19*	9	5	81	67	-4,4,L		
15	19*	4	13	22*	22	6	13*	5	1	64	53
16	22*	2	14	49*	9	7	15*	16	2	43	47
-5,6,L			15	36*	12	8	15*	5	3	19*	15
			-5,9,L			9	32*	26	4	58	59
1	58	44	1	17*	19	10	49*	3	5	55	50
2	36*	28	2	51*	25	11	17*	26	6	273	272
3	32*	34	3	30*	8	12	19*	18	7	333	333
4	119	114	4	51*	15	13	19*	12	8	60	59
5	39*	32	5	41*	2	14	19*	2	9	90	83
6	51*	43	6	56*	53	15	19*	6	10	56*	37
7	26*	2	7	30*	17	16	19*	3	11	19*	24
8	49*	17	8	19*	5	17	22*		12	19*	1
9	45*	26	9	19*	8	-4,2,L			13	58*	9
10	17*	19	10	19*	4	1	222	223	14	43*	28
11	19*	15	11	22*	25	2	328	310	15	19*	24
12	34*	14	12	22*	51	3	11*	5	16	32*	3
13	19*	8	13	22*	9	4	36*	33	-4,5,L		
14	19*	6	14	19*	4	5	34*	41	1	15*	10
15	22*	2	-5,10,L			6	26*	15	2	24*	43
16	19*	11	1	28*	2	7	62	68	3	92	91
-5,7,L			2	66	62	8	45*	31	4	207	198
1	47*	64	3	45*	44	9	15*	9	5	62	51
2	175	176	4	45*	26	10	36*	7	6	111	118
3	32*	44	5	19*	14	11	45*	52	7	156	161
4	15*	19	6	41*	22	12	56*	34	8	15*	11
5	26*	27	7	19*	9	13	19*	10	9	94	93
6	17*	21	8	36*	25	14	51*	46	10	45*	7
7	34*	29	9	39*	21	15	45*	17	11	62	47
8	32*	1	10	19*	24	16	22*	23	12	45*	43
9	19*	22	11	75	47	17	51*	5	13	19*	1
10	19*	20	12	54*	5	-4,3,L			14	19*	1
11	32*	4	13	22*	3	1	260	258	15	22*	3
12	19*	16	-4,0,L			2	403	391	16	30*	4
13	30*	5	2	17*	12	3	53	48	-4,6,L		
14	45*	6	4	335	344	4	85	92	1	32*	39
15	58*	9	6	188	192	5	70	60	2	134	138
-5,8,L			8	92	91	6	130	133	3	15*	1
1	49*	44	10	34*	39	7	15*	9	4	207	210
2	70	67	12	28*	17	8	15*	19	5	126	130
3	36*	22	14	60	20	9	60	63	6	58	47
4	60	57	16	30*	4	10	77	68	7	17*	16
5	17*	6	-4,1,L			11	34*	14	8	17*	28
6	79	65	1	503	504	12	43*	21	9	19*	26
7	43*	40	2	1231	1233	13	54*	50	10	36*	33
8	17*	12	3	41	33	14	56*	57	11	62	42
9	19*	13	4	173	169	15	54*	18	12	19*	19
10	19*	7				16	22*	13	13	19*	15
11	36*	21				17	77	6	14	19*	9
									15	22*	7

-4,6,L			14	19*	4	5	36*	32	-3,5,L		
16	22*	8	-4,10,L			6	141	146	1	13*	10
-4,7,L			1	45*	30	7	143	148	2	24*	8
1	98	84	2	17*	6	8	73	54	3	62	45
2	15*	12	3	60*	59	9	30*	15	4	26*	32
3	36*	30	4	49*	27	10	17*	9	5	134	120
4	34*	23	5	49*	10	11	58	57	6	211	223
5	17*	11	6	17*	19	12	54*	49	7	85	67
6	30*	25	7	19*	29	13	43*	19	8	17*	18
7	17*	7	8	43*	26	14	19*	37	9	43*	26
8	43*	46	9	81	66	15	41*	1	10	36*	35
9	19*	21	10	22*	7	16	19*	10	11	41*	9
10	41*	10	11	47*	30	17	22*	20	12	19*	29
11	19*	12	12	22*	0	-3,3,L			13	19*	5
12	19*	15	13	22*	7	1	75	78	14	41*	8
13	19*	4	-3,0,L			2	81	89	15	22*	8
14	22*	8	1	520	527	3	87	87	16	43*	0
15	22*	7	3	390	379	4	119	131	-3,6,L		
-4,8,L			5	147	143	5	79	82	1	39*	42
1	15*	1	7	476	486	6	51	34	2	15*	1
2	62	50	9	70	57	7	36*	15	3	15*	0
3	62	47	11	30*	10	8	28*	1	4	94	89
4	17*	34	13	19*	7	9	15*	16	5	58	49
5	111	104	15	19*	4	10	60	63	6	15*	12
6	17*	4	17	22*	2	11	19*	29	7	66	55
7	17*	5	-3,1,L			12	19*	29	8	39*	7
8	19*	13	1	659	672	13	41*	33	9	17*	17
9	68	53	3	612	601	14	32*	20	10	17*	16
10	19*	5	4	162	164	15	19*	4	11	79	58
11	34*	18	5	162	168	16	56*	25	12	32*	12
12	32*	6	6	15*	32	17	30*	19	13	28*	7
13	19*	6	7	134	127	-3,4,L			14	19*	5
14	54*	2	8	75	79	1	151	155	15	22*	5
15	22*	6	9	15*	21	2	156	154	16	22*	8
-4,9,L			10	26*	14	3	39*	34	-3,7,L		
1	55	24	11	94	77	4	43*	38	1	15*	24
2	17*	7	12	19*	28	5	132	134	2	36*	34
3	17*	19	13	19*	10	6	107	101	3	15*	9
4	47*	42	14	19*	5	7	43*	53	4	15*	25
5	19*	46	15	45*	1	8	17*	17	5	66	78
6	17*	21	16	34*	1	9	32*	38	6	68	58
7	17*	7	17	22*	5	10	17*	33	7	39*	57
8	19*	19	-3,2,L			11	17*	20	8	47*	17
9	58*	46	1	143	140	12	34*	9	9	77	58
10	19*	25	2	147	131	13	30*	4	10	79	39
11	22*	26	3	149	135	14	39*	12	11	19*	3
12	47*	14	4	13*	4	15	22*	2	12	49*	3
13	22*	4				16	22*	11	13	19*	12
									14	19*	4
									15	22*	4

-3,8,L			6	43*	32	7	198	201	4	53	50
1	15*	13	8	15*	21	8	55	53	5	77	67
2	28*	11	10	26*	25	9	15*	15	6	15*	14
3	17*	35	12	17*	1	10	17*	17	7	98	92
4	15*	13	14	19*	3	11	73	63	8	17*	3
5	66	66	16	19*	3	12	39*	42	9	79	56
6	43*	42				13	32*	6	10	51*	37
7	51*	29	-2,1,L			14	47*	30	11	36*	31
8	45*	30	1	30	26	15	36*	5	12	43*	8
9	83	81	2	1258	1262	16	22*	35	13	19*	11
10	32*	46	3	47	61	17	58*	20	14	22*	14
11	19*	0	4	105	109				15	22*	12
12	36*	13	5	160	163	-2,4,L			16	43*	5
13	22*	0	6	200	194	1	62	57	-2,7,L		
14	22*	17	7	98	93	2	117	115			
15	22*	2	8	51	54	3	55	54	1	32*	18
			9	77	76	4	32*	15	2	32*	24
-3,9,L			10	55	39	5	92	103	3	28*	9
1	17*	18	11	17*	7	6	119	123	4	15*	13
2	17*	7	12	41*	13	7	81	71	5	26*	15
3	17*	12	13	19*	3	8	28*	6	6	26*	7
4	190	170	14	19*	9	9	17*	6	7	32*	9
5	39*	2	15	19*	1	10	17*	5	8	39*	6
6	19*	1	16	19*	7	11	30*	9	9	56*	61
7	19*	11	17	28*	8	12	17*	9	10	19*	14
8	19*	7				13	19*	18	11	47*	33
9	19*	43	-2,2,L			14	30*	7	12	19*	1
10	75	60	1	232	235	15	22*	5	13	19*	2
11	19*	33	2	177	173	16	22*	23	14	22*	26
12	22*	12	3	399	394	-2,5,L			15	22*	1
13	22*	8	4	53	58	1	36*	12	-2,8,L		
14	22*	1	5	77	74	2	341	328	1	15*	4
-3,10,L			6	15*	13	3	190	189	2	49*	50
1	17*	13	7	34*	17	4	435	448	3	17*	16
2	19*	43	8	32*	24	5	363	373	4	30*	41
3	30*	29	9	15*	3	6	68	61	5	36*	21
4	45*	44	10	17*	18	7	32*	7	6	28*	36
5	19*	35	11	43*	4	8	32*	30	7	19*	39
6	32*	4	12	107	124	9	17*	18	8	32*	23
7	43*	0	13	19*	6	10	34*	24	9	58*	51
8	19*	26	14	19*	9	11	66	36	10	19*	42
9	105	93	15	56*	2	12	19*	4	11	19*	2
10	68	36	16	43*	13	13	19*	12	12	19*	2
11	22*	0	17	22*	10	14	19*	1	13	41*	7
12	22*	8	-2,3,L			15	34*	4	14	43*	13
13	22*	2	1	380	385	16	19*	3	15	22*	10
-2,0,L			2	350	340	-2,6,L			-2,9,L		
			3	13*	18	1	36*	39	1	17*	4
2	1325	1292	4	34*	29	2	49	40	2	17*	19
4	245	248	5	43*	0	3	87	79			
			6	81	81						



-2,9,L			11	17*	14	6	15*	7	4	34*	12
3	47*	40	12	70	78	7	145	149	5	15*	20
4	64	47	13	19*	18	8	105	90	6	28*	2
5	73	65	14	34*	8	9	60	66	7	68	38
6	17*	19	15	22*	9	10	17*	13	8	43*	7
7	34*	17	16	19*	12	11	28*	1	9	45*	30
8	19*	3	17	22*	9	12	19*	13	10	62	28
9	19*	9				13	19*	15	11	19*	13
10	32*	3	-1,2,L			14	19*	7	12	19*	26
11	49*	17	1	350	356	15	19*	3	13	39*	6
12	45*	2	2	70	62	16	22*	10	14	22*	14
13	30*	20	3	81	83	-1,5,L			15	22*	8
14	32*	2	4	254	263	1	51	45	-1,8,L		
-2,10,L			5	13*	6	2	286	282	1	43*	1
1	17*	12	6	81	79	3	166	169	2	34*	32
2	58	3	7	68	76	4	301	316	3	17*	22
3	17*	7	8	22*	37	5	431	430	4	17*	33
4	54*	59	9	66	75	6	30*	3	5	51*	32
5	19*	17	10	17*	39	7	15*	27	6	17*	10
6	19*	16	11	117	115	8	53	33	7	54*	49
7	17*	16	12	87	85	9	17*	16	8	19*	18
8	19*	1	13	19*	23	10	17*	4	9	54*	42
9	96	86	14	34*	12	11	30*	37	10	19*	10
10	87	45	15	19*	14	12	49*	6	11	45*	14
11	19*	3	16	58*	6	13	19*	16	12	47*	3
12	22*	4	17	22*	16	14	19*	5	13	45*	7
13	54*	1	-1,3,L			15	22*	5	14	22*	19
-1,0,L			1	70	72	16	19*	3	-1,9,L		
1	275	268	2	87	76	-1,6,L			1	17*	20
3	382	380	3	22*	4	1	24*	14	2	119	119
5	58	52	4	24*	17	2	26*	4	3	60	45
7	15*	16	5	188	188	3	102	105	4	17*	21
9	15*	18	6	15*	19	4	181	173	5	17*	26
11	17*	26	7	147	158	5	102	100	6	28*	7
13	19*	0	8	32*	45	6	85	83	7	30*	1
15	19*	5	9	28*	31	7	36*	29	8	43*	3
17	41*	15	10	70	65	8	30*	27	9	41*	7
-1,1,L			11	26*	6	9	77	70	10	19*	15
1	346	350	12	58*	67	10	98	78	11	34*	3
2	87	89	13	19*	21	11	19*	1	12	22*	12
3	75	77	14	19*	20	12	19*	10	13	30*	12
4	60	63	15	36*	16	13	19*	1	14	22*	6
5	252	235	16	22*	21	14	22*	3	-1,10,L		
6	60	54	17	22*	9	15	47*	3	1	28*	4
7	141	145	-1,4,L			16	22*	1	2	77	64
8	39*	38	1	149	145	-1,7,L			3	17*	19
9	15*	20	2	205	198	1	64	73	4	19*	29
10	49*	21	3	119	118	2	117	121	5	19*	15
			4	186	179	3	126	129			
			5	77	85						

-1,10,L			11	54*	64	6	36*	34	3	109	101
6	49*	13	12	73	66	7	30*	20	4	28*	25
7	19*	7	13	36*	10	8	30*	8	5	17*	11
8	70	43	14	19*	13	9	39*	10	6	32*	26
9	22*	14	15	19*	11	10	49*	21	7	109	108
10	34*	28	16	22*	13	11	19*	8	8	85	57
11	22*	11	17	51*	5	12	28*	15	9	19*	44
12	34*	6	0,3,L			13	47*	13	10	19*	9
13	22*	4	1	179	173	14	19*	13	11	19*	9
0,0,L			2	26*	19	15	19*	13	12	45*	7
2	205	218	3	60	55	16	22*	2	13	19*	3
4	134	132	4	22*	12	0,6,L			14	22*	22
6	79	80	5	47	41	0	171	167	0,9,L		
8	45*	15	6	75	75	1	15*	16	1	49*	40
10	34*	3	7	81	77	2	162	171	2	32*	45
12	19*	20	8	15*	16	3	15*	25	3	115	117
14	19*	2	9	17*	42	4	209	208	4	70	75
16	22*	13	10	98	91	5	15*	25	5	17*	5
0,1,L			11	17*	13	6	47*	39	6	17*	6
1	158	157	12	19*	29	7	47*	45	7	32*	0
2	196	196	13	19*	3	8	54*	49	8	43*	9
3	58	55	14	19*	35	9	17*	10	9	19*	15
4	70	63	15	30*	20	10	19*	33	10	19*	16
5	544	539	16	22*	6	11	19*	14	11	19*	9
6	39*	44	17	41*	0	12	19*	21	12	22*	20
7	98	88	0,4,L			13	19*	18	13	22*	6
8	51	41	0	122	112	14	19*	1	14	22*	6
9	58	56	1	79	72	15	22*	1	0,10,L		
10	17*	1	2	322	331	16	22*	5	0	41*	14
11	45*	25	3	34*	18	0,7,L			1	17*	2
12	54*	54	4	45*	52	1	15*	15	2	98	86
13	19*	12	5	149	156	2	134	133	3	19*	26
14	19*	13	6	43*	33	3	15*	11	4	30*	7
15	19*	4	7	49	47	4	43*	56	5	19*	4
16	19*	8	8	87	89	5	58	60	6	41*	8
17	22*	5	9	17*	15	6	77	74	7	19*	10
0,2,L			10	17*	24	7	158	162	8	62*	47
0	424	447	11	47*	2	8	94	101	9	22*	36
1	139	146	12	19*	27	9	49*	61	10	19*	5
2	85	85	13	43*	10	10	54*	26	11	22*	11
3	290	278	14	45*	50	11	19*	2	12	22*	2
4	224	222	15	22*	30	12	68	8	13	62*	3
5	171	168	16	32*	4	13	19*	7	1,0,L		
6	68	63	0,5,L			14	19*	13	1	815	824
7	66	55	1	115	115	15	22*	7	3	194	202
8	15*	12	2	49	49	0,8,L			5	250	252
9	17*	28	3	102	119	0	17*	10	7	60	65
10	113	110	4	203	209	1	15*	12	9	111	103
			5	151	161	2	94	94			

1,0,L			4	232	220	1,6,L			1,9,L		
11	17*	8	5	288	269	0	75	72	0	47*	48
13	19*	19	6	64	55	1	15*	17	1	17*	14
15	34*	11	7	41*	53	2	286	301	2	79	73
17	22*	3	8	30*	35	3	378	387	3	45*	24
			9	124	122	4	30*	39	4	17*	24
1,1,L			10	83	87	5	79	83	5	17*	21
0	378	393	11	39*	23	6	47*	52	6	17*	3
1	435	452	12	19*	13	7	66	62	7	17*	25
2	198	201	13	43*	20	8	17*	13	8	19*	8
3	181	186	14	19*	24	9	68	65	9	45*	20
4	241	228	15	34*	35	10	19*	14	10	51*	11
5	235	222	16	22*	8	11	39*	20	11	43*	2
6	115	123	17	51*	8	12	19*	26	12	34*	23
7	15*	17	1,4,L			13	19*	0	13	47*	5
8	15*	11	0	599	604	14	19*	7	14	22*	5
9	81	80	1	205	215	15	22*	3	1,10,L		
10	113	106	2	105	101	1,7,L			0	17*	1
11	17*	9	3	73	81	0	58	55	1	17*	23
12	19*	10	4	47	32	1	64	53	2	92	97
13	43*	8	5	79	66	2	15*	20	3	49*	30
14	19*	11	6	15*	14	3	73	67	4	32*	3
15	30*	8	7	36*	28	4	53	44	5	17*	6
16	30*	1	8	81	85	5	36*	41	6	43*	2
17	58*	2	9	26*	35	6	45*	15	7	36*	2
1,2,L			10	28*	6	7	70	67	8	19*	9
0	350	377	11	17*	3	8	105	98	9	22*	20
1	132	143	12	19*	18	9	54*	43	10	19*	4
2	11*	5	13	34*	14	10	28*	14	11	39*	15
3	179	177	14	56*	54	11	19*	4	12	22*	10
4	235	220	15	30*	33	12	36*	13	13	56*	0
5	36*	31	16	47*	1	13	19*	15	2,0,L		
6	117	117	1,5,L			14	54*	6	0	619	643
7	81	76	0	100	86	15	22*	5	2	269	277
8	15*	21	1	13*	5	1,8,L			4	181	180
9	107	106	2	47	43	0	94	89	6	132	134
10	62	56	3	173	175	1	70	65	8	55	53
11	17*	2	4	81	91	2	58	33	10	158	151
12	36*	21	5	41*	45	3	17*	33	12	30*	35
13	51*	3	6	15*	1	4	36*	15	14	62*	45
14	41*	8	7	30*	24	5	17*	11	16	22*	12
15	49*	5	8	15*	22	6	64	45	2,1,L		
16	22*	7	9	28*	24	7	19*	51	0	523	555
17	34*	4	10	28*	10	8	17*	2	1	644	658
1,3,L			11	19*	11	9	19*	12	2	188	188
0	17*	1	12	28*	22	10	19*	30	3	22*	9
1	115	125	13	49*	23	11	19*	2			
2	81	82	14	41*	27	12	41*	25			
3	75	79	15	41*	16	13	19*	6			
			16	51*	10	14	22*	0			



3,1,L			7	85	75	3	51	66	4	17*	4
16	19*	12	8	128	122	4	15*	17	5	19*	2
17	22*	0	9	17*	3	5	162	162	6	19*	20
3,2,L			10	47*	29	6	70	87	7	19*	8
			11	19*	3	7	36*	28	8	34*	42
			12	58*	57	8	73	68	9	19*	8
0	126	141	13	64*	51	9	28*	1	10	22*	18
1	224	230	14	47*	14	10	41*	0	11	22*	12
2	320	319	15	22*	12	11	43*	14	12	22*	1
3	279	269	16	39*	1	12	19*	3	13	22*	2
4	13*	24	3,5,L			13	19*	3	4,0,L		
5	62	51	0	36*	19	14	54*	4	0	45	62
6	34*	17	1	164	166	15	22*	1	2	318	306
7	87	101	2	115	101	3,8,L			4	346	345
8	15*	20	3	166	165	0	205	220	6	26*	9
9	28*	33	4	115	121	1	96	104	8	260	259
10	54*	42	5	47*	31	2	64	71	10	28*	37
11	19*	13	6	109	113	3	34*	21	12	43*	41
12	45*	14	7	26*	27	4	30*	15	14	19*	4
13	51*	23	8	17*	16	5	17*	15	16	19*	4
14	19*	13	9	17*	18	6	28*	32	4,1,L		
15	22*	6	10	34*	7	7	43*	39	0	11*	22
16	22*	4	11	19*	1	8	19*	26	1	107	108
17	54*	0	12	81	58	9	19*	19	2	222	221
3,3,L			13	51*	36	10	36*	1	3	456	458
0	420	421	14	22*	3	11	41*	3	4	13*	1
1	196	206	15	19*	1	12	19*	16	5	13*	22
2	262	265	16	22*	5	13	22*	14	6	36*	19
3	13*	0	3,6,L			14	22*	1	7	15*	21
4	77	82	0	32*	36	3,9,L			8	143	136
5	15*	14	1	30*	9	0	17*	4	9	17*	24
6	60	58	2	115	116	1	17*	21	10	77	81
7	109	116	3	260	271	2	17*	2	11	47*	48
8	64	82	4	81	67	3	17*	13	12	41*	20
9	17*	2	5	207	203	4	17*	24	13	32*	1
10	26*	25	6	124	125	5	41*	34	14	19*	12
11	17*	16	7	54*	49	6	17*	17	15	39*	0
12	19*	0	8	17*	3	7	36*	29	16	56*	3
13	19*	26	9	39*	8	8	19*	6	17	30*	1
14	43*	27	10	19*	35	9	19*	9	4,2,L		
15	19*	8	11	34*	13	10	19*	13	0	162	167
16	22*	2	12	19*	13	11	19*	4	1	196	199
3,4,L			13	60*	20	12	22*	31	2	90	85
0	196	201	14	19*	6	13	43*	7	3	994	962
1	73	76	15	19*	12	3,10,L			4	26*	26
2	47	42	3,7,L			0	58*	46	5	109	105
3	94	81	0	175	176	1	66	61	6	15*	30
4	15*	22	1	32*	9	2	75	56	7	41*	20
5	41*	34	2	111	111	3	19*	37			
6	41*	27									

4,2,L			1	177	179	4,8,L			5	28*	10
8	15*	12	2	98	92	0	41*	55	7	375	372
9	68	63	3	28*	6	1	17*	33	9	117	118
10	17*	30	4	15*	18	2	17*	18	11	58*	39
11	36*	10	5	98	96	3	122	122	13	81	69
12	19*	10	6	156	157	4	17*	10	15	34*	2
13	19*	25	7	28*	28	5	85	78	5,1,L		
14	19*	1	8	41*	25	6	49*	36	0	64	65
15	19*	5	9	30*	41	7	92	75	1	117	113
16	56*	5	10	70	72	8	19*	0	2	320	316
4,3,L			11	58*	44	9	19*	13	3	245	241
0	188	192	12	87	75	10	19*	4	4	290	269
1	24*	29	13	81	51	11	19*	11	5	102	94
2	24*	16	14	19*	7	12	19*	3	6	64	69
3	224	229	15	22*	10	13	51*	4	7	132	127
4	64	56	16	47*	3	14	22*	5	8	301	295
5	107	113	4,6,L			4,9,L			9	17*	34
6	15*	19	0	77	80	0	39*	27	10	17*	5
7	90	95	1	250	250	1	17*	22	11	17*	5
8	85	82	2	51	36	2	47*	45	12	51*	37
9	17*	9	3	15*	16	3	17*	1	13	19*	19
10	17*	6	4	60	53	4	56*	44	14	19*	20
11	17*	5	5	213	206	5	17*	14	15	34*	3
12	19*	11	6	132	120	6	19*	4	16	41*	7
13	19*	7	7	26*	29	7	43*	21	5,2,L		
14	56*	10	8	17*	21	8	19*	24	0	45	49
15	34*	22	9	17*	4	9	19*	21	1	160	178
16	22*	11	10	19*	19	10	43*	10	2	66	57
4,4,L			11	62	35	11	30*	16	3	299	276
0	94	89	12	19*	23	12	19*	2	4	28*	29
1	51	48	13	22*	25	13	34*	4	5	15*	12
2	149	148	14	19*	3	4,10,L			6	51	34
3	85	77	15	22*	11	0	36*	14	7	15*	37
4	68	70	4,7,L			1	17*	6	8	15*	23
5	15*	7	0	145	137	2	17*	12	9	49*	30
6	81	80	1	305	307	3	17*	15	10	66	46
7	141	149	2	90	78	4	51*	0	11	36*	29
8	62	69	3	30*	32	5	58*	10	12	54*	14
9	36*	19	4	15*	9	6	47*	24	13	19*	16
10	41*	24	5	85	83	7	19*	11	14	32*	2
11	64	42	6	17*	13	8	19*	3	15	19*	5
12	54*	35	7	17*	6	9	19*	8	16	19*	1
13	66	53	8	49*	23	10	22*	26	5,3,L		
14	22*	20	9	43*	17	11	19*	8	0	81	77
15	22*	17	10	19*	5	12	51*	10	1	119	111
16	22*	3	11	19*	5	5,0,L			2	207	203
4,5,L			12	19*	17	1	388	387	3	177	174
0	55	47	13	66	10	3	154	154	4	41*	41
			14	22*	3						
			15	22*	1						

5,3,L			5,6,L			5,9,L			8	17*	12
5	143	133	0	96	94	0	17*	20	9	17*	7
6	36*	35	1	51	72	1	17*	25	10	49*	52
7	39*	21	2	60	51	2	51*	17	11	43*	31
8	58	61	3	288	275	3	17*	24	12	19*	18
9	17*	11	4	105	104	4	64	55	13	19*	23
10	60	14	5	177	179	5	17*	25	14	19*	4
11	54*	30	6	39*	16	6	49*	53	15	19*	6
12	19*	7	7	47*	16	7	19*	17	16	19*	3
13	19*	10	8	28*	11	8	19*	17	6,2,L		
14	39*	13	9	28*	3	9	19*	1	0	346	357
15	22*	13	10	19*	2	10	75	58	1	26*	31
16	22*	5	11	19*	24	11	34*	8	2	102	99
5,4,L			12	19*	16	12	19*	1	3	171	160
0	200	191	13	19*	4	13	51*	5	4	92	85
1	41*	39	14	22*	4	5,10,L			5	22*	1
2	181	176	15	22*	11	0	45*	17	6	28*	1
3	98	94	5,7,L			1	41*	13	7	15*	18
4	87	84	0	296	299	2	19*	7	8	79	63
5	64	67	1	309	321	3	19*	25	9	36*	42
6	15*	32	2	36*	25	4	43*	16	10	43*	4
7	58	56	3	24*	32	5	32*	30	11	17*	1
8	55	60	4	55	52	6	19*	21	12	32*	10
9	17*	18	5	17*	11	7	19*	6	13	34*	31
10	34*	8	6	17*	25	8	30*	12	14	19*	3
11	19*	16	7	68	60	9	22*	1	15	54*	1
12	19*	15	8	17*	6	10	34*	18	16	30*	1
13	19*	16	9	19*	2	11	22*	16	6,3,L		
14	41*	1	10	17*	4	12	54*	1	0	209	209
15	19*	2	11	19*	9	6,0,L			1	437	444
16	47*	13	12	47*	1	0	847	850	2	94	100
5,5,L			13	34*	18	2	124	120	3	143	132
0	13*	11	14	22*	2	4	171	161	4	81	69
1	15*	23	15	22*	15	6	83	96	5	73	81
2	98	105	5,8,L			8	433	442	6	81	79
3	30*	16	0	15*	27	10	26*	0	7	24*	16
4	15*	11	1	62	61	12	32*	1	8	41*	44
5	15*	4	2	15*	1	14	47*	23	9	26*	31
6	83	81	3	26*	42	16	43*	12	10	17*	13
7	17*	0	4	17*	1	6,1,L			11	19*	1
8	36*	41	5	60	51	0	83	88	12	51*	16
9	47*	9	6	132	111	1	282	283	13	66	37
10	17*	13	7	17*	12	2	75	80	14	22*	0
11	19*	5	8	19*	19	3	288	271	15	19*	3
12	81	57	9	19*	1	4	32*	20	16	22*	1
13	47*	32	10	41*	30	5	186	184	6,4,L		
14	22*	9	11	45*	24	6	36*	33	0	68	74
15	22*	1	12	19*	0	7	66	60			
16	36*	3	13	19*	4						
			14	22*	6						





7,5,L			7,8,L			4			13*			9			5			64		
0	15*	13	0	79	71	6	653	666	6	15*	26	6	15*	26	6	15*	26	6	15*	26
1	15*	4	1	15*	7	8	73	75	7	34*	41	7	34*	41	7	34*	41	7	34*	41
2	45*	46	2	15*	31	10	49*	8	8	17*	9	8	17*	9	8	17*	9	8	17*	9
3	341	353	3	17*	27	12	34*	51	9	51*	54	9	51*	54	9	51*	54	9	51*	54
4	113	105	4	39*	15	14	41*	20	10	17*	2	10	17*	2	10	17*	2	10	17*	2
5	96	98	5	17*	9	16	22*	7	11	19*	20	11	19*	20	11	19*	20	11	19*	20
6	122	126	6	17*	12	8,1,L			12	19*	2	12	19*	2	12	19*	2	12	19*	2
7	15*	1	7	39*	19	0	279	279	13	19*	14	13	19*	14	13	19*	14	13	19*	14
8	17*	34	8	39*	31	1	599	582	14	19*	4	14	19*	4	14	19*	4	14	19*	4
9	49*	42	9	19*	28	2	250	235	15	19*	7	15	19*	7	15	19*	7	15	19*	7
10	58*	49	10	49*	16	3	24*	1	16	22*	12	16	22*	12	16	22*	12	16	22*	12
11	81	85	11	45*	13	4	22*	30	8,4,L			0	98	105	0	98	105	0	98	105
12	41*	1	12	19*	3	5	117	122	1	51	46	1	51	46	1	51	46	1	51	46
13	45*	1	13	22*	2	6	277	289	2	124	125	2	124	125	2	124	125	2	124	125
14	19*	4	14	22*	9	7	145	137	3	79	85	3	79	85	3	79	85	3	79	85
15	22*	1	7,9,L			8	39*	41	4	92	93	4	92	93	4	92	93	4	92	93
7,6,L			0	17*	8	9	17*	15	5	15*	33	5	15*	33	5	15*	33	5	15*	33
0	49*	55	1	70	60	10	34*	2	6	85	94	6	85	94	6	85	94	6	85	94
1	111	115	2	17*	6	11	70	50	7	15*	13	7	15*	13	7	15*	13	7	15*	13
2	30*	17	3	49*	49	12	19*	17	8	58	66	8	58	66	8	58	66	8	58	66
3	177	177	4	54*	48	13	19*	3	9	17*	26	9	17*	26	9	17*	26	9	17*	26
4	200	205	5	56*	33	14	19*	5	10	28*	12	10	28*	12	10	28*	12	10	28*	12
5	81	80	6	51*	16	15	19*	2	11	17*	5	11	17*	5	11	17*	5	11	17*	5
6	43*	22	7	19*	6	8,2,L			12	28*	7	12	28*	7	12	28*	7	12	28*	7
7	49*	30	8	19*	4	0	66	55	13	19*	3	13	19*	3	13	19*	3	13	19*	3
8	32*	29	9	41*	19	1	203	206	14	51*	12	14	51*	12	14	51*	12	14	51*	12
9	30*	35	10	19*	12	2	158	157	15	22*	9	15	22*	9	15	22*	9	15	22*	9
10	81	82	11	22*	1	3	117	117	8,5,L			0	39*	20	0	39*	20	0	39*	20
11	70	56	12	22*	0	4	45*	44	1	250	251	1	250	251	1	250	251	1	250	251
12	19*	16	13	22*	5	5	160	165	2	49*	40	2	49*	40	2	49*	40	2	49*	40
13	19*	15	7,10,L			6	149	148	3	181	176	3	181	176	3	181	176	3	181	176
14	19*	3	0	62	34	7	36*	39	4	126	123	4	126	123	4	126	123	4	126	123
15	22*	14	1	17*	10	8	15*	26	5	75	83	5	75	83	5	75	83	5	75	83
7,7,L			2	36*	34	9	17*	14	6	39*	48	6	39*	48	6	39*	48	6	39*	48
0	15*	1	3	30*	43	10	58	35	7	17*	8	7	17*	8	7	17*	8	7	17*	8
1	60	47	4	47*	34	11	68	44	8	17*	31	8	17*	31	8	17*	31	8	17*	31
2	15*	6	5	32*	12	12	19*	17	9	77	48	9	77	48	9	77	48	9	77	48
3	49*	52	6	19*	17	13	19*	23	10	19*	12	10	19*	12	10	19*	12	10	19*	12
4	17*	11	7	19*	23	14	19*	0	11	19*	24	11	19*	24	11	19*	24	11	19*	24
5	17*	23	8	43*	24	15	51*	2	12	41*	14	12	41*	14	12	41*	14	12	41*	14
6	54*	39	9	22*	14	16	19*	11	13	43*	11	13	43*	11	13	43*	11	13	43*	11
7	54*	36	10	56*	18	8,3,L			14	45*	21	14	45*	21	14	45*	21	14	45*	21
8	19*	10	11	22*	2	0	51	54	15	22*	2	15	22*	2	15	22*	2	15	22*	2
9	47*	34	12	22*	0	1	77	77	8,0,L											
10	19*	15				2	13*	20	0	117	112	0	117	112	0	117	112	0	117	112
11	41*	23				3	34*	29	1	124	127	1	124	127	1	124	127	1	124	127
12	30*	12				4	49	42												
13	22*	2																		
14	19*	5																		

8,6,L			8,9,L			9			4		
0	41*	26	0	49*	21	10	34*	43	4	26*	8
1	68	63	1	17*	21	11	49*	30	5	15*	3
2	62	52	2	17*	34	12	47*	52	6	134	144
3	24*	46	3	107	105	13	19*	17	7	64	64
4	105	109	4	43*	36	14	19*	2	8	47*	29
5	47*	52	5	17*	13	15	32*	9	9	26*	15
6	66	57	6	49*	14	16	19*	17	10	17*	2
7	34*	14	7	19*	2		22*	16	11	19*	20
8	17*	42	8	19*	20	9,2,L			12	19*	5
9	39*	2	9	54*	21	0	243	245	13	19*	6
10	126	99	10	19*	13	1	55	59	14	19*	7
11	19*	35	11	19*	0	2	205	202	15	22*	13
12	19*	3	12	30*	3	3	111	112	9,5,L		
13	22*	18	13	43*	1	4	124	127	0	90	84
14	22*	14				5	51	53	1	292	316
15	22*	0	8,10,L			6	39*	10	2	264	266
			0	17*	15	7	26*	22	3	241	258
8,7,L			1	17*	22	8	51*	46	4	124	138
0	53	61	2	39*	22	9	17*	3	5	55	49
1	34*	33	3	17*	6	10	85	80	6	30*	14
2	30*	27	4	64	56	11	19*	13	7	30*	27
3	15*	4	5	19*	22	12	19*	4	8	17*	19
4	83	82	6	30*	12	13	19*	12	9	41*	47
5	60	50	7	19*	7	14	34*	0	10	45*	47
6	47*	32	8	45*	36	15	22*	1	11	41*	15
7	17*	5	9	39*	16	16	22*	6	12	19*	2
8	19*	26	10	41*	5	9,3,L			13	19*	18
9	70	52	11	49*	2	0	75	81	14	22*	11
10	45*	19	12	22*	1	1	102	103	15	22*	1
11	19*	3	9,0,L			2	55	68	9,6,L		
12	19*	17	1	576	559	3	26*	14	0	51	62
13	43*	36	3	94	97	4	53	55	1	87	82
14	19*	11	5	79	73	5	15*	15	2	130	124
8,8,L			7	15*	19	6	30*	35	3	43*	57
0	15*	2	9	70	72	7	58	60	4	15*	27
1	26*	6	11	49*	4	8	49*	33	5	15*	24
2	15*	9	13	19*	7	9	17*	28	6	58	46
3	17*	10	15	22*	1	10	17*	6	7	17*	9
4	66	82				11	19*	34	8	75	82
5	17*	30	9,1,L			12	19*	2	9	81	78
6	34*	3	0	143	144	13	19*	4	10	19*	18
7	19*	25	1	397	384	14	22*	1	11	19*	11
8	45*	41	2	115	115	15	22*	14	12	43*	19
9	19*	7	3	49	30	16	22*	12	13	51*	26
10	19*	3	4	32*	27	9,4,L			14	19*	6
11	19*	25	5	262	270	0	30*	12	15	22*	12
12	19*	3	6	262	266	1	107	118	9,7,L		
13	22*	14	7	17*	23	2	85	114	0	43*	26
14	19*	3	8	17*	33	3	43*	38			

9,7,L			2	19*	18	9	17*	32	6	49*	15
1	15*	16	3	56*	33	10	17*	12	7	34*	16
2	28*	23	4	17*	8	11	96	74	8	43*	36
3	15*	14	5	19*	5	12	36*	21	9	39*	25
4	15*	11	6	51*	18	13	19*	5	10	75	69
5	34*	13	7	19*	15	14	32*	11	11	41*	18
6	36*	30	8	64*	37	15	19*	7	12	19*	16
7	30*	16	9	30*	20	16	22*	12	13	19*	6
8	96	76	10	39*	7				14	22*	9
9	19*	33	11	22*	2	10,3,L			15	22*	5
10	19*	2	12	22*	4	0	124	129	10,6,L		
11	19*	4	10,0,L			1	13*	19	0	15*	25
12	22*	3	0	19*	23	2	15*	20	1	34*	18
13	22*	3	2	34*	24	3	49*	60	2	100	93
14	32*	16	4	85	88	4	107	90	3	58	67
9,8,L			6	115	122	5	87	82	4	100	105
0	17*	16	8	17*	26	6	130	139	5	36*	37
1	58	49	10	43*	38	7	17*	6	6	17*	18
2	17*	26	12	19*	6	8	39*	8	7	17*	14
3	17*	9	14	19*	7	9	17*	3	8	75	69
4	17*	25	16	51*	33	10	32*	17	9	126	113
5	28*	22	10,1,L			11	19*	6	10	19*	15
6	17*	18	0	55	59	12	19*	5	11	47*	24
7	19*	29	1	188	191	13	19*	2	12	19*	5
8	19*	21	2	13*	16	14	19*	9	13	22*	11
9	19*	32	3	171	191	15	43*	4	14	56*	9
10	56*	16	4	271	287	10,4,L			10,7,L		
11	19*	14	5	47*	43	0	151	151	0	53	56
12	22*	6	6	188	197	1	407	423	1	49*	47
13	32*	6	7	87	81	2	160	172	2	154	155
9,9,L			8	34*	23	3	85	88	3	17*	6
0	17*	16	9	26*	15	4	43*	40	4	17*	11
1	17*	2	10	36*	25	5	15*	13	5	17*	14
2	17*	9	11	66	56	6	49*	54	6	39*	11
3	17*	20	12	19*	13	7	17*	12	7	34*	7
4	79	88	13	19*	10	8	17*	7	8	100	83
5	19*	15	14	22*	7	9	17*	24	9	60*	64
6	41*	4	15	22*	17	10	17*	3	10	34*	5
7	30*	1	16	49*	14	11	32*	4	11	34*	31
8	19*	14	10,2,L			12	28*	14	12	19*	3
9	19*	7	0	66	73	13	19*	6	13	22*	20
10	19*	1	1	51	49	14	19*	12	14	22*	19
11	22*	1	2	64	68	15	22*	13			
12	22*	7	3	53	45	10,5,L			10,8,L		
13	22*	5	4	45*	42	0	53	52	0	15*	8
9,10,L			5	58	44	1	333	334	1	34*	20
0	17*	11	6	15*	22	2	145	147	2	17*	25
1	19*	29	7	68	71	3	34*	8	3	66	58
			8	64	48	4	190	196	4	17*	6
						5	15*	6			
10,8,L			10,9,L			9	19*	18	4	19*	12
5	17*	8	0	17*	14	10	58*	6	5	34*	2
6	19*	10	1	60	63	11	22*	10	6	39*	9
7	19*	33	2	36*	42	12	41*	9	7	19*	19
8	49*	23	3	17*	15	10,10,L			8	54*	3
9	47*	15	4	17*	7	0	17*	4	9	19*	2
10	19*	3	5	41*	25	1	19*	43	10	22*	5
11	19*	18	6	34*	3	2	19*	35	11	22*	3
12	41*	4	7	49*	25	3	19*	31			
13	30*	17	8	36*	2						

VITA 2

Chinda Wongwiechintana

Candidate for the Degree of

Doctor of Philosophy

Thesis: OXIDATION-REDUCTION REACTIONS OF MORPHINE IN  
DIFFERENT MEDIA

Major Field: Chemistry

Biographical:

Personal Data: Born in Bangkok, Thailand, September  
23, 1951, the daughter of Mr. and Mrs. Wong.

Education: Graduated from Saynumpung High School,  
Bangkok, Thailand, in March, 1970; received the  
Bachelor of Science degree from Mahidol Univer-  
sity, Bangkok, Thailand, with a major in  
Chemistry in March, 1974; received the Master  
of Science degree from Oklahoma State University,  
with a major in Natural Science in July, 1977;  
completed the requirements for the Doctor of  
Philosophy degree at Oklahoma State University  
in December, 1983.

Professional Experience: Graduate teaching assistant,  
Chemistry Department, Oklahoma State University  
from 1975 to 1977 and from 1979 to 1983; taught  
at Srinakarintarawit University, Bangkok,  
Thailand from 1977 to 1979.